

Welch Allyn®

KleenSpec® Cordless Illumination System

The Welch Allyn Cordless Illumination System is a unique design with enhanced features including bright, white LED light to help improve patient exams. The compact, efficient system offers on-time of 80 minutes that will get you through the day, and recharges in 6 hours. Standardise your facility with KleenSpec Cordless Illumination Systems and KleenSpec Disposable Specula to help please both your staff and your patients

- For use with the KleenSpec Disposable Vaginal Specula
- No bulb replacement necessary, for reduced cost of ownership
- Compact LED light source with homogenous spot improves exam area visualisation
- No cord issues and on-time of 80 minutes, to help improve exam efficiency
- Automatic illumination when specula is attached / Once removed light goes off
- Produces much less heat than traditional lamps, to help improve comfort

WA80015 **\$345**

Complete KleenSpec Cordless Illuminator System with Charging Station



KleenSpec® Vaginal Specula
(590 Series Kleenspec - premium series)

Small	WA59000	Box/24	\$142.80	4+	\$129.65
Medium	WA59001	Box/24	\$142.80	4+	\$129.65
Large	WA59004	Box/18	\$129.95	4+	\$119.25



Vaginal Specula Disposable

GRAVES PLASTIC (ratchet type)

Small	ARTA300	Box/24	\$24.10	6+	\$21.25
Medium	ARTA310	Box/24	\$24.10	6+	\$21.25
Large	ARTA320	Box/20	\$22.50	6+	\$20.30



See page 384-385 for more gynaecology instruments & women's health

Proctoscopes Disposable

PLASTIC CLEAR

Small	WWZZP008	Each	\$8.65	10+	\$7.95
Medium	WWZZP009	Each	\$10.05	10+	\$8.65
Large	WWZZP010	Each	\$15.50	10+	\$13.35



Algerbrush II with 1mm burr

The worldwide standard for corneal rust ring removers. Dependable, versatile and easy to use

Foreign bodies in the eye can leave rust rings when removed. These rust rings are normally caused by a ferrous foreign body accidentally striking the cornea or sclera of the patient. When the foreign body (usually from a spark produced by a grinding wheel or welding torch) strikes the eye it sticks to the eye and eventually produces rust when combined with the moisture in the eye. The ALGERBRUSH II is used for the expressed purpose of removing these rust rings from the eye of a patient.

The ALGERBRUSH II has a very low torque motor powered by a single AA battery. The burr in this rotary instrument is used as a "brush" rather than a drill and the rust ring is thus "brushed" from the cornea or sclera leaving a smooth surface which in turn heals much faster than with other instruments which may be used for this purpose.



Algerbrush II with 1 mm burr	WWBR21	\$169.00
Spare Burr 1.0mm	AGRBU1U	\$22.00
Spare Burr 0.5mm	AGRBU5U	\$22.00
Chuck only	AGRCHUCK2	\$33.50

Medline Surgical Clipper & Blades



- Cordless, rechargeable surgical clipper features an ergonomic handle for greater comfort for the user and patient
- Helps maintain skin integrity with all moving blades away from the patients skin
- Easy to use and clean - Clippers are easy to attach and remove, and the clipper handle is fully submersible
- Nickel-metal hydride battery runs long for continuous maximum performance
- Table mount-in charger for upright charging / storing on table counter
- Wall plate / mounting included
- Three interchangeable blade attachments available



Surgical Clipper (clipper only; please order charger and blades separately)	MDL70801	Each	\$98.60
Charger & Base	MDL70802	Each	\$49.35
Blade, clipper, universal, gray Designed for all types of hair. Recommended for general surgery, cardiology and orthopaedics	MDL70880	Pkt/50	\$229.00
Blade, clipper, coarse hair, blue Designed for head hair and other very thick hair. Recommended for neurosurgery	MDL70885	Pkt/50	\$793.90
Blade, clipper, narrow blade, magenta Designed for hair removal in the groin and perineal area. Recommended for urology, labour and delivery, obstetrics and gynaecology	MDL70890	Pkt/50	\$250.00



3M™ Professional 9681 Surgical Clipper

- Blade helps prevent skin from being pulled into the blade
- Blade positioning minimises the potential for nicks or cuts
- Removes all types of hair in a single pass
- 120 minutes of use time on a 4-hour recharge



Surgical Clipper (clipper only; please order charger and blades separately)	3M9681	Each	\$215.00
Charger & Base with Cord	3M9685	Each	\$98.35
Single-Use Surgical Clipper Blades	3M9680	Box/50	\$347.00

Also available

3M™ Single-Use Surgical Clipper Blades for
3M™ Surgical Clipper with Pivoting Head - 9661L

3M9660 Box/50 \$359.90



LOOKING FOR ANOTHER PRODUCT?
GO ONLINE AND SEARCH OUR ENTIRE RANGE
www.teammed.com.au

Razors

Single Edged

A comb that facilitates the shaving of long hair

XXB006

Pkt/12 \$9.70
2+ \$9.05



Razors

Single Edged

SPMMH0812

Box/100 \$36.60
3+ \$34.05



Razors

Double Edged

SPMMH0806

Box/100 \$59.65
3+ \$55.50




Biopsy punches

- Sterile, single use
- Seamless stainless steel sharp cutting edge
- Ribbed handle for sure grip and control
- Sterile and ready to use

2mm	SPMKAI700
3mm	SPMKAI701
4mm	SPMKAI702
5mm	SPMKAI704
6mm	SPMKAI705
8mm	SPMKAI706





Box/20	\$84.70
3+	\$78.85
4+	\$71.40



Biopsy punch with plunger system

- Sterile, single use
- Seamless stainless steel sharp cutting edge
- Internal plunger system
- Ejects lodged skin specimen inside the metal lumen of the punch



Box/20	\$93.25	2mm	SPMKAI697
3+	\$86.75	3mm	SPMKAI698
		4mm	SPMKAI699




Biopsy blade

- Sterile and ready to use
- Blade with sure grip handle for precise biopsy



SPMKAI610

Box/20	\$56.70
3+	\$53.90



Dermal curettes

- Stainless steel semi sharp cutting edge for soft or weak tissue
- Size embossed on each curette for easy identification
- Used in curettage of basal cell carcinoma, verruca or other benign neoplasms of the skin



Box/20	\$120.85	2mm	SPMKAI402
3+	\$112.50	3mm	SPMKAI403
		4mm	SPMKAI404
		5mm	SPMKAI405
		7mm	SPMKAI407



Stiefel Biopsy Punches

- Biopsy Punch 2mm STI1B002
- Biopsy Punch 3mm STI1B003
- Biopsy Punch 4mm STI1B004
- Biopsy Punch 5mm STI1B005
- Biopsy Punch 6mm STI1B006



Box/10 \$68.90 3+ \$62.95



Stiefel Dermal Curette

- Dermal Curette 4mm STIIC004
- Dermal Curette 7mm STISF1C007

Box/10 \$88.40 3+ \$78.90

Acu-Punch® Biopsy Punches

- Disposable Sterile biopsy punches
- Made in the U.S.A.
- Always razor sharp and ready to use
- Handle is ribbed for maximum control, even when wearing gloves

- 2mm BYDPP250
- 3mm BYDPP350
- 4mm BYDPP450
- 5mm BYDPP550
- 6mm BYDPP650 Box/10 \$61.75
- 8mm BYDPP850 5+ \$51.05
- 10mm BYDP1050 Box/10 \$83.45
- 12mm BYDPP1250 5+ \$69.00



Each punch comes individually wrapped in an easy peel pouch

Acu-Dispo-Curette® Dermal Curettes

- Made in the U.S.A.
- Carefully weighted and ergonomically designed to fit the physician's hand
- Each curette comes individually wrapped in an easy peel pouch.

Box/10 \$57.45 3+ \$51.85

- 2mm BYDRRO250
- 3mm BYDRRO350
- 4mm BYDRRO450
- 5mm BYDRRO550
- 7mm BYDRRO750



Bionix Safe Ear Curettes

Bionix Safe Ear Curettes are flexible and make them safer than stainless steel curettes, reducing injury to the ear canal. Safely and gently remove cerumen with styles to match different types of wax and personal preference.



Adult

Yellow CeraSpoon
Perfect for softer wax, especially in patients who produce excessive amounts of cerumen. 4mm
AUR6333 Box/50 \$93.45



Adult

Orange ControLoop
The Orange ControLoop has a bent "Shapleigh" style tip that allows the wax to be pulled straight out of the ear canal. 4mm
AUR2999 Box/50 \$104.85



Adult

Red AngleLoop
The Red AngleLoop has an angled tip that facilitates reaching the hard-to-get, tough cerumen common in older patients. 4mm
AUR3444 Box/50 \$95.45



General

White FlexLoop
The White FlexLoop is an all purpose curette designed for normal curettage in a variety of patients. 4mm
AUR9555 Box/50 \$114.95



Youth

Purple VersaLoop
The Purple VersaLoop has a smooth flexible tip that gives the user exceptional control while offering the patient more comfort. 3mm
AUR4111 Box/50 \$88.70



Youth

Blue InfantScoop
The Blue InfantScoop is designed for smaller ears and may be used through an otoscope. 2mm
AUR4888 Box/50 \$94.10



Youth

Green MicroLoop
The Green MicroLoop is more rigid in construction and is designed to remove tough, impacted cerumen. 3mm
AUR1222 Box/50 \$95.45



Variety Kit

(10/each style & 15/white)
Box/75 \$152.90 AUR5777

The Lighted Ear Curette™ with Magnification

The Lighted Ear Curette with Magnification brilliantly illuminates the ear canal eliminating the practice of blindly curettaging the sensitive area.



Variety Pack
Includes all tip styles except CeraPik
Includes 1 Light Source and 1 Magnification Lens

Box/50 \$189
AUR2245





LIGHTING THROATS AROUND THE WORLD

Throat Scope all-in-one light and tongue depressor allows the clinician one handed control for oral cavity examination. The light is guided by the blade to illuminate the inside of the mouth/throat to provide the clinician with a fast, accurate and pleasant oral cavity examination.

Throat Scope is proudly designed and developed in Australia.

Throat Scope Starter Kit

Includes illuminator + 2 blades **THRTS106** **\$27.95**



See
page
67

Throat Scope Blades

All individually wrapped

Box/50 **THRTS102** **\$35.00**

Box/100 **THRTS103** **\$38.00**



Throat Scope Quick Pack

Includes 1 starter kit and 1 box of blades

THRTHSQP **\$65.00**



Tongue Depressors

Box/100 **\$4.20**

50+ **\$3.60**

AXP10107211



Box/250 **\$9.40**

20+ **\$8.75**

AXP10107212



SMI Surgical Marker Pen

- Double Sided; normal and fine tip
- Includes Ruler
- Sterile

Each **\$8.45**

BYDZSM20 6+ **\$6.90**



Symmetry Surgical Markers


includes ruler (previously Codman)

Symmetry Surgical Marker Pen

JJ431030 Each **\$10.65** 6+ **\$9.40**

Symmetry Surgical Marker Pen Micro

JJ431032 Each **\$10.15** 6+ **\$9.30**




Medi-Safe Safety Scalpel

Top Activated, Retractable, Single Use

- No passing of scalpel with exposed blade
- No attaching or removing of blade from scalpel handle
- Audible "clicks" indicate blade position (retracted, extended and final lock for disposal)
- Peace of mind – no compromise of your professional skills and technique
- Simple, convenient, effective, affordable

SAFETY + SCALPEL™



Sz10 **BYDMSS10** Sz22 **BYDMSS22**

Sz11 **BYDMSS11** Sz23 **BYDMSS23**

Sz15 **BYDMSS15** Sz24 **BYDMSS24**

Box/10 **\$26.65**

3+ **\$23.90**



BP Medical (WA) | www.teammed.com.au
PHONE 1300 33 55 60 | info@bpmedical.com.au

All prices are subject to 10% GST unless marked with * (GST Free)



Team Medical | www.teammed.com.au
PHONE 1300 22 44 50 | info@teammed.com.au



Surgical Blades

- Quality surgical blades
- Made in Japan

Sz10	SPMN010
Sz11	SPMN011
Sz12	SPMN012
Sz15	SPMN015
Sz23	SPMN023



Box/100	\$36.05
3+	\$34.30



Disposable scalpels

- Stainless steel blade
- Easy dispenser box

Sz10	SPMSCALP10
Sz11	SPMSCALP11
Sz15	SPMSCALP15

Box/20	\$28.40
5+	\$26.40



Stitch Cutters

- Quality stitch cutter
- Made in Japan

SPMKAI6999

Box/20	\$37.85
3+	\$36.30



Surgical blades



Sz10	BYD263	Sz21	BYD268
Sz11	BYD264	Sz22	BYD269
Sz12	BYD265	Sz23	BYD270
Sz15	BYD266	Sz24	BYD271
Sz20	BYD267		

Box/100	\$21.45
3+	\$18.70



Disposable scalpels



Sz10	BYD27510
Sz11	BYD27511
Sz15	BYD27515

Box/10	\$12.90
5+	\$11.80



Stitch Cutters



Box/100	\$41.95
3+	\$35.80


BYD274

Scalpel Handles

Handle No 3	BYDA4950	\$14.00	suitable for blades 10-15
Handle No 7	BYDS4954	\$10.95	suitable for blades 10-15
Handle No 4	BYDA4951	\$12.65	suitable for blades 20-24



Swann-Morton®



Blades 10 - 15
suitable for
Handle no 3 or 7

Blades 20 - 24
suitable for
Handle no 4

Surgical Blades

The comprehensive Swann-Morton range of surgical blades and scalpels is the result of dedication to quality of the highest standard and the latest technology.

- Unsurpassed sharpness
- Superior feel
- Exceptional durability and quality



Sz10	SNSM0201
Sz11	SNSM0203
Sz12	SNSM0204
Sz15	SNSM0205
Sz20	SNSM0206
Sz21	SNSM0207
Sz22	SNSM0208
Sz23	SNSM0210
Sz24	SNSM0211

Box/100	\$46.10
3+	\$42.55



Stitch Cutters



	SNSM0420
Box/100	\$46.90
3+	\$44.15



Scalpel Handles



No 3	SNSM0803	\$18.20
No 4	SNSM0804	\$18.20



Disposable Scalpels

Sz10	SNSM0501
Sz11	SNSM0503
Sz12	SNSM0504
Sz15	SNSM0505
Sz20	SNSM0506
Sz21	SNSM0507
Sz22	SNSM0508
Sz23	SNSM0510

Box/10	\$23.00
5+	\$21.20

