

ZOLL AED 3™ Defibrillator

★★★★★ PREMIUM PRODUCT

SMART. READY. TRUSTED.


The AED that leads the way

The ZOLL AED 3™ defibrillator guides rescuers through the process of performing high-quality CPR and can deliver a potentially lifesaving shock to the heart when needed. With a ZOLL AED 3 and quick action, tragedy could be averted.

ZOLL's proprietary Real CPR Help® technology helps guide rescuers in performing high-quality CPR in accordance with the most recent guidelines for depth and rate of compression. The innovative, intuitive design and enhanced features of the ZOLL AED 3 will give unexpected heroes both the confidence and knowledge needed to treat sudden cardiac arrest (SCA).

The ZOLL AED 3 defibrillator provides:

- Enhanced Real CPR Help, featuring a full-colour display with vivid rescue images, CPR cycle timer, and a large colour bar gauge that shows CPR compression depth
- CPR Uni-padz™ electrodes that can be used on both adult and paediatric patients, so there is no need for two separate sets of pads
- Cloud connectivity that can automatically report AED readiness.
- 5-year shelf life for pads and batteries increases readiness and decreases frequency of maintenance.



Leading Technology. Reasonable Cost.

Once installed, the AED 3 costs significantly less than almost any other AED to maintain. Because the CPR Uni-padz can defibrillate both children and adults, there are no special pads to buy and regularly replace to be ready to rescue a child. And because the CPR Uni-padz and the battery both last five years, replacing disposables is easier and less expensive over time. Longer lasting, fewer replacements, lower costs.

No other AED's paediatric pads last 5 years. No other AED lets rescuers rescue a child at the push of a button, and replace pads and battery only once in the first 10 years of service.



AED 3™ Fully automatic

ZOL100120213 **\$3295***

AED 3™ Semi-automatic

ZOL100120113 **\$3195***

AED 3™ BLS Fully automatic

Provides on-screen display of the patient's ECG

ZOL8502120313 **\$3695***

The ZOLL AED™ with RapidShock™ analysis enables the shortest rhythm analysis for more continuous life-saving CPR. RapidShock analysis allows the ZOLL AED 3 to deliver a shock, if needed, in as little as five seconds. Minimising the pre-shock pause and providing more CPR can improve patient outcomes.

ZOLL AED PLUS® Defibrillator

SMART. READY. TRUSTED.

Not all AEDs are created equal

The ZOLL AED Plus® was the first AED with Real CPR Help® to help rescuers know when they are performing high-quality CPR. Today, the ZOLL AED Plus® continues to set the benchmark for design, dependability, durability, and performance.

- Easily recognisable with its distinctive green colour – it won't be missed in an emergency
- Large graphics to guide rescuers through the rescue process
- Fully automatic operation – when a shock is required the device will deliver a shock on its own without the need for the rescuer to press a button
- Clear voice and text prompts to assist you
- Real-time CPR feedback and a metronome to help you provide high-quality CPR
- Longest electrode pad life on the market at 5 years
- Convenient and cost-effective replacement batteries which last for five years even with the device performing a weekly self test
- Integrated Paediatric Capability – when the paediatric pads are connected, the defib automatically proceeds with a paediatric ECG and adjusts energy to paediatric levels
- Built in automatic weekly self-test function making sure the unit is always rescue ready
- Comes complete with a soft carry case, CPR-D padz Adult Electrodes (includes Rescue-Ready Kit that includes CPR barrier mask, gloves, razor, dry towel, antiseptic wipe and scissors), batteries, operators manual, A3 Basic Life Support poster and instructional DVD
- 7-year warranty
- Fully Automatic model also available

★★★★★ PREMIUM PRODUCT



When it's time for CPR, the AED Plus® with Real CPR Help® provides the best support to help save a life.


\$2795*

ZOL102011050

7YR
warranty



Unique one piece
Electrode →



AED Plus features Real CPR Help®, a tool to help the rescuer achieve the correct rate and depth of chest compressions during CPR through visual and voice prompts. Feedback such as "PUSH HARDER", "GOOD COMPRESSIONS" and an adaptive metronome to guide compression rate to help accomplish quality compressions

ZOLL AED Pro® Defibrillator

The ZOLL AED Pro® is designed to perform in any environment. The AED Pro provides the ruggedness, portability and advanced functionality that professional rescuers and services require from an AED.

- Real CPR Help for rate and depth of chest compressions during CPR using ZOLL's unique one-piece CPR-D-padz®
- ECG monitoring with 3-lead cable and high-quality display
- See-thru CPR allows a rescuer to view a patient's underlying ECG rhythm while performing CPR
- Semi-automatic with manual override
- Built-in internal memory and external USB memory technology allows expanded data storage capacity and maximum flexibility for data management and transfer
- Only AED in the industry to withstand a 1.5 metre drop test, making it suitable to withstand rough treatment typically found in an emergency medical environment
- Exceptionally resistant to dust and water jets
- Unrivalled battery options (long-life disposable or rechargeable)
- 5-year warranty



ZOL800AEDPRO

\$5070*

ZOLL Spare Parts and Accessories

Defibrillator Pads for AED Plus & AED Pro


Adult ZOL890008000
\$230.95*

Child ZOL8900081
\$145.95*



Defibrillator Pads for AED 3 CPR Uni-padz™ Universal (adult/paediatric)

ZOL8900000270
\$259.05



Batteries Lithium 123A (AED Plus)

ZOL8000080701
Pkt/10 \$125.10



Battery Pack Lithium AED 3

ZOL800000696
\$261.25



Battery Pack Lithium AED PRO

ZOL8000086001
\$295.00



ILCOR (AED) Wall Signs

Flush Wall Sign
ZOL93100737
\$35.65*

3D Wall Sign
ZOL800001261
\$63.00*



Wall Mount AED Plus

ZOL80000809
\$178.10



Wall Mount Bracket AED 3

ZOL00001266
\$116.75



Wall Cabinet AED 3 Standard Surface

Includes: alarm, set-up instructions and ILCOR label set
ZOL800001256
\$440.00



Intelligent Public Access Defibrillator

IPAD CU-SP1
Intelligent PUBLIC ACCESS DEFIBRILLATOR

The iPAD CU-SP1 - Advanced Features - Advanced Performance

Auto volume adjust

The iPAD CU-SP1 listens to the ambient noise, and automatically adjusts the volume of its prompts to make them clear and easily heard over the background noise.

Smart pads - smart storage!

The pre-connected 'Smart' electrode pads are stored in a clear storage compartment on the underside of the unit. This ensures that the iPAD SP1 is always ready to be used in the shortest amount of time possible. The unique technology built in to the iPAD CU-SP1 and the 'Smart' electrode pads allows the unit to detect, when connected, or if a set of pads has already been used. The visual indicator on the front of the iPAD CU-SP1 indicates if the connected pads are within their life expectancy.

CPR detection

The iPAD CU-SP1 detects if CPR is being performed. If CPR is not being performed, the iPAD CU-SP1 prompts the first responder to 'perform CPR'. If CPR is already being performed, the iPAD CU-SP1 encourages the first responder by prompting them to 'continue CPR'.

Easy communication with cu-ex1 software

Installing CU-EX1 software onto a computer allows you to see and analyse the usage history of the unit. It can store up to 5 events with up to 3 hours of ECG analysis on an SD memory card.

- High Capacity Battery (5 year life span)
- Unit includes Adult Pads
- 5-year warranty

BIOCUSP1 **\$2155***



5YR
warranty



Spare Parts and Accessories IPAD SP1

Defib Pads for IPAD CU-SP1 Adult (Pkt/2)	BIOCUA1007S	\$155.30
Defib Pads for IPAD CU-SP1 Child (Set/2)	BIOCUA1102S	\$223
Battery for IPAD SP1	BIOCUSA1103BB	\$444.50
Wall Bracket for IPAD SP1	BIOCUSP1A12	\$209.80
Carry Case for IPAD SP1	BIOCUSP1BAG	\$282.40

Defibrillator Pads

Have you checked the expiry date of your defibrillator pads lately?



Zoll AED Plus, PRO	Adult	ZOL890008000	Pkt
Zoll AED Plus, PRO	Paediatric	ZOL8900081	Pkt
Zoll AED 3	Adult/Paediatric	ZOL8900000270	Pkt
IPAD NF1200	Adult	BIOCUA0512	Pkt/2
IPAD NF1200	Paediatric	BIOCUA0512P	Pkt/2
IPAD SP1	Adult	BIOCUA1007S	Pkt/2
IPAD SP1	Paediatric	BIOCUA1102S	Pkt/2
Heartstart FR2	Paediatric	LAEM3870A	Pkt/2
Heartstart First Aid	Adult	LAEM5071A	Pkt /2
Heartstart First Aid	Paediatric	LAEM5072A	Pkt /2
Powerheart AED	Adult	CDS9131	Pkt /2
Powerheart AED	Paediatric	CDS9730	Pkt /2
3M Defib pads	11 x 15cm	3M2346	Pkt /2



View Defibrillator Pads prices online