

## Bovie® BANTAM PRO High Frequency Desiccator

The Bovie Bantam PRO is a high frequency desiccator that also has cutting ability. The unit provides 50 watts of cutting power and coagulation modes, plenty of power for most procedures performed in the doctor's office. Power is adjustable in 1/10th watt increments from 0 to 10 watts, for the most accurate output available.



- Large illuminated display
- Self test circuits
- Monopolar & bipolar outputs
- Rotary power dial
- Isolated circuitry
- Digital error detection
- 4-year warranty

Package Includes


- A952 electrosurgical generator
- A902 power control pencil
- A910 disposable sheath x 2
- A910ST sterile disposable sheath x 2
- A804 non-sterile single use sharp dermal tips x 5
- A805 sterile single use sharp dermal tips x 2
- A806 non-sterile single use blunt dermal tips x 5
- A807 sterile single use blunt dermal tips x 2
- A806DE non-sterile single use DERM-Elite premium blunt dermal tips x 2
- A807DE sterile single use DERM-Elite premium blunt dermal tips x 2
- ES01 sterile single use blade electrode x 3
- ES02 sterile single use needle electrode
- ES20 sterile single use ball electrode
- ESRE single use split grounding pad x 5
- A1252C reusable grounding cord
- A837 wall mount kit
- AC power cord
- User guide



**\$4635**

DALA952

Dimensions: 18.8 x 22.8 x 10.5cm



## Bovie® Bipolar Accessories

### Adson Tissue Forceps

- Serrated Tip 127mm [DAL381250A](#) **\$565.00**
- Smooth Tip 127mm [DAL381250B](#) **\$534.30**

### Jewellers Forceps

- 4mm tip 101mm [DAL0800011](#) **\$314.65**

- Bipolar Cord** [DALA827V](#) **\$155.00**

- Footswitch** [DALA803](#) **\$425.00**  
(required with bipolar forceps)




## Bipolar vs. Monopolar Electrosurgical Procedures

**Sharp Electrode Tips:** are most commonly used for procedures requiring pinpoint precision and control.

**Blunt Electrode Tips:** are most commonly used for general blood based coagulation procedures.

**Bipolar Forceps:** are designed to grasp, manipulate and coagulate selected tissue. The electrical current alternates between the two tips, reducing the length of the current path producing a precise, controlled therapeutic effect.



## Bovie® DERM 942 High Frequency Desiccator

The DERM 942 and its handpiece are designed to be a reliable and durable desiccator.

- 40 watt high frequency desiccator
- Adjustable in 1/10th watt increments in the first 10 watts, in monopolar and bipolar
- More accurate power with Universal Power Control 100-240 VAC
- No patient plate required
- Made in the USA
- 4-year warranty

### Package Includes

- A942 electro-surgical generator
- A902 power control pencil
- A910 disposable sheath x 2
- A804 non-sterile single use sharp dermal tips x 5
- A805 sterile single use sharp dermal tips x 2
- A806 non-sterile single use blunt dermal tips x 5
- A807 sterile single use blunt dermal tips x 2
- A806DE non-sterile single use DERM-Elite premium blunt dermal tips x 2
- A807DE sterile single use DERM-Elite premium blunt dermal tips x 2
- A837 wall mount kit
- AC power cord
- User guide

4YR  
warranty

Bovie®



\$2475

DALA942

Dimensions: 18.8 x 22.8 x 10.5cm

## Bovie® Disposable Electrodes



### Sharp Tips

Non-Sterile	DALA804	Box/100	\$140.80
Sterile	DALA805	Box/50	\$140.80



### Blunt Tips

Non-Sterile	DALA806	Box/100	\$140.80
Sterile	DALA807	Box/50	\$140.80

## Bovie® Reusable Electrodes

a	Epilation Needle	DALA834	\$175.50
b	Extended Ball	DALES07R	\$155.00
c	Angled Ball	DALA831	\$108.50



## Bovie® Autoclavable 3-Button Pencil

DALA901


\$455.00



## Bovie® Disposable Handpiece Sheaths

Non-sterile DALA910 Box/100 \$91.55  
Protects the pencil body and cord from external contamination during procedures.

Sterile DALA910ST Box/25 \$256.00  
Protects the handpiece body from external contamination during procedures.



## Bovie® Mounts & Stands

Desk Mount	DALA813	\$140.00
Wall Mount	DALA837	\$73.50
Mobile Stand	DALA812	\$905



# Hyfrecator® 2000



The Hyfrecator has been the #1 product-of-choice by physicians in office-based electrosurgery for almost 70 years

**Features only available in the Hyfrecator® system:**

- High or Low Output Ports; allows you convenient switching back and forth between settings depending on your procedure. Save your most common procedural settings in both high and low mode. Precise output at low and high settings.
- Dual Microprocessors; one microprocessor controls the output and functionality of the system while the second terminates activation if a safety hazard is detected. This makes the Hyfrecator one of the safest office-based electrosurgery units in the world.

1-year warranty

**CONTENTS**

- Power up/down handswitching pencil with 3m cable
- Wall mount kit
- 5 sharp and 5 blunt disposable tips
- Operator's manual, warranty card

The "Common Dermatologic Procedures" video is available online.  
<http://www.vimeo.com/14252858> Password: Electrosurgery



**\$2595**  
 CML7900230

## Electrolase disposable electrodes sharp tips

Sharp Non-Sterile	CML710012BX	Box/100	<b>\$163.40</b>
Sharp Sterile	CML71008BX	Box/50	<b>\$163.40</b>



Sharp Tips - Ideal for pinpoint coagulation

## Electrolase disposable electrodes blunt tips

Blunt Non-Sterile	CML710112BX	Box/100	<b>\$163.40</b>
Blunt Sterile	CML71018BX	Box/50	<b>\$163.40</b>



Blunt Tips - Ideal for broad based coagulation

## Disposable pencil sheaths

Pencil Sheath Non-sterile	CML779618BX	Box/100	<b>\$103.20</b>
Protects the pencil body from external contamination during procedures. 20cm			
Handpiece Sheath Sterile	CML779619CS	Pkt/25	<b>\$292.40</b>
Protects the handpiece body from external contamination during procedures. 92cm			



## Reusable Hyfrecator Accessories



a	CML78097	Bipolar Adson serrated tip forceps	h	CML138019	Short ball electrode (box/5)
b	CML78094	Bipolar Adson smooth tip forceps	i	CML7222L	Long ball electrode
c	CML78091	Bipolar jewellers forceps	j	CML79005	Hand Control Pencil, Autoclavable
d	CML780911	Bipolar cord 12ft	k	CML79006	Footcontrol Pencil (non-switching)
e	CML711B	Electrode general purpose	l	CML79004	Foot Switch (required with bipolar forceps)
f	CML714	Epilation needle (pkt/6)	m	CML79001	Mobile base for hyfrecator
g	CML7222A	Angled ball electrode	n	CML779620	Wall mount kit



Bipolar Starter Pack; Most common starter combination = a + d + i



## Bovie® Smoke Shark™ II Smoke Evacuation System

Bovie

Combines customisable (variable) filter life, a compact shelf design, 3-port filter with filter life tracking technology, while exceeding expectations for value and efficacy.

2YR warranty

- **Procedure Safety**  
AORN recommends the use of smoke evacuation systems during any medical procedure producing surgical smoke
- **Reduced Risk**  
minimises exposure to these toxic smoke plumes during electrosurgery, laser, and harmonic scalpel procedures
- **Variable Filter Life (up to 35 hours)**  
Delivers unsurpassed value and reduces per procedure costs
- **Simple LED Display**  
Two button design (motor speed and on/standby) makes this system extremely easy to operate
- **Quiet Operation**



\$4310

DALSE02220

Dimensions: 27.9 x 39.4 x 15.2cm

Package Includes


- Smoke Shark™ II Smoke Evacuator
- Pneumatic Footswitch
- SP35 Filter
- AC power cord
- User guide

## Bovie® Smoke Shark™ III Smoke Evacuation System

Bovie

- Uniquely designed for the physician's office
- Specifically intended for shorter (<15 minutes) procedures
- Reduces biological, physical and chemical hazards present in electrosurgical smoke plume
- Filter life tracking technology
- Hand-operated activation for maximising filter life

DESIGNED FOR GENERAL PRACTICE



2YR warranty

\$3855

DALSE03220

Dimensions: 23.9 x 15.5 x 18cm

Package Includes

- Smoke Shark™ III Smoke Evacuator
- Smart Shark Automatic Activation Device
- SP06 Filter
- Wall mounting kit
- AC power cord
- User guide

### Bovie® Smoke Shark™ II 3 Port Filter

Four layers of proven filtration gross particulate, ULPA media, activated carbon, and post filter.

3 port 6.4mm, 9.5mm and 22mm.

Up to 35 hours of operations.



\$858.15

DALSF35

### Smoke Shark™ III Filter

Four layers of proven filtration gross particulate, ULPA media, activated carbon, and post filter.

Up to 6 hours of operation

Pkt/2




\$354.00

DALSF06

### 4-Stage ULPA Filtration

99.999%

efficient in trapping 0.1-0.2 µm particulates and micro-organisms






## What is surgical plume and what are the risks?

Surgical plume is a by-product from the use of electrosurgical instruments when cutting, vaporising or coagulating tissue during surgical procedures. It is generated when thermal damage of tissue causes a release of smoke. The smoke is often not visible and has an unpleasant odour.

Surgical plume can contain a mixture of potentially hazardous components including ultra-fine particulates, bacteria, viruses, cellular debris, noxious and toxic aerosols, gases, vapours and fumes. Human papillomavirus (HPV) and the human immunodeficiency virus (HIV) have also been found in the plume, although it has not been confirmed if that includes infective virus.

Whilst there is insufficient scientific research on adverse health effects from exposure to surgical plume, there is evidence that smoke from surgical plume is a potential health hazard, at least similar to smoking cigarettes, and research emphasises the basic principle that exposure to smoke of any kind should be minimised.

Reference: WorkSafe Victoria, updated on 22/06/2021, accessed 17/05/2022 <<https://www.worksafe.vic.gov.au/managing-surgical-plume-exposure-healthcare>>



### Bovie® Smart Shark Automatic Activation Device

Automatic activation device, can be used to remotely turn the evacuator on and off, saving filter life.

**\$975** DALSERS2



### Bovie® Smoke Evacuator Articulated Armstand

Enables user to place smoke evacuation tubing in place at the surgery site.


DALSEAS  
**\$527.05**



### Bovie® Shark Skin Electrosurgical Pencil Adaptor

Ergonomically designed smoke evacuation electrosurgical pencil adaptor. Transforms most electrosurgery hand-switch pencils and electrosurgical pencils into a surgical smoke evacuation device to capture smoke plume.

DALSS95 Box/20 **\$857.40**



### Bovie® Smoke Evacuation Tubing



22 x 152mm Non-Sterile  
DAL786T Box/24  
**\$448.40**

## Disposable Fine Tip Cautery

- A highly portable and easy to use cauterising tool.
- Fully disposable including the outer plastic handle, the batteries inside and the tip
- Individually sterile packed
- Includes 2 x AA batteries

Low Temp Box/10 DALAA00 **\$318.45**  
High Temp Box/10 DALAA01 **\$378.70**



## Silver Nitrate Applicators

With well-established antiseptic and antimicrobial properties, silver nitrate from an applicator or pencil is applied topically usually in aqueous solution and used as a haemostatic, for pain relief in certain conditions and for the removal of small quantities of unwanted tissue.

- Simple, Effective, Economic
- 75% Silver Nitrate
- Single use 15cm Bendable Applicator
- 100 Applicators in a resealable tube
- Removal of non-genital warts (including verruca's)
- Removal of unwanted tissue
- Cauterisation



SUTSNA7484  
Pkt/100 **\$98.10**