


Waterproof Drapes (sterile)


60 x 60cm	MGT29220	\$2.55	100+	\$2.30
90 x 90cm	MGT29223	\$4.70	100+	\$4.25
45 x 60cm with 7.5 x 6cm Fenestration	MGT29612	\$1.85	120+	\$1.65
45 x 75cm with 7.5 x 6cm Fenestration	MGT29614	\$1.95	120+	\$1.75
60 x 60cm with 7.5 x 6cm Fenestration	MGT29620	\$2.20	100+	\$1.85
90 x 90cm with 7.5 x 6cm Fenestration	MGT29634	\$3.70	100+	\$3.30
60 x 60cm with 7.5 x 6cm Fenestration & 10 x 10cm Adhesive Surround	MGT29780	\$3.45	100+	\$3.05




Plain Drape



Drape with fenestration



Drape with fenestration and adhesive surround




Plastic waterproof sheets (sterile)


60 x 45cm	MGT22200	\$0.50	100+	\$0.35
61 x 91cm	MGT22201	\$0.60	100+	\$0.50
90 x 120cm	MGT22203	\$1.00	100+	\$0.75
120 x 140cm	MGT22205	\$1.35	100+	\$1.05



Convertor Drapes



45 x 65cm
with 7cm Fenestration
CAR3523
Each \$5.95
50+ \$5.35
100+ \$4.65




45 x 65cm
CAR3520
Each \$9.15
50+ \$7.60
100+ \$6.90

Minor Procedure Surgical Drape with Circular Aperture

- Non-woven bilaminate drape with 6cm circular aperture and adhesive around the fenestration
- 61 x 66 cm

Box/25 \$62.05
2+ \$56.75

MDLDYNJSD1092



Steri-Drape™




Small Drape with Adhesive Aperture

Overall size 40 x 40cm
Adhesive area 10 x 12.5cm
Circle aperture 6.3cm

Ctn/40 \$176.65

3M1020




Drape with Small Adhesive Aperture

Overall size 57 x 56cm
Adhesive area 8.5 x 8.5cm
Circle aperture 6.3cm

Box/25 \$146.95

3M1092



Drape Sheet/Table Cover

Overall size 114 x 150cm


Box/20 \$215.30

3M9072

Looking for a different drape?
Visit www.teammed.com.au

Dressing Towels (sterile)

Paper Dressing Towel	40 x 40cm	MGT30666	\$0.55	100+	\$0.45
MediClean Low Lint Towel	35 x 30cm	MGT47680	\$0.80	400+	\$0.70
MediClean Low Lint Towel	35 x 60cm	MGT47681	\$1.25	250+	\$1.10
Huck Towel - Peel Pack	42 x 66cm	MGT47683	\$2.15	100+	\$1.85



Compro Standard Gowns



Compro gowns are made from Softesse, a fluid repellent spunlace fabric from Dupont.

- Provides uncompromising performance with a splash resistant barrier that is fluid and bacteria repellent
- Cloth-like texture that is low linting, very drapeable and designed for comfort
- Compro Standard gowns provide AAMI level 2 protection for light to moderate fluid exposure

Small	MGT28000	Ctn/20	\$215.70
Medium	MGT28001	Ctn/20	\$199.35
Large	MGT28002	Ctn/20	\$212.75
X-Large	MGT28003	Ctn/20	\$217.65




Eclipse Surgical Gowns



Light and breathable, for maximum comfort

Eclipse gowns offer outstanding softness and enhanced breathability. This light weight line has great moisture evaporation, leaving the wearer cool and dry. Those who spend much of the day in a surgical gown appreciate the comfort of Eclipse.

- Wide cut in the chest and sleeves
- AAMI Level 2
- Non-reinforced
- Sterile



Small (104cm)	MDLJP2005	Ctn/30	\$171.60	2+	\$158.85
Large (109cm)	MDLJP2001	Ctn/30	\$166.90	2+	\$152.65
X-Large (119cm)	MDLJP2002	Ctn/30	\$159.75	2+	\$146.65

Surgical Scrub suits



Single-use scrub suits contribute to the hospital's infection control. BARRIER® scrub suits are produced in a controlled environment from low-linting materials and are only worn once.

Scrub Shirt			
Small	MOL21610	Box/48	\$127.05
Medium	MOL21620	Box/48	\$131.95
Large	MOL21630	Box/48	\$136.35

Scrub Pants			
Small	MOL 21710	Box/48	\$168.65
Medium	MOL21720	Box/48	\$179.95
Large	MOL21730	Box/48	\$187.70



Disposable Economy Patient Gown

Medium (110cm long x 125cm around)	HAIDSMSEGM45	Ctn/100	\$135.70
X-Large (120cm long x 140cm around)	HAIDSMSEGL45	Ctn/100	\$129.05
XX-Large (130cm long x 155cm around)	HAIDSMSEGXXL45	Ctn/100	\$129.05



X-Ray Gown with Tie

- 76 x 106cm
- Light blue

MDLNON24354
Ctn/50 \$88.20



Aprons Plastic

Unsterile	
710mm x 1320mm	
Pkt/100	\$27.55
3+	\$26.35
MIAAP100	



MEET OUR TEAM: THERESA
CUSTOMER SERVICE - WA

Theresa is the newest member of our Perth customer service team. We are looking forward to seeing all the wonderful things Theresa will achieve in her first year working with us.

Theresa started sewing at the age of 6 and spends a lot of time with her daughter working on exciting projects. She also enjoys board games and camping.

What's your favourite takeaway food? Chinese

What series have you been watching recently? Blown away





Thumb-Hook Impervious Open Back Gown

Get hooked on convenience, comfort and cover

Reorder Code	Size
AGO30RR AGO30XL	Regular Extra Large

- Thumb-hook design for easy application
- Individually wrapped
- Impervious material
- Open-back for maximum comfort and breathability for longer wear
- Blue colour

Regulatory Compliance

TGA ARTG 266929

15 gowns per box
5 boxes per carton

Contact us to place an order



A brand by



Mun Australia
Toll free: 1800 456 837
munglobal.com.au

