

Walmay Cryoflask WM227 with 4 Tips



- New and improved model
- Strong, ergonomic design
- Increased travel in regulator handle allows for better control in the flow of liquid nitrogen.
- Fine gauze strainer on the end of siphon tube reduces contamination.
- Safety valve ensures consistent pressure in the unit.
- High centre of gravity makes for easy use and balance.
- Wide base stand provides stability when setting down device.

WALWM227

\$965



Walmay Medical
PRODUCTS



Our two Nicks are here to help you source all of your medical Nick knack needs, just in the Nick of time.



SALES - VIC

Meet Our Team: Nick L

With over 25 years of experience in the medical industry our East Melbourne Sales Rep Nick is a true asset to our team. A people person, Nick loves engaging with customers and sharing his expertise and detailed product knowledge.

For fun, he enjoys fishing, beer brewing and cooking, often choosing recipes that go well with finds from his wine cellar.

What sport do you love watching? AFL (Essendon's Mighty Bombers – established in 1872). Our VIC team is keeping up the friendly rivalry!

What is your biggest pet peeve? Teabags in the sink.

A memorable show you have been to? The Point at Portsea - like Big Day Out for old people.



SALES - QLD

Meet Our Team: Nick G



Another new addition to the Team Medical family, Nick brings seven years of dedicated experience in the medical supplies industry, most notably in the Day Hospital, GP and retail markets. His passion lies in networking and building lasting relationships, which has built him a strong network of loyal clients over the years. His ability to thrive in environments of collaboration and teamwork makes him a perfect fit for Team Medical.

When he's not hard at work, Nick carries a certain passion for all things health and fitness and he loves playing sports.

What TV shows are you streaming at the moment?

"Yellowstone" and "Jack Reacher" whenever I get a chance.



Dewar

Liquid Nitrogen Container

- Made of high-strength aluminium alloy
- Adopting advanced powder coating technology
- Each unit include leather bag and ladle
- Designed for large capacity storage and low liquid nitrogen consumption
- Safety lock lid
- 5 year vacuum warranty
- 2 year container warranty

| Model | LNC-10-50 | LNC-20-50 | LNC-30-50 |
|-------------------------------|-----------|-----------|-----------|
| Capacity | 10L | 20L | 30L |
| Weight empty (kg) | 6.0 | 9.5 | 13 |
| Caliber (mm) | 50 | 50 | 50 |
| Outer diameter (mm) | 300 | 394 | 462 |
| Height (mm) | 552 | 672 | 706 |
| Static evaporation loss (L/d) | 0.12 | 0.12 | 0.12 |
| Static holding time (days) | 86 | 168 | 254 |
| Code | PACLNC10L | PACLNC20L | PACLNC30L |



10 ltr **PACLNC10L** **\$1525**
 20 ltr **PACLNC20L** **\$1735**
 30 ltr **PACLNC30L** **\$1945**



Each dewar includes ladle

Gloves Liquid Nitrogen

- Gloves for cryogenic applications
- Length of glove is approx. 35cm
- Can bare temperatures from -180°C to -250°C
- External waterproof material, inner material has a high quality cold-proof sponge and lining



PACLNCGLV **\$52.50**

ActivArm® 23-700 Glove

- Gloves that shed water, resist salts and chemicals
- Rough-finish PVC coating provides gripping power on wet and slippery surfaces



AN23700
\$38.00

Extractor 20/30 Liter Pacific Dewar

Liquid Nitrogen extractor for filling cryospray units.

PACLNCEXT **\$215**



We recommend having a ladle as backup when purchasing an extractor.

Ladles

10 Ltr
PACLNL10L **\$140**
46cm for 10 ltr ladle

20 Ltr
PACLNL20L **\$145**
56cm for 20/30 ltr ladle



Measuring rod

The measuring rod will clearly show the LN2 level of a dewar in inches and centimetres to 92cm/36".

TAYR0338C11 **\$40**



Meet Our Team: Dee

When it comes to customer service in the medical supplies industry, it would be hard to find someone with more experience than Dee. Dee joined the team 6 years ago and has been in the industry for a total of 20 years. Her favourite part of working at Team Medical is "The people, Team Medical is a great place to work with plenty of support."

Like her favourite animal the tiger, Dee is tenacious yet calm by nature.

Pet Peeve? Toilet paper hanging wrong way (paper should face outward)

Favourite Food? Chinese




CUSTOMER SERVICE - QLD

CRYOSUCCESS Liquidfreezing

CryoSuccess with 1 and 4 mm tip

The unique method with direct contact between the liquid and the skin allowing a constant penetrating temperature of -89°C

- Medical Grade Gas
- Low treatment costs
- Simple and secure handling
- Autoclavable tips
- Extremely precise application, with a constant penetrating temperature of -89° C
- Approximately 300 seconds of treatment time per cartridge
- Cartridge can be removed and reattached without wasting gas
- Gas Cartridges have no shelf life or evaporation
- Excellent cosmetic effects; almost no scar formation
- Non bleeding treatment, no anaesthesia required
- Practically painless, therefore suitable for children and sensitive patients



Visit our website to watch the product video

Kit includes:

- 1 x CryoSuccess
- 2 x 23.5g N2O Cartridge
- 1 x 1mm Tip
- 1x 4mm Tip

\$1305

STMCSUC01

Cryosuccess 1mm Tip

- For fine work
- Autoclavable

\$360

STMSTS01



Cryosuccess 4mm Tip

- For larger areas
- Autoclavable

\$420

STMDS04



Cryosuccess Refill 23.5G Cartridge

\$112.75

STMCA2310



HistoFreezer® Portable Cryosurgical system

A professional-grade product that offers medical professionals an easy and effective method to treat nine different types of skin lesions, including:

- Acrochordon (Skin Tags)
- Actinic Keratosis (Facial & Non-Facial)
- Condyloma Acuminata (Genital Warts)
- Lentigo (Facial & Non-Facial)
- Molluscum Contagiosum (MCV)
- Seborrheic Keratosis
- Verruca Plana (Flat Warts)
- Verruca Plantaris (Plantar Warts)
- Verruca Vulgaris (Common Warts)

ECO KIT Contents:
2 x 80ml Canisters of HistoFreezer Cryogen
24 x 2mm Applicators
36 x 5mm Applicators

\$210

ECO10010412

Cryo Professional 170ml Portable Cryosurgical system

The Cryo Professional can be used to treat various types of benign skin lesions on patients over 4 years old. Examples of the types of skin lesions treatable with the Cryo Professional include:

- Genital warts
- Molluscum Contagiosum
- Fibroma, Acrochordon (skin tags)
- Verruca (Plantaris, Vulgaris, Plana)
- Seborrheic Keratosis (SKs)
- Aktinica Keratosis (AKs), both in the face and on other body parts
- Lentigos, on the face and other body parts

\$140

CRUJTM0171

Contents:
1 x 170ml aerosol dispenser mix
30 x 2mm foam sticks
30 x 5mm foam sticks

