

bodichek®
Bodichek Spots
Ø23mm

AXP13010005

Box/50 \$2.95
10+ \$2.70



Curity Spots
Ø22mm

COV6403LF

Box/100 \$5.20



Protect-A-Plast
Plastic Spots

MIAM39126

Box/100 \$6.05

MIAM39127

Box/1000 \$59.25



bodichek®
Bodichek Plastic Strips
72 x 20mm

AXP13010003

Box/50 \$4.35
10+ \$3.80

AXP13010004

Box/100 \$5.30
10+ \$4.30



Curaplast® Sensitive
Injection Dressing

20 x 40mm
Easy Dispenser box

LOH30625

Roll/250 \$20.60



Curaplast® Sensitive
Adhesive Round Dressing

Ø16mm
Easy Dispenser Box

LOH30626

Roll/100 \$9.65



Nichiban IV Pressure Pad

- Overall Product Ø36mm
- IV Pressure Pad Ø16mm
- White colour

DEF929

Box/100 \$57.00
2+ \$52.90



IV Pressure Pad

- Overall Product Ø36mm
- IV Pressure Pad Ø16mm
- Tan colour

MIDAWIVPAD

Box/100 \$44.75
2+ \$42.90



Band-Aid - Advanced Healing

Plaster seals the wound to act like an instant scab, and absorbs wound fluid to create an optimal moist healing environment for fast wound healing.

- Hydrocolloid promotes fast healing
- Acts like an instant scab
- 100% waterproof
- Hypoallergenic

Regular 25 x 66mm
Box/10 \$13.85

XXB035

Large 44 x 70mm
Box/6 \$13.85

XXB036



Bandaid Strips

JJ37206156

Box/100 \$10.60
6+ \$9.60



Elastoplast Sensitive

Dressing lengths
(6 x 10cm)

BDF46040

Box/10 \$8.45



Nexcare Tattoo Strips

"Cool" Collection

Waterproof bandage
stays on skin in
the shower or pool

3M5942

Box/20 \$8.90*



Leukoplast Elastic Spots

Latex Free



2.2cm diameter

BSN7645523


Box/100 \$9.05 12+ \$8.60

Leukoplast Universal

Plastic First Aid Dressing

Leukoplast Standard dressings are high quality, perforated plastic film strips ideal for first aid or minor wound management.

Latex-free



7.2cm x 1.9cm

BSN7646201

Box/100 \$6.65 12+ \$6.15



Leukoplast Strong Fabric Dressings

Fabric First Aid Dressing

- Ideal for use with cuts, lacerations and skin loss wounds
- Absorbent, low-adherent pad backed by a high quality stretch fabric
- Practical and highly conformable
- Strong, porous adhesive ensures a strong and lasting fixation

First aid dressing strips

2.2cm x 7.2cm

BSN7646007

Box/100 \$9.45

12+ \$8.70

BSN7646008

Box/50 \$5.10

6+ \$4.75



Finger Tip Dressing Small
45 x 50mm

BSN7646012

Box/50 \$17.80

3+ \$16.65



Finger Tip Dressing Large
50 x 63mm

BSN7646011

Box/100 \$51.05

3+ \$47.35



Knuckle Dressing
Strong Sterile 38 x 75mm Tan

BSN7646010

Box/100 \$69.40

3+ \$64.00



Leukoplast Strong Dressing Lengths

4cm x 1m	BSN7645907	Roll	\$3.75	5+	\$3.40
6cm x 1m	BSN7645905	Roll	\$4.35	5+	\$3.95
6cm x 5m	BSN7645901	Roll	\$18.20	5+	\$16.50
8cm x 5m	BSN7645902	Roll	\$23.45	5+	\$21.25

Leukoplast Strong Dressing rolls can be cut to the exact length required, making it highly cost-effective.



Leukoplast Barrier

Waterproof Barrier Dressing

2.2cm x 6.3cm BSN7645201

Box/100 \$41.55

3+ \$37.80



2.2 x 7.2cm BSN7645208

Box/100 \$43.30

3+ \$39.20



3.8cm x 3.8cm BSN7645202

Box/100 \$47.80

3+ \$44.05



5cm x 7.2cm BSN7645205

Box/100 \$98.25

3+ \$90..60



Leukoplast Barrier is a waterproof, microporous first aid dressing that breathes and is an effective bacterial and viral barrier.



Leukoplast Detectable

X-Ray and Metal Detectable, Waterproof Dressing

The waterproof film and its unique viral and bacterial barrier properties offer the best prevention of food contamination.

7.5cm x 2.2m

SN7645303 Box/100 \$79.70

7.5 x 5cm

BSN7645304 Box/100 \$132.25

5cm x 4.4cm (finger tip)

BSN7645305 Box/100 \$54.80



Cotton Tips and Applicators

7.5cm Double ended	MGT22666A	Pkt/100	\$1.25	24+	\$0.90
15cm Single end	AXP10102121	Pkt/100	\$2.45	24+	\$2.00
18cm Jumbo Swabsticks	AXP10102214	Pkt/20	\$3.00	6+	\$2.80
7.5 single end sterile - pkt/5	MGT121893	Box/100	\$24.40	3+	\$22.05
15cm single end sterile - pkt/2	MGT121891	Box/100	\$19.75	3+	\$17.85



Cotton balls

Cotton wool balls sterile	MGT02222P	Pkt/5	\$0.25	300+	\$0.20
Non-woven balls sterile	MGT00005	Pkt/5	\$0.30	500+	\$0.25
Cotton wool 0.35g/ball S&M	AXP10102311	Pkt/60	\$2.30	6+	\$1.80
Cotton wool 0.6g/ball S&M	AXP10102312	Pkt/120	\$4.35	6+	\$3.60
Cotton wool balls Swisspers	SN12260	Pkt/160	\$3.55	6+	\$3.30
Cotton wool balls large	SN2902021	Ctn/3000	\$67.20	3+	\$61.95
Cotton wool balls small	SN2903141	Ctn/4000	\$48.15	3+	\$44.40
Cotton wool roll 375g	AXP10102411	Each	\$14.05	6+	\$13.10



Eye pad

Eye Pad with Crimped Edges
(Individually wrapped in
Dispenser Box) - Sterile

MGT19111D	
Box/25	\$7.15
3+	\$6.25



Nexcare Opticlude

Nexcare Eye Patch Junior

3M1537	
Box/20	\$21.75

Nexcare Eye Patch Regular

3M1539	
Box/20	\$23.25



Combine pieces



10 x 10cm	AXP11013010	Each	\$0.50	100+	\$0.45
10 x 20cm	AXP11013012	Each	\$0.50	100+	\$0.45



10 x 12cm	MGT09888P	Each	\$0.35	100+	\$0.30
10 x 20cm	MGT09889P	Each	\$0.50	100+	\$0.35
20 x 20cm	MGT09890P	Each	\$0.60	100+	\$0.55



Combine rolls

9cm x 10m

AXP10133133	
Roll	\$12.90

20cm x 10m

AXP10133134	
Roll	\$27.00



Packing Gauze



1.25cm x 1m	MGT11111			
Each	\$0.90	100+	\$0.65	
2.5cm x 1m	MGT11112			
Each	\$0.75	100+	\$0.60	




X-Ray Gauze

10 x 10cm	MGT15666B	Pkt/5	\$0.80	
5cm x 1m	MGT15112	Each	\$1.00	



Green Gauze

10 x 10cm (non-sterile)				
MGT18888B	Pkt/100	\$5.25		
10 x 10cm (sterile)				
MGT17888	Pkt/5	\$0.90		



Ear Wick

MEROCEL Ear Wicks, Pope Oto Wick

- Ideal for treating otitis externa
- Expands to conform to ear canal when hydrated with otic solutions and local anesthesia

XOM1425008	15mm
Pkt/5	\$41.15





Gauze Swabs
non-sterile (8ply)

This soft and highly absorbent gauze can be used to clean a wound before it is dressed.



5cm	BP010805BA		
Pkt/100	\$2.15	200+	\$1.95
7.5cm	BP010875BA		
Pkt/100	\$3.25	100+	\$2.95
10cm	BP010810BA		
Pkt/100	\$4.15	50+	\$3.75



Gauze non-sterile
Low-Lint

7.5cm	AXP10101562		
Pkt/100	\$4.10		
100+	\$3.50		



Gauze non-sterile

7.5cm	SN36361349		
Pkt/100	\$3.80		
100+	\$2.95		



Sterile Gauze

7.5 x 7.5cm	MGT11600A	Pkt/3	\$0.30	500+	\$0.20
7.5 x 7.5cm	MGT11666	Pkt/5	\$0.50	500+	\$0.35
10 x 10cm	MGT11559B	Pkt/3	\$0.45	500+	\$0.40
10 x 10cm	MGT11668	Pkt/5	\$0.45	500+	\$0.40



Sterile Gauze in dispenser pack

7.5cm x 7.5cm (pkt/5) dispenser pack		
Box/50	\$18.10	MGT11666D
6+	\$15.45	
10cm x 10cm (pkt/3) dispenser pack		
Pkt/50	\$27.60	MGT11559D
8+	\$22.45	



Non-woven Gauze sterile

7.5 x 7.5cm	MGT47600A	Pkt/3	\$0.25	500+	\$0.20
7.5 x 7.5cm	MGT47601P	Pkt/5	\$0.35	500+	\$0.30
7.5 x 7.5cm (pkt/3) dispenser pack					
Box/30	\$7.80	MGT47600C			
6+	\$6.65				



Sterile Gauze Handy

5 x 5cm	SN7111705		
Pkt/100	\$5.45		
10+	\$5.10		
7.5 x 7.5cm	SN5673		
Pkt/100	\$9.40		
10+	\$8.65		
10 x 10cm	SN7111706		
Pkt/100	\$16.05		
10+	\$14.80		



Topper 8 Gauze

Topper 8 is composed of nonwoven viscose rayon and polyester. These cost effective low-linting swabs with high absorbency are ideal for use in all general wiping and cleansing procedures associated with wound management.


- Ribbed structure provides a textured surface to facilitate debris removal
- The non-woven fabric permits high absorbency without linting, unlike gauze swabs

Non Sterile

5 x 5cm	JJM50405	Pkt/100	\$5.50
7.5 x 7.5cm	JJM53407N	Pkt/200	\$17.40

Sterile

5 x 5cm (5's)	JJTS8055	Pkt/25	\$19.00
7.5 x 7.5cm (2's)	JJTS8072	Box/50	\$22.30
10 x 10cm	JJTS8102	Box/75	\$49.75





Interpose®-Lite

Non-adherent dressing, light absorption

7.5 x 5cm	MGT30501	Box/100	\$16.10	3+	\$14.20
7.5 x 10cm	MGT30503	Box/100	\$23.95	3+	\$21.35
7.5 x 20cm	MGT30504	Box/100	\$42.40	3+	\$31.05



Interpose®

Non-adherent dressing, high absorption

5 x 5cm	MGT30381	Box/100	\$20.20	3+	\$17.10
10 x 10cm	MGT30382	Box/100	\$51.80	3+	\$46.70
10 x 20cm	MGT30383	Box/100	\$86.00	3+	\$76.50



Bodichek Non Adherent Dressing Light

- All in one dressing suitable for lightly draining wounds
- Sterile and ideal for first and second degree wounds
- Highly absorbent to keep wound dry
- Ideal for minor cuts, grazes and abrasions
- Non-stick for pain free dressing removal
- Hypoallergenic and gentle so suitable for sensitive skin



5 x 7.5cm	AXP13070026	Box/50	\$8.20	3+	\$7.45
7.5 x 10cm	AXP13070027	Box/50	\$15.45	3+	\$13.85
10 x 10cm	AXP13070028	Box/50	\$17.40	3+	\$15.90
10 x 20cm	AXP13070029	Box/25	\$16.35	3+	\$14.70

Cutilin

Cutilin is a low-adherent dressing for exuding wounds that can be cut to suit, without linting.

- Step 1: Clean the wound with saline solution. Do not dry the wound bed.
- Step 2: Cutilin can be secured with an adhesive tape or bandage.
- Step 3: Cutilin can be left in place until the capacity of the pad is exhausted, which can be seen when exudate reaches the edges of the dressing

5 x 5cm	SN76240	Box/100	\$36.05	3+	\$33.25
10 x 10cm	SN76242	Box/100	\$88.75	3+	\$81.90




SmithNephew

Melolite

Lightly absorbent low-adherent dressing

- Double sided low-adherent wound dressing
- Absorbent pad
- Highly conformable to all wound sites
- Either side of the dressing can be placed directly onto the wound site

7.5 x 5cm	SN4811	Box/100	\$23.15	3+	\$21.35
7.5 x 10cm	SN4812	Box/100	\$36.60	3+	\$33.75
7.5 x 20cm	SN4814	Box/100	\$59.55	3+	\$55.60



SmithNephew

Melolin

Highly absorbent low-adherent dressing

Melolin consists of a highly absorbent cotton and acrylic fibre pad which is heat bonded on one side to a very thin perforated polyester film. The film side of the dressing is placed next to the wound.

- Rapid drainage of exudate, reducing trauma to healing tissue
- Comfortable and minimises pain on removal
- Easy to cut to shape

5 x 5cm	SN66974940	Box/100	\$38.90	3+	\$35.90
10 x 10cm	SN66974941	Box/100	\$96.20	3+	\$88.75
10 x 20cm	SN66974939	Box/100	\$170.40	3+	\$157.20



Smith+Nephew

Zetuvit

Non-adherent absorbent dressing

Zetuvit is a sterile single use absorbent dressing used for the treatment of exuding wounds, acute or chronic, and for upholstery and cushioning of wounds. It's high absorbency reduces the risk of maceration to the wound and surrounding tissue.

10 x 10cm	HTM413701	Box/25	\$58.45	3+	\$53.40
10 x 20cm	HTM413702	Box/25	\$66.00	3+	\$60.30
20 x 20cm	HTM413703	Box/15	\$61.60	3+	\$56.25



Zetuvit® Plus

Super absorbent dressing

Zetuvit Plus is a highly absorbent dressing suitable for the management of heavily exuding wounds. Zetuvit Plus has outstanding absorbency capacity due to super absorbent polymers within the dressing. This cotton based dressing is soft and conformable for maximum patient comfort. The special backing material prevents strike-through of exudate which promotes extended wear time.

- Outstanding absorbency
- Superior patient comfort
- Excellent cost-effectiveness

10 x 10cm	HTM413710	Box/10	\$34.90	3+	\$31.85
10 x 20cm	HTM4137113	Box/10	\$43.55	3+	\$39.80
15 x 20cm	HTM413712	Box/10	\$57.45	3+	\$52.50
20 x 25cm	HTM4137133	Box/10	\$78.40	3+	\$71.65
20 x 40cm	HTM4137150	Box/10	\$113.20	3+	\$103.40



Zetuvit®

Plus Silicone Border

Zetuvit Plus Silicone Border is a super absorbent dressing with a silicone adhesive contact layer and border. The absorbent core draws fluid away from the wound and traps it in the dressing via a unique combination of loose cellulose and super absorbent polymers (SAPs). The silicone interface and border enable atraumatic removal and repeated lifting/ repositioning. The backing film is shower-proof and breathable.

- Effective
- Versatile
- Patient friendly

10 x 10cm	HTM4139100	Box/10	\$115.40	2+	\$109.70
12.5 x 12.5cm	HTM4139200	Box/10	\$125.05	2+	\$118.85
17.5 x 17.5cm	HTM4139300	Box/10	\$183.45	2+	\$172.00
15 x 25cm	HTM4139400	Box/10	\$352.80	2+	\$330.75
20 x 25cm	HTM4139500	Box/10	\$423.35	2+	\$396.90



Telfa® Ouchless™ Non-Adherent Dressings

Superior “Ouchless” dressing won’t disrupt healing tissue by sticking to wound. Ideal as a primary dressing for lightly draining wounds.

5 x 7.5cm	COV1961	Box/100	\$22.65	3+	\$19.80
7.5 x 10cm	COV2132	Box/100	\$27.80	3+	\$24.25

Telfa® Ouchless™ Adhesive Dressings

Highly absorbent cotton fabric bonded on both sides with perforated non-adherent film. Won’t disrupt healing tissue by sticking to wound.

5 x 7.5cm	COV6017	Box/100	\$22.35	3+	\$18.80
7.5 x 10cm	COV7643	Box/100	\$32.00	3+	\$27.95



Mesorb®



Highly Absorbent dressing with a fluid-repellent backing

Mesorb is a soft highly absorbent dressing that can retain large amounts of wound exudate. The fluid-repellent strikethrough barrier helps prevent exudate from soiling clothing and bed linen.

- High absorption and retention capacity
- Conformable and comfortable for patients

10 x 13cm	MOL677001	Box/50	\$68.90
15 x 23cm	MOL677601	Box/50	\$121.95



Maximum exudate level

Jelonet

Paraffin Gauze Dressing

JELONET is a soft paraffin dressing made from open-weave gauze. The gauze has interlocking threads, which minimise fraying when the dressing is cut to shape.

JELONET is soothing and low-adherent, and allows the wound to drain freely into an absorbent secondary dressing.

Indications: Burns, scalds, lacerations, abrasions and leg ulcers.

5cm x 5cm	SN7403	Box/50	\$79.95	3+	\$73.75
10cm x 10cm	SN66007450	Box/50	\$82.15	3+	\$75.80
10cm x 10cm	SN7404	Box/10	\$16.30	3+	\$15.05



Smith+Nephew



Bactigras


Medicated Tulle Gras

BACTIGRAS is a gauze of leno weave impregnated with White Soft Paraffin BP containing 0.5% Chlorhexidine Acetate BP in a sterile presentation.

BACTIGRAS is effective against a wide range of gram-positive and gram-negative wound pathogens, including MRSA.

BACTIGRAS is for use as a topical treatment only. It is indicated for a wide range of wounds where there is a risk of infection. For already infected wounds, BACTIGRAS may be used in conjunction with systemic antibiotics.

The BACTIGRAS tulle gras dressing, made from an open-weave gauze, has low adherence and allows the wound to drain freely into an absorbent secondary dressing.



Indications: Suitable for covering wounds such as minor burns, lacerations, abrasions, graft sites and leg ulcers.

5 x 5cm	SN7456	Box/50	\$85.75	3+	\$79.10
10 x 10cm	SN7457	Box/10	\$21.25	3+	\$19.60
10 x 10cm	SN66003650	Box/50	\$104.90	3+	\$96.75

Smith+Nephew



Cuticerin

Low-Adherent, Impregnated Gauze Dressing

An impregnated, non-aqueous, non-linting dressing ideal as a wound contact layer for burns, abrasions and other wounds with superficial skin loss. It reduces the risk of granulation tissue growing through the dressing, decreasing trauma on removal.

Indications: For treatment of extensive, exuding wounds – burns, abrasions, radiation injuries, skin graft and donor sites.



Smith+Nephew

7.5cm x 7.5cm	SN66045560	Box/50	\$76.45
		3+	\$70.50

Adaptic
Non-Adherent Dressing

Indicated for dry to highly exuding wounds where adherence of dressing and exudate is to be prevented.
A primary dressing made of knitted cellulose acetate fabric and impregnated with a specially formulated petrolatum emulsion.

When to use it:
1st and 2nd degree burns, Lacerations and abrasions, Grafts, Venous ulcers, Pressure ulcers, Nail extractions, Eczema, Staples and suture lines, Surgical incisions and reconstructive procedures

7.6 x 7.6cm	JJ2012	Box/50	\$68.95	3+	\$63.85
7.6 x 20cm	JJ2015	Box/24	\$78.45	3+	\$73.45

PROTECTS THE WOUND
MINIMISES PAIN
AND MACERATION



Adaptic Touch
Non-Adhering Silicone Dressing

A flexible, open-mesh silicone wound contact layer

- Designed to facilitate fluid transfer to secondary dressing and minimise both adherence and pain at dressing change
- Protects the wound and works as a primary dressing
- Helps prevent adherence of the secondary dressing to the wound
- Minimises patient pain and trauma at dressing change
- Allows free passage of exudate, reduces the risk of exudate pooling and maceration

5 x 7.6cm	KCITCH501	Each	\$8.10	3+	\$7.90
7.6 x 11cm	KCITCH502	Box/10	\$143.35	3+	\$137.50

ALL YOU COULD WISH FOR
IN A SILICONE CONTACT LAYER



Adaptic Digit
Non-Adhering Digit Dressing

Made of a unique, tacky, silicone-coated fabric covered with a tubular bandage. It is specifically designed for dressing fingers and toes.

ADAPTIC™ DIGIT does not adhere to the wound. It is conformable and non-bulky, therefore allowing maximum digit flexibility. It does not adhere to the wound. Easy to apply and remove.

Indicated to manage a range of mild to severe digit injuries, such as:
Lacerations and abrasions, traumatic injuries, nail extractions, suture lines

Small	Ø2cm	KCIMAD003	Box/10	\$94.65
Medium	Ø2.4cm	JJMAD013	Box/10	\$103.80
Large	Ø2.8cm	JJMAD023	Box/10	\$106.90
Toe	Ø2.8cm	KCIMAD062	Box/10	\$95.90


COMFORTABLE
EASY TO APPLY
AND REMOVE



GOING ABOVE AND BEYOND

Our entire team is committed to providing a positive customer experience, every time you order. Everything from the way you place your order to the way our warehouse team packs your order is optimised to provide the ultimate customer experience.

By listening to our customers' feedback and empowering our employees, we continually improve the customer experience through innovation.






ACTIVHEAL® ACADEMY

FREE ONLINE WOUND CARE EDUCATION

The ActivHeal® Academy is a multi-tiered wound care education programme designed to support clinicians of all levels with easily accessible, online education resources.

What’s great about the ActivHeal® Academy? You learn at your pace. Access and complete courses at a time that suits you, no start date or deadlines.



ACTIVHEAL ACADEMY SIMPLIFIED MODULES

Access to the simplified modules is free. Each module takes 30 minutes to complete.

Why Wounds Fail to Heal Simplified

Dressing Selection Simplified

Wound Assessment Simplified

Moisture Lesions Simplified

Pressure Ulcers Simplified

Pressure Ulcer Prevention Simplified

Skin Tears Simplified

Infection Simplified

Biofilm, Bioburden and Antimicrobials Simplified

ACTIVHEAL ACADEMY ADVANCED MODULES

Access to the advanced modules is available exclusively for ActivHeal Customers. Each module takes 60-90 minutes to complete.

Anatomy and Physiology of the Skin

Wound Healing

Patient and Wound Assessment

Wound Classification and Wound Management


Dressing Selection

- 
SIMPLIFIED GUIDES
- 
COMPREHENSIVE MODULES
- 
ONLINE TESTS
- 
PDF CERTIFICATES
- 
REFLECTIVE LEARNING

Sign up today at

academy.activheal.com





ActivHeal® SILICONE WOUND CONTACT LAYER

MINIMISES
EPIDERMAL
STRIPPING

REDUCES PAIN
AT DRESSING
CHANGES


PROTECTS THE
PERIWOUND

14 DAYS
IN SITU

ActivHeal® Silicone Wound Contact Layer
prevents the dressing from sticking to the wound
and reduces pain during dressing changes.



- ✓ Highly conformable
- ✓ Allows passage of exudate
- ✓ Gentle yet secure adhesion
- ✓ Can remain in situ for 14 days
- ✓ Minimises epidermal stripping¹
- ✓ Minimises pain at dressing changes¹
- ✓ Protects the periwound
- ✓ Maintains a moist wound environment



BEST
BUY

Indications for use:

ActivHeal® Silicone Wound Contact Layer dressing is indicated for use on nil to heavily exuding chronic and acute wounds (with appropriate secondary dressing). The dressing may be used throughout the healing process on the following wounds:

- ✓ Surgical Incisions
- ✓ Skin Tears
- ✓ Pressure Ulcers
- ✓ Venous Leg Ulcers
- ✓ Arterial Leg Ulcers
- ✓ Partial Thickness Burns
- ✓ Partial Thickness Burns Skin Grafts
- ✓ Can be used in conjunction with Negative Pressure Wound Therapy (NPWT)

Nil > Low > Moderate > Heavy Exudate



Reference: 1. Barrett, S.(2012) BJN 21 (21):1271-7.
ActivHeal is a registered trademark of Advanced Medical Solutions Limited.

5 x 7.5cm	ADM9001241	Box/10	\$60.10*	3+	\$55.40*
10 x 10cm	ADM9001258	Box/10	\$94.55*	3+	\$87.25*
10 x 20cm	ADM9001265	Box/10	\$170.20*	3+	\$154.20*
15 x 15cm	ADM9001272	Box/10	\$139.30*	3+	\$126.20*

Mepitel®




Wound contact layer dressing with Safetac® on both sides



Mepitel® is a gentle, effective wound contact layer using Safetac®. The open mesh design allows exudate to pass into an outer absorbent dressing. The Safetac layer prevents the outer dressing from sticking to the wound and ensures atraumatic dressing changes. The Safetac layer seals around the wound edges, preventing exudate leaking onto the surrounding skin, thus minimising the risk of maceration.

- Minimises pain and trauma at dressing changes¹
- Can remain in place for up to 14 days
- Maintained product properties over time - leaves no residues and does not dry out
- Transparent for easy wound inspection during application and during wear
- Low potential for skin irritation and allergy

5 x 7.5cm	MOL290510	Box/10	\$113.75*	3+	\$103.50*
7.5 x 10cm	MOL290710	Box/10	\$186.75*	3+	\$169.90*



Mepitel® One




Wound contact layer dressing with Safetac® on one side



An innovative one-sided Safetac® wound contact layer, Mepitel One lets you give the healing process the time it needs as it can be left in place for up to 14 days. Avoiding the need for dressing removal helps minimise patient pain and unnecessary disturbance to the healing process.

- Minimises trauma to the wound and pain at dressing changes¹
- Adheres gently to the surrounding skin on removal and minimises the risk of maceration
- Mepitel One supports optimal undisturbed healing
- Built-in directional arrows, perfect to guide nurses on the direction it should be removed

6 x 7cm	MOL289170	Box/5	\$65.60*	3+	\$59.65*
9 x 10cm	MOL289270	Box/5	\$113.75*	3+	\$103.45*



What is the technical difference between Mepitel One and Mepitel?
Mepitel One consists of a polyurethane film coated on one side with Safetac technology. Mepitel is a polyamide-net, coated on both sides with Safetac technology.

Atrauman®


Atraumatic wound contact layer

Atrauman is a triglyceride impregnated wound contact layer. It provides a soothing protective layer to wounds such as skin tears, abrasions, grafts, and minor burns. The small mesh weave reduces the risk of tissue growing into the dressing but allows the passage of exudate to pass freely. Atrauman leaves less residue in the wound than traditional paraffin gauze dressings.

5 x 5cm	HTM449550	Box/50	\$59.00
7.5 x 10cm	HTM499553	Box/50	\$82.15
10 x 20cm	HTM499536	Box/30	\$88.10



APPLICATION:



Atrauman® Silicone


Silicone wound contact layer

Atrauman Silicone is a soft and conformable wound contact layer ideal for wounds such as skin tears and other superficial wounds. It allows pain free dressing changes for those patients with fragile skin. Atrauman Silicone can stay in place for up to 7 days and often only the secondary dressing needs to be changed. The small weave structure prevents the tissue growing into the dressing. Ideal for challenging to dress wounds, such as fingers and toes.

5 x 7cm	HTM499568	Box/10	\$75.25
7.5 x 10cm	HTM499561	Box/5	\$64.80



APPLICATION:





Non-Woven Island Dressing

- Flexible and soft fabric dressing
- Easily conforms to the body maintaining high comfortability
- Secure, light stretch fabric with superior elasticity
- Hypoallergenic and latex-free acrylic adhesive
- Absorbent non-adherent pad absorbs more than 8 times its weight

5 x 4cm	AXP13070001	Box/50	\$19.65	3+	\$18.30
5 x 7cm	AXP13070002	Box/50	\$22.25	3+	\$20.75
7 x 9.5cm	AXP13070003	Box/25	\$15.30	3+	\$14.25
9 x 12cm	AXP13070004	Box/25	\$22.40	3+	\$20.85
9 x 15cm	AXP13070005	Box/20	\$22.95	3+	\$21.40
9 x 20cm	AXP13070006	Box/20	\$27.85	3+	\$25.90



BEST
SELLER

Primapore

Adhesive Wound Dressing

Primapore combines an absorbent pad with a soft and conformable fixative layer

- Soft breathable cover
- Highly absorbent (melolin) pad
- Low allergy adhesive

7.2 x 5cm	SN66003634	Box/50	\$23.85	3+	\$22.05
8.3 x 6cm	SN7135	Box/50	\$36.75	3+	\$33.90
10 x 8cm	SN66000317	Box/20	\$39.10	3+	\$36.05
15 x 8cm	SN66000318	Box/20	\$55.85	3+	\$51.55
10 x 25cm	SN66000320	Box/20	\$70.75	3+	\$65.30



Smith+Nephew



Cutiplast Steril

Conformable Fabric Island Dressing

Cutiplast Steril is a dressing for primary wound management. The wound pad is low-adherent so dressing changes do not disturb or interfere with the wound healing process.

7.2 x 5cm	SN66001478	Box/100	\$33.95	2+	\$31.35
10 x 8cm	SN66001473	Box/50	\$57.10	2+	\$52.65
15 x 8cm	SN66001474	Box/50	\$81.60	2+	\$75.30



Smith+Nephew

Cutiplast

Conformable Fabric Dressing Strips

Cutiplast conformable fabric island dressings are ideal for the treatment of minor wounds. Economical; dressing lengths can be cut to size.

4cm x 5m	SN66001464	Roll	\$15.65	2+	\$14.45
6cm x 5m	SN66001466	Roll	\$19.30	2+	\$17.80
8cm x 5m	SN66001468	Roll	\$23.00	2+	\$21.20



Smith+Nephew

Compose® Island Dressing

Adhesive Island Dressing

7.5 x 5cm	MGT30611	Box/50	\$9.45	3+	\$9.05
8.25 x 6cm	MGT30612	Box/50	\$16.95	3+	\$14.45
12 x 8.25cm	MGT30613	Box/20	\$11.70	3+	\$10.20
8.25 x 18cm	MGT30614	Box/20	\$15.85	3+	\$14.80



MULTIGATE

Bodichek® Waterproof Island Film Dressing

- Highly absorbent non-adherent pad, won't stick to wounds
- Waterproof but breathable
- Hypo-allergenic adhesive

5 x 4cm	AXP13070009	Box/50	\$45.30	3+	\$42.15
5 x 7cm	AXP13070010	Box/50	\$60.40	3+	\$56.25
7 x 9.5cm	AXP13070011	Box/25	\$46.85	3+	\$41.85
9 x 12cm	AXP13070012	Box/25	\$56.30	3+	\$50.30
9 x 15cm	AXP13070013	Box/20	\$55.50	3+	\$49.65
9 x 20cm	AXP13070014	Box/20	\$67.60	3+	\$60.45



Claripose® Island Dressing

Waterproof Island Dressing

6.5 x 5cm	MGT30651	Box/100	\$66.10*	3+	\$61.40*
9.5 x 8.5cm	MGT30652	Box/20	\$29.65*	3+	\$27.55*
10 x 12cm	MGT30653	Box/10	\$19.15*	3+	\$16.90*
15 x 8.5cm	MGT30654	Box/20	\$38.10*	3+	\$34.55*



Leukomed® Control

Leukomed Control combines a highly innovative hydropolymer technology with the protective and fully transparent benefits of the Leukomed film. This advanced Post-Op dressing makes early detection of wound infection possible without lifting or changing it. Leukomed Control provides ideal healing conditions by creating an optimal moisture level and adapting to skin's elasticity. The dressing puts medical caregivers in control of the healing progress and allows for quick mobilisation of the patient.

Leukomed Control comes with:

- Transparent, visual control to assess therapeutic wound progress
- Unique application via simple step-by-step technique allows for precise placement and quick wound protection
- Thin and absorbent hydropolymer pad keeps the wound clean and ready for epithelialisation
- Ultra-thin, highly breathable and conformable film for a second skin feeling



5 x 7cm	BSN7323000	Box/10	\$54.90*	3+	\$53.10*
7 x 10cm	BSN7323001	Box/10	\$57.20*	3+	\$55.30*
8 x 15cm	BSN7323002	Box/10	\$78.25*	3+	\$75.10*



Looking for more information on a product? Visit our website or call our team. We're always happy to help.



3M™ Tegaderm™ + Pad Film Dressing with Non-Adherent Pad



A waterproof, absorbent dressing.

- All-in-one sterile dressing – reduces application time
- A bacterial and viral barrier – helps prevent external bacteria and contaminants from entering the site
- A waterproof barrier – patient may shower without removing dressing
- Conforms to body contours, providing comfort
- Non-adherent absorbent pad
- Breathable film allows skin to function normally with good exchange of moisture vapour and oxygen
- Provides moist wound environment to enhance healing and lessen pain
- Frame delivery system ensures precise placement without applying tension
- Hypoallergenic

Suggested Applications:

- Post-operative sites
- Superficial and partial thickness wounds
- Lacerations, abrasions and burns
- Light to moderately draining wounds
- Secondary wound dressing

5 x 7cm	3M3582	Box/50	\$54.85*	2+	\$50.95*
6 x 10cm	3M3584	Box/50	\$105.65*	2+	\$98.15*
9 x 10.5cm Oval	3M3587	Box/25	\$63.25*	2+	\$58.75*
9 x 10cm	3M3586	Box/25	\$59.75*	2+	\$55.50*
9 x 15cm	3M3589	Box/25	\$74.45*	2+	\$69.15*
9 x 25cm	3M3591	Box/25	\$99.90*	2+	\$92.15*



Picture frame delivery for easy application

3M™ Tegaderm™ Transparent Film Dressing



- Framed delivery system allows for one handed application and prevents dressing sticking to itself
- Transparent dressings allow easy monitoring of wounds and I.V. sites without the need to remove the dressing
- Secure adhesion, yet gentle to the skin

4.4 x 4.4cm	3M1622W	Box/100	\$56.20*	2+	\$51.85*
6 x 7cm	3M1624W	Box/100	\$57.60*	2+	\$53.50*
10 x 12cm	3M1626W	Box/50	\$68.95*	2+	\$64.05*
15 x 20cm	3M1628	Box/10	\$51.05*	2+	\$48.35*



3M™ Tegaderm™ IV Transparent Film Dressing with Border



- Notched design provide a better seal around the catheter
- Sterile tape strips are provided to assist with improved catheter securement
- Breathable film, allows oxygen in and moisture out
- Transparent dressings allow easy monitoring of I.V. sites without the need to remove the dressing
- Framed delivery system allows for one handed application and prevents dressing sticking to itself

7 x 8.5cm	3M1633	Box/100	\$114.95*	2+	\$106.80*
-----------	--------	---------	-----------	----	-----------



3M™ Tegaderm™ Absorbent Clear Acrylic Dressing



Transparent plus absorbent for optimal clinical performance and wear time.

A transparent wound dressing with fluid handling capacity that allows for wound monitoring without changing the dressing. For low to moderate exudate. It is waterproof, washable and breathable.

- Easy monitoring of the wound, exudate and periwound skin
- No upper wear time - eliminates unnecessary dressing changes
- Eliminates hydrocolloid-related odour
- No dressing breakdown in the wound

7.6 x 9.5cm	3M90800	Box/5	\$49.65*	3+	\$46.65*
11.1 x 12cm	3M90801	Box/5	\$73.50*	3+	\$69.60*
14.9 x 15.2cm	3M90802	Box/5	\$124.90*	3+	\$117.35*



Suggested Applications:

- Partial and full thickness dermal ulcers
- Pressure ulcers
- Abrasions, superficial wounds, donor sites and skin tears

+ When you value quality and patient satisfaction

The unique design of OPSITE POST-OP VISIBLE allows continual monitoring of the incision site without the need to disrupt the healing process 'providing both visibility and absorbency to promote longer wear times and fewer unnecessary dressing changes.'²

Getting you closer to zero delays in wound healing.

OPSITE POST-OP VISIBLE

When you need continual monitoring.

- REACTIC[®] film
- High MVTR film
- Lattice hydrocellular foam
- Bacterial barrier
- Waterproof film

OPSITE POST-OP

Offers a significantly longer wear time when compared with non-woven or passive dressings.¹

- REACTIC[®] film
- High MVTR film
- Low-allergy adhesive
- Bacterial barrier
- Waterproof film



Smith+Nephew

OPSITE[®] POST-OP

Absorbent, Waterproof
and Bacteria-proof Film Dressing

OPSITE[®] POST-OP VISIBLE

Waterproof, Bacteria-proof Dressing
with See-through Absorbent Pad



1. Barker A, et al. A ward-based in-practice evaluation of a post-operative dressing. Poster presented at EWMA conference, Geneva; 2010. 2. O'Brien G, et al. A prospective evaluation to assess performance of a new post-operative dressing, Poster presented at Infection Prevention Society Conference, Harrogate; 2009. Smith & Nephew Pty Ltd (Australia) www.smith-nephew.com/australia. Smith & Nephew Ltd (New Zealand) www.smith-nephew.com/new-zealand. [®]Trademark of Smith+Nephew. All trademarks acknowledged. SN13505 V1 03/20. For detailed product information, including indications for use, contraindications, precautions and warnings, please consult the product's applicable Instructions for Use (IFU) prior to use.

CUTIFILM Plus

CUTIFILM Plus is a waterproof dressing ideal for use on surgical wounds, lacerations and abrasions.

5 x 7.2cm	SN76248	Box/50	\$78.00*	3+	\$72.00*
8 x 10cm	SN76244	Box/25	\$84.00*	3+	\$77.50*
8 x 15cm	SN76245	Box/25	\$111.45*	3+	\$102.85*
10 x 12.5cm	SN76246	Box/25	\$115.40*	3+	\$106.45*
10 x 20cm	SN76247	Box/25	\$150.90*	3+	\$139.20*



- Sterile, bacteria proof and waterproof
- Hypo-allergenic adhesive
- Secure adhesion
- Low adherent pad
- Protects the skin
- Allows the skin to breathe

OPSITE Post-Op

Post-operative dressing combines the advantages of OPSITE film with an absorbent dressing pad to provide a superior dressing for many types of wound.

- Waterproof, breathable and conformable
- Bacterial barrier helps to reduce the risk of secondary infection
- Absorbent, non-adherent pad minimises the number of dressing changes

Indications: Post-operative wounds, skin tears, minor cuts, abrasions and lacerations

6.5 x 5cm	SN66000708	Box/100	\$176.35*	3+	\$162.70*
9.5 x 8.5cm	SN66000709	Box/20	\$81.15*	3+	\$74.85*
12 x 10cm	SN66000710	Box/10	\$57.05*	3+	\$52.65*
15.5 x 8.5cm	SN66000712	Box/20	\$107.70*	3+	\$99.35*

OpSite® Post-Op is a revolutionary post-operative dressing combining a high Moisture Vapour Transmission Rate (MVTR) film and the cushioning and absorbency of a low adherent pad together with an easy application carrier system. The hydrophilic polyurethane film on OpSite Post-Op has REACTIC® properties that 'switch on' in the presence of moisture allowing any excess to transpire into the atmosphere.

Smith+Nephew



BEST
SELLER

★★★★★ PREMIUM PRODUCT

OPSITE Post-Op Visible


Waterproof, bacteria-proof dressing with see-through absorbent pad

OpSite Post-Op Visible is a transparent film with a unique see-through absorbent pad. The highly breathable film is waterproof and bacteria-resistant, and the unique lattice structure foam pad enables regular wound assessment without the need to lift or remove the dressing.

- Absorbent lattice structured foam pad
- Highly breathable waterproof outer film
- Perforated low-adherent wound contact layer
- Low allergy adhesive
- Easy application flexible carrier

Indications: ideal for moderately exuding wounds including surgical wounds, lacerations and abrasions. The low allergy adhesive reduces the risk of irritation and keeps the dressing in place, even over awkward contours.

8 x 10cm	SN66800136	Box/20	\$110.65*	3+	\$102.05*
15 x 10cm	SN66800137	Box/20	\$162.45*	3+	\$149.90*
20 x 10cm	SN66800138	Box/20	\$198.40*	3+	\$183.05*



Smith+Nephew

For complete listings of Smith & Nephew products, treatment options, education and event calendar visit <https://www.smith-nephew.com/australia/>



Smith+Nephew

Claripose® Transparent Dressing

Transparent Dressing

6 x 7cm	MGT30632			
Box/100	\$44.35*	3+	\$40.20*	
10 x 12cm	MGT30633			
Box/50	\$56.55*	3+	\$50.15*	



Bodichek® Waterproof Film Dressing

- Ideal for minor cuts, grazes, abrasions, cover for wounds and IV / Catheters
- Waterproof protective film that is also breathable
- Hypo-allergenic adhesive
- Release paper frame allows easy and effective fixation without adherence to itself

5 x 4cm	AXP13070017	Box/50	\$20.50	3+	\$18.35
5 x 7cm	AXP13070018	Box/50	\$34.70	3+	\$29.60
9.5 x 7cm	AXP13070019	Box/25	\$33.15	3+	\$28.00
12 x 9cm	AXP13070020	Box/25	\$34.20	3+	\$31.40



Mepitel® Film

Breathable, transparent film dressing with Safetac® Technology

Mepitel® film provides a flexible, transparent covering to protect the skin or wound from microbial contamination and fluid strikethrough. It maintains a moist environment even though the vapour permeability allows excess moisture to pass away from the skin. It provides instant tack adhesion that minimises the requirement of extra fixation.

Benefits:

- Minimises pain and trauma at dressing changes¹
- Limits the risk for damage of newly formed tissue
- Transparent for easy inspection
- Conforms well to body contours, promoting patient comfort during wear
- Bacteria and viral barrier where microbes larger than 25nm
- Low potential for skin irritation and allergy
- Leaves no residues
- Fluid impermeable
- Showerproof



6 x 7cm	MOL296170	Box/10	\$23.55*	2+	\$21.30*
10 x 12cm	MOL296270	Box/10	\$36.90*	2+	\$33.65*

MEET OUR TEAM: PIETER MANAGEMENT - SYDNEY

Our CEO Pieter decided the corporate world wasn't for him and joined his sister Marijke and son Matt at Team Medical in 2012. His son Tom soon followed in 2013 truly making Team Medical a family company.

With extensive industry experience, Pieter has worked tirelessly towards ensuring continued growth by being the best at what we do – whilst having fun along the way.

Pieter consistently embeds his family values into our company culture. "At Team Medical, everyone is valued and respected. Seeing our people grow and thrive alongside our company and customers is very rewarding."

Born in the Netherlands, Pieter has recently been back there for a family wedding. He enjoys gardening, playing tennis and watching soccer – no game is on too early or too late for Pieter, and always better watched with a crowd.

What is your biggest pet peeve? Cleaning dog poo, I have 3 dogs.

What app do you use every day? Fitbit

What show have you been watching lately? Breaking Bad



OPSITE Flexifix

Transparent, Non-Sterile, Adhesive Film Roll

- Waterproof and bacteria proof
- Breathable
- Highly conformable: easy application to awkward areas of the body
- Ideal over primary dressings, fixation of tubes

5cm x 10m	SN66000040	
Roll	\$59.95*	3+ \$55.35*
10cm x 10m	SN66000041	
Roll	\$107.80*	3+ \$99.45*



The unique and novel roll format enables rapid and economical use of the OpSite film to create comfortable, waterproof fixation for dressings.

SmithNephew

OPSITE spray

Spray Film Dressing

- Transparent and quick drying film
- Waterproof
- Easy and convenient to apply
- Reduces the risk of secondary infection

Indications: The film provides protection for dry minor surgical and surface wounds.

SN66004978	
110ml	\$39.15*
3+	\$36.15*




SmithNephew

OPSITE Flexigrid

Moisture Vapour Permeable Adhesive Film Dressing

- Transparent, adhesive film
- Unique wound measurement grid
- Very conformable
- Clinically proven to provide moist wound healing for superficial wounds

6 x 7cm	SN4628	Box/100	\$79.25*	3+	\$73.10*
10 x 12cm	SN4629	Box/10	\$30.75*	3+	\$28.70*
10 x 12cm	SN4630	Box/50	\$86.90*	3+	\$80.15*
12 x 25cm	SN4632	Box/20	\$108.70*	3+	\$101.00*



Indications: Minor wounds, cuts, grazes, superficial minor burns and skin tears. Ideal for use as a secondary dressing and to help prevent pressure injuries


SmithNephew

IV3000

Moisture Responsive IV Catheter Dressing

IV3000* is a moisture responsive, transparent film dressing, specifically designed to meet the needs of IV catheter fixation.

	7 x 9cm	Peripheral	SN4006	Box/100	\$109.05*	3+	\$101.85*
	6 x 7cm	Peripheral/Ported	SN4007	Box/100	\$100.55*	3+	\$92.75*
	6 x 8.5cm	Peripheral	SN4924	Box/100	\$108.05*	3+	\$99.70*
	10 x 14cm	Central/Jugular	SN4925	Box/10	\$40.60*	3+	\$37.90*



Orange-handle Incorporates tape for security of fixation



SmithNephew

OPSITE Flexifix Gentle

Transparent silicone adhesive film roll

OPSITE Flexifix Gentle combines breathability, transparency and up to a seven day wear time to offer a versatile silicone film roll for patients with fragile or at-risk skin. It removes cleanly, without disrupting fragile skin layers or causing patients undue pain.

Indication: to retain primary dressings, fixation of tubing, protect skin from friction and external contamination, for prophylaxis against pressure sores, relieve pain in painful diabetic neuropathy, protect the skin around stoma and under leg bags.

2.5cm x 5m	SN66801195	Roll	\$41.60*
5cm x 5m	SN66801196	Roll	\$72.90*
10cm x 5m	SN66801197	Roll	\$113.95*



SmithNephew