

Setting up a new practice, upgrading or expanding?

We're the right team for the job. With the largest GP Sales Team and our extensive experience, we look after your team from beginning to end. No set-up is too big or too small.

Need a Quote Today?

Give us a call or fill out the form through the QR code to be contacted by one of our specialists. Let us give you a hand, no matter where you are in the process.



team medical supplies

HVF128 TENS Therapy Device

Engage TENS technology for drug-free treatment.

- Block the pain messages to the brain
- Trigger the release of endorphins (natural pain killers)
- Improve blood circulation (as result of repeated muscle contraction and muscle relaxation)

Adapt your treatment according to your preference and pain symptoms.

- Hyper Mode - For relief of stiffness with gradual increase of intensity during treatment
- Soft Mode - For relief with soft stimulation
- Repeat Mode - Repeats your choice of stimulation
- Point and Wide Mode - Treats a small or wide area respectively

Washable & Durable Long-Life Pads, reusable up to 150 times



OMRON

\$270*

DAVHVF128

Your patients can buy all these items via our online shop www.teammed.com.au

OMRON HV-F013 Electronic Nerve Stimulator

The OMRON HV-F013 is a pocket-sized therapy device that targets and alleviates muscle and joint pains with 5 massage modes. It has a compact design and simple one-hand operation – perfect for TENS therapy on the go.

Features:

- 3 body area modes: Arm/Shoulder, Lower Back, and Leg/Foot
- Switch between standard high-frequency and low-frequency treatment as needed
- 2 massage-like modes: Knead and Steady
- Washable & Durable Long-Life Pads, reusable up to 150 times

\$170*

DAVHVF013



Pro Tens Machine

A regular TENS with adjustable pulse width, pulse rate, intensity, modes and an accurate timer of 15, 30 and 60 minutes, will stop output at presenting time.

Channel: Dual, isolated between channels
 Pulse Amplitude: 0-80mA peak to peak into 500 ohm load
 Pulse Rate: Adjustable, from 2Hz ~ 150Hz
 Pulse Width: Adjustable, from 30uS ~ 260uS
 Mode: B, N, M Timer: Selectable, 15, 30, 60 minutes (or 15, 30, constant)
 Wave Form: Asymmetrical Bi-Phasic Square Pulse
 Voltage: 0-40V peak into 500ohm load




ALLPROTENS \$70*

Supplied with Carry Case, 2 x TENS Leads, Operating Manual, 1 x 9V Battery, 1 x Pack of Self-Adhesive Electrodes



Tens Electrodes

5 x 5cm square ALLF35050 Pkt/4 \$14.55*
 5 x 9cm rectangle ALLF35090 Pkt/4 \$11.50*
 5cm round ALLF350R Pkt/4 \$7.25*



Bodisure Massage Gun Mini

Relax your tired muscles on the go with this compact, yet powerful deep tissue massage gun.

The BodiSure Massage Gun Mini has a 2500rpm brushless high torque motor with 4 interchangeable massage heads. Featuring 4 speed adjustable intensity coupled with a cordless and lightweight design with comfortable grip handle, this easy-to-use muscle massager can be used to relax tired muscles and is perfect for post-exercise.



EACH **DAVBMGMINI** \$120

Quick Ultrasound Gel

- Acoustically correct; for the broad range of frequencies used
- Completely aqueous; will not stain clothing or damage transducers
- Water-soluble and odourless
- Hypoallergenic, bacteriostatic, non-sensitising and non-irritating



Quick-Eco Gel Blue

blue colour

250ml **BIO9980020** \$9.80
5Ltr **BIO9980021** \$37.75



Quick-Cristal Gel Clear

no colour, totally transparent

250ml **BIO9980025** \$9.80
5Ltr **BIO9980022** \$37.25

Aquasonic 100 Ultrasound Gel

Most widely used gel for diagnostic and therapeutic ultrasound

250ml **BIBUS0108** \$20.15
5Ltr **BIBUS0150** \$104.65



Hwato Acupuncture Needles (with Guide Tubes)

Hwato Acupuncture Needles are China's finest quality traditional needles. Established 140 years ago, the Hwato factory is the world's oldest, largest and most respected manufacturer of hand-made needles (metal handle needles).

20 x 30mm tube	ACUHT2030	Box/100	\$16.80
20 x 40mm tube	ACUHT2040	5+	\$15.40
22 x 13mm tube	ACUHT2213	10+	\$13.95
22 x 30mm tube	ACUHT2230		
25 x 30mm tube	ACUHT2530		
25 x 40mm tube	ACUHT2540		
25 x 50mm tube	ACUHT2550		
25 x 75mm tube	ACUHT2575		
30 x 30mm tube	ACUHT3030		
30 x 40mm tube	ACUHT3040		
30 x 50mm tube	ACUHT3050		
30 x 60mm tube	ACUHT3060		



Seirin Acupuncture Needles (with Guide Tubes)

The SEIRIN J Type is made with a new technology that allows for smooth painless insertion using less pressure. Features an easy grip, colour coded, lightweight 2.0 x 20mm plastic handle with a soft-touch top.

20 x 30mm tube	ACUSJ2030
20 x 40mm tube	ACUSJ2040
25 x 30mm tube	ACUSJ2530
25 x 40mm tube	ACUSJ2540
25 x 50mm tube	ACUSJ2550
30 x 30mm tube	ACUSJ3030
30 x 40mm tube	ACUSJ3040

Box/100 \$39.00
5+ \$35.65
10+ \$34.20



S+M Instant Ice Pack

230 x 130mm



AXP10130331
Each \$3.65
12+ \$2.60

Instant Ice Pack

230 x 150mm



LRINTI1069
Each \$2.70 40+ \$2.45

S+M Hot/Cold Pack

Medium (110 x 220mm)

AXP10130311
\$6.60

Large (135 x 260mm)

AXP10130312
\$7.25

Heat: Microwavable or hot water



Bodichek Premium Hot/Cold Pack
(with towel bag)

Medium (125 x 270mm)

AXP62630353
\$13.50

Large (175 x 270mm)

AXP62630352
\$22.75

Heat: Microwavable or hot water



Bodichek Reusable Hot/Cold Gel Pack
(clear plastic)

200 x 115mm

Heat: hot water only

AXP62630311
\$8.30



Hot/Cold Gel Pack Reusable

PVC

100 x 150mm

Heat: Microwavable or hot water

LRINTGEL3057
\$3.70



Hotteeze Hand Warmers

- Keeps hands and fingers warm
- Self-heating - air activated
- Non-stick
- Designed to hold in your hand, keep in your pockets, or tuck into your gloves
- Up to 10 hours of heat
- Max temp 63° C

Pad size – 100mm x 70mm

DAVHHW10
Pkt/10 \$13.25



Hotteeze for Feet

- Keeps your toes nice and warm
- Self-heating - air activated
- Non-stick
- Designed to fit into men's and women's shoes
- Up to 6 hours of heat
- Max temp 39° C

Pad size – 96mm x 70mm

DAVHFP5
Pkt/5 \$13.25



Hotteeze Heat Pads

- Keeps hands and fingers warm
- Self-heating - air activated
- Non-stick
- Designed to stick over clothes or to a bodigrip tube
- Up to 10 hours of heat
- Max temp 63° C

Pad size – 132mm x 100mm

DAVHHP10
Pkt/10 \$31.00



ALL-Band 1.5m Exercise Band

All-Band Exercise Band is used in rehabilitation and exercise programs. Available in a handy 1.5 metre length single pack or 25m rolls.

Resistance		1.5m x 14.5cm		25m x 14.5cm	
Light	Yellow	ALL04002401	\$12.10	ALLACBAND25YE	\$48.65
Medium	Red	ALL04002402	\$12.10	ALLACBAND25R	\$53.25
Firm	Green	ALL04002403	\$13.30	ALLACBAND25G	\$67.50
Extra Firm	Blue	ALL04002404	\$13.30	ALLACBANDBL	\$65.95
X-Extra Firm	Black	ALL04002405	\$13.30	ALLACL BANDBK	\$80.55




Exercise and Posture Ball

A great tool to enhance:

Posture, Performance, Flexibility, Strength and Muscle tone

- A durable, flexible and firm ball
- Feature a non slip surface with ribbing for extra security
- An anti-burst ball rated at 550kg

● Red	55cm	ALL03001512	\$40.60
● Orange	60cm	ALL03001527	\$36.70
● Green	65cm	ALL03001533	\$38.65
● Blue	75cm	ALL03001554	\$46.80
● Pearl	85cm	ALL03001568	\$48.75
● Red	90cm	ALL03001569	\$62.60



Foam Rollers

- Made of high density eva foam
- Provides firm support for all body types
- Excellent for balance and motor skill development
- Also use for self massage, de-activating trigger points
- Use for soft tissue work, stretching, massage, stability exercises, core strengthening, injury prevention
- Strong enough to be used by large adults



Small Round	30L x 15Dcm	ALL700461R	\$26.15
Large Half Round	90L x 7.5Dcm	ALL700463H	\$34.85
Small Half Round	30L x 7.5Dcm	ALL700461H	\$18.15

Crutches

Forearm Medium	890-1200mm	AUS4010
Forearm Large	890-1370mm	AUS4020

\$62.50*

Underarm Small	1400-1600mm	AUS4050
Underarm Medium	1570-1780mm	AUS4051
Underarm Large	1570-1980mm	AUS4052

\$64.25*



Wheelchair Patient mover

An essential in every clinic

- 110kg user weight capacity
- Seat 46cm W x 41cm D
- Frame: Single cross brace, powder coating
- Front wheel: 8", PVC tires
- Rear wheel: 12" PU tires
- Armrest: Flip back with desk length
- Footrest: Swing away footrest

TMF18WCM **\$385***



Wheelchair - Standard

A versatile and practical wheelchair

- 110kg user weight capacity
- Seat 45cm W x 40cm D
- Frame: Single cross brace, powder coating
- Front wheel: 8"x6", rubber tires
- Rear wheel: 24" spoke wheel, pneumatic tires
- Armrest: Flip back with desk length
- Footrest: Swing away footrest

TMF18WCS **\$365***



Wheelchair - Bariatric

High Capacity wheelchair

- 160kg user weight capacity
- Seat 55cm W x 40cm D
- Frame: Single cross brace, powder coating
- Front wheel: 8"x2", PU tires
- Rear wheel: 24" spoke wheel, PU tires
- Armrest: Flip back with desk length
- Footrest: Swing away footrest
- Seat and backrest: Nylon with sponge

TMF22WCB **\$495***

