

Soffban Natural
Orthopaedic Padding



- Conforms well for easy shaping
- Natural fleece absorbs sweat and exudate
- Ideal for plaster casts
- 2.7m



5cm	SN7147200	Pkt/12	\$31.80	2+	\$29.70
7.5cm	SN7147300	Pkt/12	\$34.40	2+	\$31.75
10cm	SN7146900	Pkt/12	\$43.85	2+	\$40.45
15cm	SN7147000	Pkt/12	\$51.20	2+	\$47.20

Soffban Synthetic
Water-Resistant Orthopaedic Padding



- Water-resistant
- Ideal for synthetic casts
- 2.7m



5cm	SN36000025	Pkt/12	\$31.40	2+	\$29.10
7.5cm	SN36000026	Pkt/12	\$34.40	2+	\$31.75
10cm	SN36000027	Pkt/12	\$43.85	2+	\$40.45
15cm	SN36000028	Pkt/12	\$51.20	2+	\$47.20



Delta-Dry
Waterproof Cast Liner



- Water resistant cast padding (2.4m)
- Allows patients to shower, swim and bathe like normal
- Better long term water resistance
- Offers more under cast protection
- Eliminates the need for stockinette



5cm	SN7344300	Box/12	\$272.65
7.5cm	SN7344301	Box/12	\$328.55
10cm	SN7344302	Box/12	\$379.55

Velband

Absorbent Orthopaedic Padding (2.7m)

Velband is a hypoallergenic underbandage with a unique three-layer structure for easy application and long-term comfort. It provides excellent cushioning and moisture management under plaster of Paris or synthetic casts , avoiding irritation and pressure sores by keeping the skin dry and comfortable.



7.5cm	SN7231501	Pkt/12	\$33.30	2+	\$30.70
10cm	SN7231502	Pkt/12	\$40.05	2+	\$36.90
15cm	SN7231504	Pkt/12	\$57.80	2+	\$54.00

Pure Cotton

Orthopaedic under-bandages

Made from high quality, long staple cotton fibres. Designed and made to have some self-adhesion property and appropriate tensile strength, resistant to abrasion and plaster strike through. 2.7m



5cm	MGT30700	Pkt/12	\$15.25	3+	\$13.95
7.5cm	MGT30701	Pkt/12	\$17.65	3+	\$15.90
10cm	MGT30702	Pkt/12	\$29.95	3+	\$26.85
15cm	MGT30703	Pkt/12	\$41.60	3+	\$37.35

Protouch

Synthetic and Cotton Stockinette

- Soft and comfortable providing good cushioning and comfort to the patient.
- Stretches to three times its normal width to conform to an arm, leg or body.

Protouch Synthetic Stockinette

- Made of synthetic fibres to react moisture away from the skin for rapid evaporation

5cm	SN301002	Roll	\$34.85
7.5cm	SN301003	Roll	\$41.45
10cm	SN301004	Roll	\$69.15



Protouch Cotton Stockinette

- Ideal material for an initial undercast layer
- Made of a high-grade knitted cotton

5cm	SN307002	Roll	\$34.85
7.5cm	SN307003	Roll	\$31.60
10cm	SN307004	Roll	\$37.50



Gypsona

Plaster of Paris Bandage

Gypsona bandages combine fast setting, high quality plaster of Paris with the conformable and stable interlocked weave of leno cloth. The bandage demonstrates minimum plaster loss whilst the rich creamy texture allows easy handling.

5cm	SN87060527	Box/12	\$66.25	2+	\$61.10
7.5cm	SN87067535	Box/24	\$142.20	2+	\$131.15
10cm	SN87061035	Box/24	\$184.10	2+	\$169.75
15cm	SN87061535	Box/12	\$129.90	2+	\$119.85




Delta-Lite Plus

Fibreglass Casting Bandage

- Easy to apply
- Greater conformability
- High lamination: smooth finished cast which is strong, long lasting and lightweight

5cm	SN7345801	Roll	\$13.05	10+	\$12.20
7.5cm	SN7345802	Roll	\$14.75	10+	\$13.75
10cm	SN7345803	Roll	\$18.75	10+	\$17.45
12.5cm	SN7345804	Roll	\$21.25	10+	\$19.85



Also available in following colours: red, blue, purple, black, pink, green, mixed
Note that we only sell them in full cartons of 10. Price is equivalent to 10+ price of standard white. Colours are not a standard stock line

Delta-Cast Elite

Polyester Casting Tape

Delta-Cast Elite has been designed with ease of application and patient comfort in mind. Our patented fibreglass-free casting material combines strength and rigidity with soft edges for enhanced patient comfort.

5cm	SN7227900	Roll	\$14.75	10+	\$13.70
7.5cm	SN7227901	Roll	\$16.70	10+	\$15.50
10cm	SN7227902	Roll	\$20.35	10+	\$18.85
12.5cm	SN7227903	Roll	\$24.75	10+	\$22.95

- Smooth cast with soft edges
- Fibreglass free
- Easy to apply and remove
- Highly conformable



Also available in following colours: red, blue, purple, black, pink, green, mixed
Note that we only sell them in full cartons of 10. Price is equivalent to 10+ price of standard white. Colours are not a standard stock line



3M™ Scotchcast™ Plus

Synthetic Casting Tape

Consists of a knitted fiberglass fabric impregnated with polyurethane resin. Exposure to moisture or water initiates a chemical reaction which causes the tape to become rigid. Finished casts produced are lightweight, strong, radiolucent, and porous. Scotchcast Plus casting tape contains a lubricating agent in the resin to provide a slippery feel when exposed to water.

5cm	3M82002	Roll	\$13.50	10+	\$12.30	10cm	3M82004	Roll	\$17.80	10+	\$16.40
7.5cm	3M82003	Roll	\$15.90	10+	\$14.50	12.5cm	3M82005	Roll	\$20.50	10+	\$18.75



Finger Stax splints



Sz1	MIAS01	Sz5	MIAS05
Sz2	MIAS02	Sz6	MIAS06
Sz3	MIAS03	Sz7	MIAS07
Sz4	MIAS04	Each	\$4.25

Finger splints

Aluminium



13mm x 480mm	ROW4813	\$6.00	2+	\$4.70
19mm x 480mm	ROW4819	\$6.45	2+	\$5.10
24mm x 480mm	ROW4824	\$6.85	2+	\$5.40
32mm x 480mm	ROW4832	\$7.80	2+	\$6.25

Sam Splint

Lightweight, reusable and very portable

The Sam Splint Original is an extremely pliable multi-purpose splint and is built from a thin core of aluminium alloy, sandwiched between two layers of closed cell foam. Bent into various shapes and curves, it becomes extremely strong and supportive for any fractured or injured limb.




9.5cm x 4cm	FERSVPSSF	\$4.40
91cm x 10mm	FERSVPSS36	\$26.50



Collar & Cuff Arm Sling

6m ROWCC06W
\$19.10



Slings

Disposable Slings

AXP10106311
110 x 110 x 155cm
Each \$1.20
12+ \$1.00



Cotton Slings


AXP10106211
110 x 110 x 155cm
Each \$1.55
12+ \$1.40



Safety Pins
Assorted Sizes

LIVSP12L

Pkt/12 \$0.60



Actimove Sling

Collar and Cuff Support

The padded strap and Y-tab closure system provide effective vertical support to the collarbone and upper limb area

- Extra patient comfort with padded strap
- Y-tab hook system facilitates secure clinical positioning
- Convenient dispenser box presentation
- 28 Y-tabs are included in each box



essity

Actimove
Feel good. Move better.

5cm x 12m
SN7285918
Box/2rolls \$66.95



Actimove Mitella

Comfort Arm Sling

The envelope sling and neck strap of the Actimove Mitella provide effective vertical support to the arm.

- Extra patient comfort with neck pad
- Adjustable neck strap, with Y-tab hook system facilitates secure clinical positioning
- Soft thumb loop maintains sling placement

Actimove
Feel good. Move better.

essity

Small SN7281961
Medium SN7281962
Large SN7281963
X-Large SN7281964

\$43.00

Jobst® Travel Socks

Medical Compression Stockings



See page 412



Actimove

Lower Leg Walker

The Actimove® Lower Leg Walker combines comfort and support. Its low profile contoured construction ensures a superior anatomical fit.

x-sml SN7627035
sml SN7627036
med SN7627037
lrg SN7627038
x-lrg SN7627039

\$63.30

Actimove
Feel good. Move better.



Standard			UK	US	EU
76270-35	XS	M	<3.5	<4.5	<36.5
		F	<3	<6	<36
76270-36	S	M	4-6.5	4.5-7	36-39.5
		F	3.5-5.5	6-8	36.5-38.5
76270-37	M	M	7-9.5	7.5-10	40-44
		F	6-8.5	8.5-11	39-43
76270-38	L	M	10-12	10.5-12.5	44.5-45.5
		F	9-11	11.5-13.5	44-46
76270-39	XL	M	>12	>12.5	>46.5
		F	>11	>13.5	>46

Cervical Collar

Adjustable Extraction Collar

Moulded from high density polyethylene with foam padding and a patented chin support. This collar restricts cervical spine flexion, extension and rotation to promote patient recovery.

The one-piece design enables efficient storage where space is limited. The multi-height adjustment fits almost all adults (four positions) and paediatrics (three positions).

The one-piece design can easily be used; simply measure patient, adjust and lock selected height and apply. For either a supine or seated position, application is quick and easy.

ADTMSEX700 \$44.50*



Cervical Collar

All 48cm length

Small 7.5cm SUTSM400
Medium 9cm SUTSM401
Large 10cm SUTSM402

\$12.50*

