

## Dressing Strips and Spots non-sterile (8ply)

This dressing strip is designed to protect minor wounds. It is water resistant and helps protect the wound from water and germs.



### Plaster Strips 70mm x 17mm

BP0113100BA

Box/100 \$9.20



### Spots 20mm

BP0114100BA

Box/100 \$7.50



### Fabric Strips 70mm x 17mm

BP0119100BA

Box/100 \$10.80



### Curity Spots Ø22mm

COV6403LF

Box/100 \$5.60



### Curaplast® Sensitive Injection Dressing

20 x 40mm

Easy Dispenser box

LOH30625

Roll/250 \$22.25



### Curaplast® Sensitive Adhesive Round Dressing

Ø16mm

Easy Dispenser Box

LOH30626

Roll/100 \$11.85



### Nichiban IV Pressure Pad

- Overall Product Ø36mm
- IV Pressure Pad Ø16mm
- White colour

DEF929

Box/100 \$32.35



### IV Pressure Pad

- Overall Product Ø36mm
- IV Pressure Pad Ø16mm
- Tan colour

MIDAWIVPAD

Box/100 \$48.30



### Band-Aid - Advanced Healing

Plaster seals the wound to act like an instant scab, and absorbs wound fluid to create an optimal moist healing environment for fast wound healing.

- Hydrocolloid promotes fast healing
- Acts like an instant scab
- 100% waterproof
- Hypoallergenic

Regular 25 x 66mm **XXB035**  
Box/10 \$15.10

Large 44 x 70mm **XXB036**  
Box/6 \$15.90



### Band-Aid Strips

JJ37206156

Box/100 \$11.40



### Elastoplast Sensitive

Dressing lengths  
(6 x 10cm)

BDF46040

Box/10 \$12.10



### Nexcare Tattoo Strips "Cool" Collection

Waterproof bandage  
stays on skin in  
the shower or pool

3M5942

Box/20 \$11.60



## Leukoplast Elastic Spots

Latex Free



**Leukoplast®**  
Wound care in best hands

2.2cm diameter

BSN7645523


Box/100 \$11.00 12+ \$10.45

## Leukoplast Universal

Plastic First Aid Dressing

Leukoplast Standard dressings are high quality, perforated plastic film strips ideal for first aid or minor wound management.

Latex-free



7.2cm x 1.9cm

BSN7646201

Box/100 \$8.05 12+ \$7.40



## Leukoplast Strong Fabric Dressings

Fabric First Aid Dressing

- Ideal for use with cuts, lacerations and skin loss wounds
- Absorbent, low-adherent pad backed by a high quality stretch fabric
- Practical and highly conformable
- Strong, porous adhesive ensures a strong and lasting fixation

First aid dressing strips

2.2cm x 7.2cm

BSN7646007

Box/100 \$11.35

12+ \$10.55

BSN7646008

Box/50 \$6.20

6+ \$5.75



Finger Tip Dressing Small  
45 x 50mm

BSN7646012

Box/50 \$24.55

3+ \$22.56




Finger Tip Dressing Large  
50 x 63mm

BSN7646011

Box/100 \$61.70

3+ \$57.50



Knuckle Dressing  
Strong Sterile 38 x 75mm Tan

BSN7646010

Box/100 \$81.35

3+ \$74.90



## Leukoplast Strong Dressing Lengths

4cm x 1m BSN7645907 Roll \$5.15 5+ \$4.65

6cm x 1m BSN7645905 Roll \$6.00 5+ \$5.45

6cm x 5m BSN7645901 Roll \$25.00 3+ \$22.80

8cm x 5m BSN7645902 Roll \$32.10 3+ \$29.25

Leukoplast Strong Dressing rolls can be cut to the exact length required, making it highly cost-effective.



## Leukoplast Barrier

Waterproof Barrier Dressing

Leukoplast Barrier is a waterproof, microporous first aid dressing that breathes and is an effective bacterial and viral barrier.

2.2cm x 6.3cm BSN7645201

Box/100 \$50.30

3+ \$45.70



2.2 x 7.2cm BSN7645208

Box/100 \$50.60

3+ \$45.85



3.8cm x 3.8cm BSN7645202

Box/100 \$57.95

3+ \$53.35



5cm x 7.2cm BSN7645205

Box/100 \$115.20

3+ \$106.05



## Tasmania Now also Covered by Mediport

At Team Medical we now also supply the Tasmanian region with Mediport. This allows us to ship vaccines in Mediport's validated refrigeration vehicles, which means that instead of polystyrene eskies, our Tasmanian customers now also receive their vaccines in recycled boxes. More and more regions are now being serviced this way. Unfortunately, given the vastness of Australia, this service is not available Australia-wide.

We will keep looking for other ways to do our bit for the environment.

Thank you to all our customers for adopting this new way of shipping.

78% of our vaccine orders in 2023 have been shipped via Mediport



## Cotton Tips and Applicators

7.5 single end sterile - pkt/5	MGT21893	Box/100	\$26.30	3+	\$23.80
15cm single end sterile - pkt/2	MGT21891	Box/100	\$20.10	3+	\$18.20



## Cotton balls

Cotton wool balls sterile	MGT02222P	Pkt/5	\$0.25	300+	\$0.20
Non-woven balls sterile	MGT00005	Pkt/5	\$0.30	500+	\$0.25
Cotton wool balls Swisspers	SN12260	Pkt/160	\$4.20	6+	\$4.00
Cotton wool balls large	SN2902021	Ctn/3000	\$81.50	3+	\$75.05
Cotton wool balls small	SN2903141	Ctn/4000	\$61.70	3+	\$56.80
Cotton wool roll 375g	AXP10102411	Each	\$15.15	6+	\$14.15



### Eye pad

Eye Pad with Crimped Edges  
(Individually wrapped in  
Dispenser Box) - Sterile

MGT19111D

Box/25	\$7.20
3+	\$6.30



### Nexcare Opticlude

Nexcare Eye Patch Junior  
3M1537

Box/20	\$48.25
--------	---------

Nexcare Eye Patch Regular  
3M1539

Box/20	\$33.55
--------	---------



### Combine pieces



10 x 12cm	BP013410BA	Each	\$0.40	100+	\$0.30
10 x 22cm	BP013422BA	Each	\$0.50	100+	\$0.40
20 x 22cm	BP013420BA	Each	\$0.65	100+	\$0.55



10 x 12cm	MGT09888P	Each	\$0.40	250+	\$0.30
10 x 20cm	MGT09889P	Each	\$0.50	150+	\$0.40
20 x 20cm	MGT09890P	Each	\$0.65	75+	\$0.55



### Packing Gauze



1.25cm x 1m	MGT11111	Each	\$1.00	100+	\$0.70
2.5cm x 1m	MGT11112	Each	\$0.80	100+	\$0.60



### Green Gauze

10 x 10cm (non-sterile)	MGT18888B	Pkt/100	\$7.00
10 x 10cm (sterile)	MGT17888	Pkt/5	\$1.15



### Ear Wick

MEROCEL Ear Wicks, Pope Oto Wick

- Ideal for treating otitis externa
- Expands to conform to ear canal when hydrated with otic solutions and local anesthesia

XOM1425008	15mm
Pkt/5	\$49.25



### X-Ray Gauze

10 x 10cm	MGT15666B	Pkt/5	\$1.00
5cm x 1m	MGT15112	Each	\$1.20








**Gauze Swabs non-sterile (8ply)**

This soft and highly absorbent gauze can be used to clean a wound before it is dressed.

<b>5cm</b>	<b>BP010805BA</b>			
Pkt/100	\$2.30	200+	\$2.10	
<b>7.5cm</b>	<b>BP010875BA</b>			
Pkt/100	\$3.50	100+	\$3.20	
<b>10cm</b>	<b>BP010810BA</b>			
Pkt/100	\$4.50	50+	\$4.05	



**Gauze non-sterile**

**7.5cm**  
**SN36361349**

Pkt/100	\$4.60
100+	\$3.55



**Gauze Squares Sterile**



<b>7.5cm x 7.5cm</b>	<b>BP010975DA</b>	Box/50	\$22.60	6+	\$20.45
<b>10cm x 10cm</b>	<b>BP010910BA</b>	Pkt/5	\$0.50	500+	\$0.45
<b>7.5cm x 7.5cm</b>	<b>BP010975YA</b>	Pkt/3	\$0.30	500+	\$0.25
<b>7.5cm x 7.5cm</b>	<b>BP010975XA</b>	Pkt/5	\$0.45	500+	\$0.30



**Sterile Gauze in dispenser pack**

<b>10cm x 10cm</b>	(pkt/3) dispenser pack		
Box/50	\$29.80	<b>MGT11559D</b>	
8+	\$24.30		



**Non-woven Gauze sterile**

<b>7.5 x 7.5cm</b>	<b>MGT47600A</b>	Pkt/3	\$0.25	500+	\$0.20
--------------------	------------------	-------	--------	------	--------

<b>7.5 x 7.5cm</b>	(pkt/3) dispenser pack		
Box/30	\$8.90	<b>MGT47600C</b>	
6+	\$7.60		



**Sterile Gauze Handy**

<b>5 x 5cm</b>	<b>SN7111705</b>	Pkt/100	\$6.55	10+	\$6.10
<b>7.5 x 7.5cm</b>	<b>SN5673</b>	Pkt/100	\$11.35	10+	\$10.45
<b>10 x 10cm</b>	<b>SN7111706</b>	Pkt/100	\$19.45	10+	\$17.90



**Topper 8 Gauze**

Topper 8 is composed of nonwoven viscose rayon and polyester. These cost effective low-linting swabs with high absorbency are ideal for use in all general wiping and cleansing procedures associated with wound management.

- Ribbed structure provides a textured surface to facilitate debris removal
- The non-woven fabric permits high absorbency without linting, unlike gauze swabs

**Non Sterile**

<b>5 x 5cm</b>	<b>JJM50405</b>	Pkt/100	\$6.15
<b>7.5 x 7.5cm</b>	<b>JJM53407N</b>	Pkt/200	\$19.70

**Sterile**

<b>5 x 5cm (5's)</b>	<b>JJTS8055</b>	Pkt/25	\$21.40
<b>7.5 x 7.5cm (2's)</b>	<b>JJTS8072</b>	Box/50	\$24.15
<b>10 x 10cm</b>	<b>JJTS8102</b>	Box/75	\$56.15



### Non-adhesive Absorbent Dressing

This non-adherent dressing helps wounds heal with a specially designed absorbent pad. Keep it in place with a bandage, adhesive tape or waterproof film dressing.

5cm x 5cm	BP010105BA	Box/100	\$24.80	3+	\$22.60
10cm x 10cm	BP010110BA	Box/100	\$58.25	3+	\$55.00
10cm x 20cm	BP010120BA	Box/50	\$52.85	3+	\$49.60



### Non-adhesive Lite Absorbent Dressing

This non-adherent dressing is a lighter version with approximately 20% less cotton. It is typically used for wounds with less wound fluid to manage.

7.5cm x 10cm	BP010210BA	Box/50	\$18.30	3+	\$16.15
7.5cm x 10cm	BP010205BA	Box/100	\$19.40	3+	\$17.20
7.5cm x 10cm	BP010220BA	Box/50	\$24.80	3+	\$22.60



### Interpose®

Non-adherent dressing, high absorption

5 x 5cm	MGT30381	Box/100	\$21.85	4+	\$18.50
10 x 10cm	MGT30382	Box/100	\$55.95	4+	\$50.50
10 x 20cm	MGT30383	Box/100	\$94.65	2+	\$83.25



### Cutilin

Cutilin is a low-adherent dressing for exuding wounds that can be cut to suit, without linting.

- Step 1: Clean the wound with saline solution. Do not dry the wound bed.
- Step 2: Cutilin can be secured with an adhesive tape or bandage.
- Step 3: Cutilin can be left in place until the capacity of the pad is exhausted, which can be seen when exudate reaches the edges of the dressing

5 x 5cm	SN76240	Box/100	\$42.50	3+	\$39.15
10 x 10cm	SN76242	Box/100	\$104.70	3+	\$96.45



Smith+Nephew

### Melolite


Lightly absorbent low-adherent dressing

- Double sided low-adherent wound dressing
- Absorbent pad
- Highly conformable to all wound sites
- Either side of the dressing can be placed directly onto the wound site

7.5 x 5cm	SN4811	Box/100	\$27.50	3+	\$25.30
7.5 x 10cm	SN4812	Box/100	\$43.45	3+	\$40.00
7.5 x 20cm	SN4814	Box/100	\$69.85	3+	\$65.50



Smith+Nephew



## Melolin

### Highly absorbent low-adherent dressing

Melolin consists of a highly absorbent cotton and acrylic fibre pad which is heat bonded on one side to a very thin perforated polyester film. The film side of the dressing is placed next to the wound.

- Rapid drainage of exudate, reducing trauma to healing tissue
- Comfortable and minimises pain on removal
- Easy to cut to shape

5 x 5cm	SN66974940	Box/100	\$45.85	3+	\$42.25
10 x 10cm	SN66974941	Box/100	\$113.45	3+	\$104.50
10 x 20cm	SN66974939	Box/100	\$201.00	3+	\$185.20



Smith+Nephew

## Zetuvit

### Non-adherent absorbent dressing

Zetuvit is a sterile single use absorbent dressing used for the treatment of exuding wounds, acute or chronic, and for upholstery and cushioning of wounds. It's high absorbency reduces the risk of maceration to the wound and surrounding tissue.

10 x 10cm	HTM413701	Box/25	\$74.55	3+	\$68.05
10 x 20cm	HTM413702	Box/25	\$79.40	3+	\$72.45
20 x 20cm	HTM413703	Box/15	\$73.85	3+	\$67.40



HARTMANN

## Zetuvit® Plus

### Super absorbent dressing

For highly exuding wounds, offering soft cushioning, excellent absorption and retention capacity.

Zetuvit® Plus is the no. 1 superabsorbent dressing. Suitable for the treatment of superficial, high to very highly exuding wounds such as acute wounds (traumatic wounds, post-operative wounds, lymphatic wounds, arterial leg ulcers) and chronic wounds (pressure injuries, leg ulcers, tumour wounds). Can also be used under compression bandages.

10 x 10cm	HTM413710	Box/10	\$41.90	3+	\$38.25
10 x 20cm	HTM4137113	Pkt10	\$55.35	3+	\$50.50
15 x 20cm	HTM413712	Box/10	\$68.95	3+	\$62.95
20 x 25cm	HTM4137133	Pkt10	\$94.10	3+	\$85.90
20 x 40cm	HTM4137150	Box/10	\$135.65	3+	\$123.80



TIP: White side to wound


APPLICATION:



HARTMANN

## Zetuvit®

### Plus Silicone Border



10 x 10cm	HTM4139100	Box/10	\$138.30	2+	\$131.85
12.5 x 12.5cm	HTM4139200	Box/10	\$149.95	2+	\$143.00
17.5 x 17.5cm	HTM4139300	Box/10	\$220.15	2+	\$205.90
16 x 26cm	HTM413904	Box/10	\$417.20	2+	\$390.20
20 x 20cm	HTM413905	Box/10	\$500.65	2+	\$468.30

### Siliconised Super Absorbent Dressing, Effectiveness Redefined

Zetuvit® Plus Silicone Border is a silicone SAP dressing suitable for the treatment of a wide range of exuding wounds.

- Versatile - Simple to choose sizes and shapes, simple to use on a variety of wounds
- Effective - Improves healing conditions by balancing the microclimate, reducing inhibitors to wound healing and protecting the wound.
- Patient-friendly - Improves patients' quality of life by reducing risks, staying in place for up to 7 days, and being comfortable.



HARTMANN



## Telfa "Ouchless" Non-Adherent Dressings

Superior "Ouchless" dressing won't disrupt healing tissue by sticking to wound. Ideal as a primary dressing for lightly draining wounds.

<b>5 x 7.5cm</b>	<b>COV1961</b>	Box/100	<b>\$25.10</b>	3+	<b>\$21.95</b>
<b>7.5 x 10cm</b>	<b>COV2132</b>	Box/100	<b>\$29.95</b>	3+	<b>\$26.20</b>



## Telfa "Ouchless" Adhesive Dressings

Highly absorbent cotton fabric bonded on both sides with perforated non-adherent film. Won't disrupt healing tissue by sticking to wound.

<b>5 x 7.5cm</b>	<b>COV6017</b>	Box/100	<b>\$24.70</b>	3+	<b>\$20.85</b>
<b>7.5 x 10cm</b>	<b>COV7643</b>	Box/100	<b>\$41.90</b>	3+	<b>\$36.65</b>

## Mesorb®



### Highly Absorbent dressing with a fluid-repellent backing

Mesorb is a soft highly absorbent dressing that can retain large amounts of wound exudate. The fluid-repellent strikethrough barrier helps prevent exudate from soiling clothing and bed linen.

- High absorption and retention capacity
- Conformable and comfortable for patients

<b>10 x 13cm</b>	<b>MOL677001</b>	Box/50	<b>\$74.55</b>
<b>15 x 23cm</b>	<b>MOL677601</b>	Box/50	<b>\$131.90</b>



## Jelonet

### Paraffin Gauze Dressing

JELONET is a soft paraffin dressing made from open-weave gauze. The gauze has interlocking threads, which minimise fraying when the dressing is cut to shape. JELONET is soothing and low-adherent, and allows the wound to drain freely into an absorbent secondary dressing.

Indications: Burns, scalds, lacerations, abrasions and leg ulcers.

<b>5cm x 5cm</b>	<b>SN7403</b>	Box/50	<b>\$96.10</b>	3+	<b>\$88.55</b>
<b>10cm x 10cm</b>	<b>SN66007450</b>	Box/50	<b>\$98.75</b>	3+	<b>\$90.95</b>
<b>10cm x 10cm</b>	<b>SN7404</b>	Box/10	<b>\$19.60</b>	3+	<b>\$18.15</b>



## Smith+Nephew



## Bactigras


### Medicated Tulle Gras

BACTIGRAS is a gauze of leno weave impregnated with White Soft Paraffin BP containing 0.5% Chlorhexidine Acetate BP in a sterile presentation.

BACTIGRAS is effective against a wide range of gram-positive and gram-negative wound pathogens, including MRSA.

BACTIGRAS is for use as a topical treatment only. It is indicated for a wide range of wounds where there is a risk of infection. For already infected wounds, BACTIGRAS may be used in conjunction with systemic antibiotics.

The BACTIGRAS tulle gras dressing, made from an open-weave gauze, has low adherence and allows the wound to drain freely into an absorbent secondary dressing.



Indications: Suitable for covering wounds such as minor burns, lacerations, abrasions, graft sites and leg ulcers.

<b>5 x 5cm</b>	<b>SN7456</b>	Box/50	<b>\$101.60</b>	3+	<b>\$93.60</b>
<b>10 x 10cm</b>	<b>SN7457</b>	Box/10	<b>\$25.15</b>	3+	<b>\$23.20</b>
<b>10 x 10cm</b>	<b>SN66003650</b>	Box/50	<b>\$124.30</b>	3+	<b>\$114.50</b>

## Smith+Nephew



## Cuticerin

### Low-Adherent, Impregnated Gauze Dressing

An impregnated, non-aqueous, non-linting dressing ideal as a wound contact layer for burns, abrasions and other wounds with superficial skin loss. It reduces the risk of granulation tissue growing through the dressing, decreasing trauma on removal.

Indications: For treatment of extensive, exuding wounds – burns, abrasions, radiation injuries, skin graft and donor sites.



## Smith+Nephew

<b>7.5cm x 7.5cm</b>	<b>SN66045560</b>	Box/50	<b>\$102.30</b>
		3+	<b>\$94.20</b>