

## Mepitel®



### Wound contact layer dressing with Safetac® on both sides



Mepitel® is a gentle, effective wound contact layer using Safetac®. The open mesh design allows exudate to pass into an outer absorbent dressing. The Safetac layer prevents the outer dressing from sticking to the wound and ensures atraumatic dressing changes. The Safetac layer seals around the wound edges, preventing exudate leaking onto the surrounding skin, thus minimising the risk of maceration.

- Minimises pain and trauma at dressing changes<sup>1</sup>
- Can remain in place for up to 14 days
- Maintained product properties over time - leaves no residues and does not dry out
- Transparent for easy wound inspection during application and during wear
- Low potential for skin irritation and allergy

5 x 7.5cm	MOL290510	Box/10	\$123.10*	3+	\$111.25*
7.5 x 10cm	MOL290710	Box/10	\$202.10*	3+	\$182.70*



## Mepitel® One



### Wound contact layer dressing with Safetac® on one side



An innovative one-sided Safetac® wound contact layer, Mepitel One lets you give the healing process the time it needs as it can be left in place for up to 14 days. Avoiding the need for dressing removal helps minimise patient pain and unnecessary disturbance to the healing process.

- Minimises trauma to the wound and pain at dressing changes<sup>1</sup>
- Adheres gently to the surrounding skin on removal and minimises the risk of maceration
- Mepitel One supports optimal undisturbed healing
- Built-in directional arrows, perfect to guide nurses on the direction it should be removed

What is the technical difference between Mepitel One and Mepitel? Mepitel One consists of a polyurethane film coated on one side with Safetac technology. Mepitel is a polyamide-net, coated on both sides with Safetac technology.

6 x 7cm	MOL289170	Box/5	\$70.95*	3+	\$64.15*
9 x 10cm	MOL289270	Box/5	\$123.05*	3+	\$111.25*




## Atrauman®


### Atraumatic wound contact layer

Atrauman is a non-adherent wound contact layer suitable for the treatment of superficial acute and chronic wounds of any type. It is particularly suitable for preventing the secondary dressing from sticking to the wound bed.

5 x 5cm	HTM449550	Box/50	\$70.70
7.5 x 10cm	HTM499553	Box/50	\$98.60
10 x 20cm	HTM499536	Box/30	\$105.80



#### APPLICATION:



## Atrauman® Silicone

### Silicone wound contact layer


Atrauman® Silicone is a single-use silicone wound contact layer for the treatment of superficial, acute and chronic wounds.


The dressing can be used as protection for the wound bed and edges which can stay in place for up to 7 days.s.

5 x 7cm	HTM499568	Box/10	\$90.05
7.5 x 10cm	HTM499561	Box/5	\$77.45



#### APPLICATION:






## ActivHeal® SILICONE WOUND CONTACT LAYER



**ActivHeal® Silicone Wound Contact Layer** prevents the dressing from sticking to the wound and reduces pain during dressing changes.



- ✓ Highly conformable
- ✓ Allows passage of exudate
- ✓ Gentle yet secure adhesion
- ✓ Can remain in situ for 14 days
- ✓ Minimises epidermal stripping<sup>1</sup>
- ✓ Minimises pain at dressing changes<sup>1</sup>
- ✓ Protects the periwound
- ✓ Maintains a moist wound environment



**BEST BUY**

**Indications for use:**

ActivHeal® Silicone Wound Contact Layer dressing is indicated for use on nil to heavily exuding chronic and acute wounds (with appropriate secondary dressing). The dressing may be used throughout the healing process on the following wounds:

- ✓ Surgical Incisions
- ✓ Skin Tears
- ✓ Pressure Ulcers
- ✓ Venous Leg Ulcers
- ✓ Arterial Leg Ulcers
- ✓ Partial Thickness Burns
- ✓ Partial Thickness Burns Skin Grafts
- ✓ Can be used in conjunction with Negative Pressure Wound Therapy (NPWT)

Nil > Low > Moderate > Heavy Exudate



5 x 7.5cm	ADM9001241 Box/10	\$64.80*	3+	\$59.80*
10 x 10cm	ADM9001258 Box/10	\$102.05*	3+	\$94.10*
10 x 20cm	ADM9001265 Box/10	\$183.85*	3+	\$166.75*
15 x 15cm	ADM9001272 Box/10	\$150.45*	3+	\$136.45*

Reference: 1. Barrett,S.(2012) BJN 21 (21):1271-7.  
ActivHeal is a registered Trademark of Advanced Medical Solutions Limited.

## Adaptic

### Non-Adherent Dressing

Indicated for dry to highly exuding wounds where adherence of dressing and exudate is to be prevented.


A primary dressing made of knitted cellulose acetate fabric and impregnated with a specially formulated petrolatum emulsion.

When to use it:

1st and 2nd degree burns, Lacerations and abrasions, Grafts, Venous ulcers, Pressure ulcers, Nail extractions, Eczema, Staples and suture lines, Surgical incisions and reconstructive procedures

<b>7.6 x 7.6cm</b>	<b>JJ2012</b>	Box/50	<b>\$69.40</b>	3+	<b>\$64.25</b>
<b>7.6 x 20cm</b>	<b>JJ2015</b>	Box/24	<b>\$84.70</b>	3+	<b>\$79.10</b>

**PROTECTS THE WOUND  
MINIMISES PAIN  
AND MACERATION**



## Adaptic Touch

### Non-Adhering Silicone Dressing

A flexible, open-mesh silicone wound contact layer

- Designed to facilitate fluid transfer to secondary dressing and minimise both adherence and pain at dressing change
- Protects the wound and works as a primary dressing
- Helps prevent adherence of the secondary dressing to the wound
- Minimises patient pain and trauma at dressing change
- Allows free passage of exudate, reduces the risk of exudate pooling and maceration

<b>5 x 7.6cm</b>	<b>KCITCH501</b>	Each	<b>\$9.10</b>	3+	<b>\$8.90</b>
<b>7.6 x 11cm</b>	<b>KCITCH502</b>	Box/10	<b>\$162.10</b>	3+	<b>\$154.65</b>

**ALL YOU COULD WISH FOR  
IN A SILICONE CONTACT LAYER**



## Adaptic Digit

### Non-Adhering Digit Dressing

Made of a unique, tacky, silicone-coated fabric covered with a tubular bandage. It is specifically designed for dressing fingers and toes.

ADAPTIC™ DIGIT does not adhere to the wound. It is conformable and non-bulky, therefore allowing maximum digit flexibility. It does not adhere to the wound. Easy to apply and remove.

Indicated to manage a range of mild to severe digit injuries, such as:  
Lacerations and abrasions, traumatic injuries, nail extractions, suture lines

<b>Small</b>	<b>∅2cm</b>	<b>KCIMAD003</b>	Box/10	<b>\$116.20</b>
<b>Medium</b>	<b>∅2.4cm</b>	<b>JJMAD013</b>	Box/10	<b>\$117.30</b>
<b>Large</b>	<b>∅2.8cm</b>	<b>JJMAD023</b>	Box/10	<b>\$120.55</b>
<b>Toe</b>	<b>∅2.8cm</b>	<b>KCIMAD062</b>	Box/10	<b>\$108.25</b>

**COMFORTABLE  
EASY TO APPLY  
AND REMOVE**



# GOING ABOVE AND BEYOND

Our entire team is committed to providing a positive customer experience, every time you order. Everything from the way you place your order to the way our warehouse team packs your order is optimised to provide the ultimate customer experience.

By listening to our customers' feedback and empowering our employees, we continually improve the customer experience through innovation.



## Soft & Flexible Fabric Island Dressing

For cuts, grazes and covering sutures

This soft and flexible dressing pad is really comfortable and conforms to your skin as it is made from a soft medical grade fabric. It has a specially designed absorbent pad which has a cushioning effect to protect the wound. The pad absorbs excess wound fluid which is important for wound healing. The nature of this dressing makes it ideal for most areas of the body, particularly near joints or areas of high movement.

- Gentle to Skin
- Absorbs and cushions
- Helps protect against germs
- Hospital grade


BP011550BA Box/50 \$14.00



## Primapore Adhesive Wound Dressing

Primapore combines an absorbent pad with a soft and conformable fixative layer

- Soft breathable cover
- Highly absorbent (melolin) pad
- Low allergy adhesive



Smith+Nephew



7.2 x 5cm	SN66003634	Box/50	\$28.15	3+	\$25.90
8.3 x 6cm	SN7135	Box/50	\$43.30	3+	\$39.90
10 x 8cm	SN66000317	Box/20	\$46.10	3+	\$42.45
15 x 8cm	SN66000318	Box/20	\$65.90	3+	\$60.70
10 x 25cm	SN66000320	Box/20	\$83.85	3+	\$77.25

## Cutiplast Steril Conformable Fabric Island Dressing

Cutiplast Steril is a dressing for primary wound management. The wound pad is low-adherent so dressing changes do not disturb or interfere with the wound healing process.

Smith+Nephew

7.2 x 5cm	SN66001478	Box/100	\$40.05	2+	\$36.90
10 x 8cm	SN66001473	Box/50	\$67.30	2+	\$62.00
15 x 8cm	SN66001474	Box/50	\$96.25	2+	\$88.65



## Cutiplast Conformable Fabric Dressing Strips

Cutiplast conformable fabric island dressings are ideal for the treatment of minor wounds. Economical; dressing lengths can be cut to size.

Smith+Nephew

4cm x 5m	SN66001464	Roll	\$18.50	2+	\$17.05
6cm x 5m	SN66001466	Roll	\$22.85	2+	\$21.05
8cm x 5m	SN66001468	Roll	\$27.25	2+	\$25.10



## Compose® Island Dressing Adhesive Island Dressing

7.5 x 5cm	MGT30611	Box/50	\$10.05	3+	\$9.06
8.25 x 6cm	MGT30612	Box/50	\$16.35	3+	\$13.85
12 x 8.25cm	MGT30613	Box/20	\$11.05	3+	\$9.60
8.25 x 18cm	MGT30614	Box/20	\$15.05	3+	\$14.10



MULTIGATE



# bpositive

## Our film feels and acts like your skin

### Water & Bacteria Proof Dressing + Pad

Our matt film feels and acts like the skin.  
The thicker pad absorbs and protects the wound.  
Yet optimise patient comfort.

### Water & Bacterial Proof Roll

Same film technology.  
Helps dress difficult wounds.  
Great secondary or primary dressing.

We also have a **Water & Bacteria Proof Dressing** which is commonly used for minor wounds and surgery as a primary or secondary dressing.



Keep up to date  
by visiting our  
website.

[bpositivehealthcare.com](http://bpositivehealthcare.com)

### Follow us on social media

- Instagram: @bpositivehealth
- Facebook: bpositivehealthcare
- Twitter: bpositivehealth
- LinkedIn: bpositivehealthcare



## Waterproof Island Dressing

For cuts, stitches and grazes

This water and bacteria proof film dressing pad helps wounds to heal with a specialty designed absorbent pad which has a soft cushioning effect to protect the wound. The pad has an extra non-adherent layer which minimise pain when removing the dressing. In addition, the water and bacteria proof protective film helps reduce the risk of infection. As it is waterproof, you can shower with it.

- Waterproof
- Helps protect against germs
- Absorbs and cushions

4 x 5cm	BP020304BA	Box/50	\$42.05*	3+	\$37.75*
5 x 7cm	BP020305BA	Box/50	\$48.55*	3+	\$44.20*
8 x 10cm	BP020308BA	Box/25	\$44.20*	3+	\$39.90*
10 x 12cm	BP020312BA	Box/25	\$50.70*	3+	\$46.35*
10 x 15cm	BP020315BA	Box/20	\$48.55*	3+	\$44.20*



## Claripose® Island Dressing

Waterproof Island Dressing

6.5 x 5cm	MGT30651	Box/100	\$66.55*	3+	\$62.15*
9.5 x 8.5cm	MGT30652	Box/20	\$32.55*	3+	\$30.40*
10 x 12cm	MGT30653	Box/10	\$20.30*	3+	\$17.95*
15 x 8.5cm	MGT30654	Box/20	\$32.75*	3+	\$29.75*



MULTIGATE



## Mepilex® Border Post-Op

Mepilex® Border Post-Op is an all-in-one post-operative dressing designed to absorb and retain surgical wound fluid.

The design of Mepilex® Border Post-Op provides high flexibility and conformability over joints, such as knees or hips, promoting early patient mobilisation.

6 x 8cm	MOL496100	Box/10	\$73.75*	3+	\$66.70*
9 x 10cm	MOL496200	Box/10	\$111.85*	3+	\$101.15*
10 x 15cm	MOL496300	Box/10	\$164.40*	3+	\$148.65*
10 x 20cm	MOL496400	Box/10	\$236.70*	3+	\$215.70*



Designed with Safetac® technology to minimise pain and peri-wound skin damage at dressing removal.!

Safetac TECHNOLOGY



## Leukomed® T Plus

Leukomed T plus is a sterile dressing consisting of a film and a wound pad for small to medium-sized wounds. Well suited for surgical wounds and post-op care. Breathable and skin friendly yet secure, it can be worn for up to 7 days.

5 x 7.2cm	BSN7238200	Box/50	\$70.85*	3+	\$68.20*
8 x 10cm	BSN7238201	Box/50	\$201.65*	3+	\$185.60*
10 x 15cm	MOL496300	Box/10	\$164.40*	3+	\$148.65*



essity

# + When you value quality and patient satisfaction

The unique design of OPSITE POST-OP VISIBLE allows continual monitoring of the incision site without the need to disrupt the healing process 'providing both visibility and absorbency to promote longer wear times and fewer unnecessary dressing changes.'<sup>2</sup>

Getting you closer to zero delays in wound healing.

## OPSITE POST-OP VISIBLE


When you need continual monitoring.

- REACTIC<sup>®</sup> film
- High MVTR film
- Lattice hydrocellular foam
- Bacterial barrier
- Waterproof film

## OPSITE POST-OP

Offers a significantly longer wear time when compared with non-woven or passive dressings.<sup>1</sup>

- REACTIC<sup>®</sup> film
- High MVTR film
- Low-allergy adhesive
- Bacterial barrier
- Waterproof film



## Smith+Nephew

### OPSITE<sup>®</sup> POST-OP

Absorbent, Waterproof  
and Bacteria-proof Film Dressing

### OPSITE<sup>®</sup> POST-OP VISIBLE



Waterproof, Bacteria-proof Dressing  
with See-through Absorbent Pad



1. Barker A, et al. A ward-based in-practice evaluation of a post-operative dressing. Poster presented at EWMA conference, Geneva; 2010. 2. O'Brien G, et al. A prospective evaluation to assess performance of a new post-operative dressing. Poster presented at Infection Prevention Society Conference, Harrogate; 2009. Smith & Nephew Pty Ltd (Australia) [www.smith-nephew.com/australia](http://www.smith-nephew.com/australia). Smith & Nephew Ltd (New Zealand) [www.smith-nephew.com/new-zealand](http://www.smith-nephew.com/new-zealand). <sup>®</sup>Trademark of Smith+Nephew. All trademarks acknowledged. SN13505 V1 03/20. For detailed product information, including indications for use, contraindications, precautions and warnings, please consult the product's applicable Instructions for Use (IFU) prior to use.

## CUTIFILM Plus

CUTIFILM Plus is a waterproof dressing ideal for use on surgical wounds, lacerations and abrasions.



5 x 7.2cm	SN76248	Box/50	\$92.30*	3+	\$85.00*
8 x 10cm	SN76244	Box/25	\$99.35*	3+	\$91.50*
8 x 15cm	SN76245	Box/25	\$131.85*	3+	\$121.45*
10 x 12.5cm	SN76246	Box/25	\$136.50*	3+	\$125.75*
10 x 20cm	SN76247	Box/25	\$178.50*	3+	\$164.45*

- Sterile, bacteria proof and waterproof
- Hypo-allergenic adhesive
- Secure adhesion
- Low adherent pad
- Protects the skin
- Allows the skin to breathe

## OPSITE Post-Op

Post-operative dressing combines the advantages of OPSITE film with an absorbent dressing pad to provide a superior dressing for many types of wound.

- Waterproof, breathable and conformable
- Bacterial barrier helps to reduce the risk of secondary infection
- Absorbent, non-adherent pad minimises the number of dressing changes

Indications: Post-operative wounds, skin tears, minor cuts, abrasions and lacerations

6.5 x 5cm	SN66000708	Box/100	\$205.45*	3+	\$189.25*
9.5 x 8.5cm	SN66000709	Box/20	\$94.80*	3+	\$87.30*
12 x 10cm	SN66000710	Box/10	\$67.65*	3+	\$62.30*
15.5 x 8.5cm	SN66000712	Box/20	\$125.80*	3+	\$115.90*

OpSite® Post-Op is a revolutionary post-operative dressing combining a high Moisture Vapour Transmission Rate (MVTR) film and the cushioning and absorbency of a low adherent pad together with an easy application carrier system. The hydrophilic polyurethane film on OpSite Post-Op has REACTIC™ properties that 'switch on' in the presence of moisture allowing any excess to transpire into the atmosphere.



## OPSITE Post-Op Visible

Waterproof, bacteria-proof dressing with see-through absorbent pad


OpSite Post-Op Visible is a transparent film with a unique see-through absorbent pad. The highly breathable film is waterproof and bacteria-resistant, and the unique lattice structure foam pad enables regular wound assessment without the need to lift or remove the dressing.

- Absorbent lattice structured foam pad
- Highly breathable waterproof outer film
- Perforated low-adherent wound contact layer
- Low allergy adhesive
- Easy application flexible carrier

Indications: ideal for moderately exuding wounds including surgical wounds, lacerations and abrasions. The low allergy adhesive reduces the risk of irritation and keeps the dressing in place, even over awkward contours.

8 x 10cm	SN66800136	Box/20	\$133.00*	3+	\$122.55*
15 x 10cm	SN66800137	Box/20	\$193.35*	3+	\$179.95*
20 x 10cm	SN66800138	Box/20	\$238.60*	3+	\$219.80*

## Smith+Nephew



For complete listings of Smith & Nephew products, treatment options, education and event calendar visit <https://www.smith-nephew.com/australia/>



Smith+Nephew



### 3M™ Tegaderm™ + Pad Film Dressing with Non-Adherent Pad




A waterproof, absorbent dressing.

- All-in-one sterile dressing – reduces application time
- A bacterial and viral barrier – helps prevent external bacteria and contaminants from entering the site
- A waterproof barrier – patient may shower without removing dressing
- Conforms to body contours, providing comfort
- Non-adherent absorbent pad
- Breathable film allows skin to function normally with good exchange of moisture vapour and oxygen
- Provides moist wound environment to enhance healing and lessen pain
- Frame delivery system ensures precise placement without applying tension
- Hypoallergenic

Suggested Applications:

- Post-operative sites
- Superficial and partial thickness wounds
- Lacerations, abrasions and burns
- Light to moderately draining wounds
- Secondary wound dressing



5 x 7cm	3M3582	Box/50	\$62.00*	2+	\$57.50*
6 x 10cm	3M3584	Box/50	\$119.40*	2+	\$110.75*
9 x 10.5cm Oval	3M3587	Box/25	\$71.50*	2+	\$66.25*
9 x 10cm	3M3586	Box/25	\$67.50*	2+	\$62.60*
9 x 15cm	3M3589	Box/25	\$84.10*	2+	\$78.00*
9 x 25cm	3M3591	Box/25	\$112.85*	2+	\$103.95*

Picture frame delivery for easy application

### 3M™ Tegaderm™ Transparent Film Dressing



- Framed delivery system allows for one handed application and prevents dressing sticking to itself
- Transparent dressings allow easy monitoring of wounds and I.V. sites without the need to remove the dressing
- Secure adhesion, yet gentle to the skin

4.4 x 4.4cm	3M1622W	Box/100	\$63.55*	2+	\$58.55*
6 x 7cm	3M1624W	Box/100	\$65.10*	2+	\$60.35*
10 x 12cm	3M1626W	Box/50	\$77.90*	2+	\$72.25*
15 x 20cm	3M1628	Box/10	\$57.45*	2+	\$54.55*



### 3M™ Tegaderm™ IV Transparent Film Dressing with Border



- Notched design provide a better seal around the catheter
- Sterile tape strips are provided to assist with improved catheter securement
- Breathable film, allows oxygen in and moisture out
- Transparent dressings allow easy monitoring of I.V. sites without the need to remove the dressing
- Framed delivery system allows for one handed application and prevents dressing sticking to itself

7 x 8.5cm	3M1633	Box/100	\$129.90*	2+	\$120.45*
-----------	--------	---------	-----------	----	-----------



### 3M™ Tegaderm™ Absorbent Clear Acrylic Dressing



Transparent plus absorbent for optimal clinical performance and wear time.

A transparent wound dressing with fluid handling capacity that allows for wound monitoring without changing the dressing. For low to moderate exudate. It is waterproof, washable and breathable.

- Easy monitoring of the wound, exudate and periwound skin
- No upper wear time - eliminates unnecessary dressing changes
- Eliminates hydrocolloid-related odour
- No dressing breakdown in the wound

7.6 x 9.5cm	3M90800	Box/5	\$74.50*	3+	\$70.25*
11.1 x 12cm	3M90801	Box/5	\$165.35*	3+	\$157.10*
14.9 x 15.2cm	3M90802	Box/5	\$187.45*	3+	\$176.75*



Suggested Applications:

- Partial and full thickness dermal ulcers
- Pressure ulcers
- Abrasions, superficial wounds, donor sites and skin tears

## Claripose® Transparent Dressing

### Transparent Dressing

<b>6 x 7cm</b>	<b>MGT30632</b>			
Box/100	\$49.95*	3+	\$45.95*	
<b>10 x 12cm</b>	<b>MGT30633</b>			
Box/50	\$60.95*	3+	\$54.15*	



## Mepitel® Film

### Breathable, transparent film dressing with Safetac® Technology

Mepitel® film provides a flexible, transparent covering to protect the skin or wound from microbial contamination and fluid strikethrough. It maintains a moist environment even though the vapour permeability allows excess moisture to pass away from the skin. It provides instant tack adhesion that minimises the requirement of extra fixation.

Benefits:

- Minimises pain and trauma at dressing changes<sup>1</sup>
- Limits the risk for damage of newly formed tissue
- Transparent for easy inspection
- Conforms well to body contours, promoting patient comfort during wear
- Bacteria and viral barrier where microbes larger than 25nm
- Low potential for skin irritation and allergy
- Leaves no residues
- Fluid impermeable
- Showerproof



<b>6 x 7cm</b>	<b>MOL296170</b>	Box/10	\$25.45*	2+	\$23.00*
<b>10 x 12cm</b>	<b>MOL296270</b>	Box/10	\$39.80*	2+	\$36.25*

## WE ARE A BUSINESS FOR GOOD



Adrian from marketing passed his 8-year anniversary with Team Medical last year. He used his B1G1 credits to provide **45 hot meals for those experiencing homelessness in Hong Kong**. Tai Kok Tsui is one of the most impoverished areas of Hong Kong, where many homeless people are going without multiple nutritious meals every day.

This project goes directly towards providing hot meals for those who need them most.

Learn more about our giving impacts at [www.teammed.com.au/esg-giving-back](http://www.teammed.com.au/esg-giving-back)



## Meet our Team: Pieter

Our CEO Pieter decided the corporate world wasn't for him and joined his sister Marijke and son Matt at Team Medical in 2012. His son Tom soon followed in 2013 truly making Team Medical a family company.

With extensive industry experience, Pieter has worked tirelessly towards ensuring continued growth by being the best at what we do – whilst having fun along the way.

Pieter consistently embeds his family values into our company culture. *"At Team Medical, everyone is valued and respected. Seeing our people grow and thrive alongside our company and customers is very rewarding."*

Born in the Netherlands, Pieter is looking forward to his mother (86) coming to Australia later this year. He enjoys driving around in his Bobcat, playing tennis and watching soccer – no game is on too early or too late for Pieter, and always better watched with a crowd.

**What is your biggest pet peeve?** Cleaning dog poo, I have 3 dogs.

**What app do you use every day?** Fitbit

**What show have you been watching lately?** Yellowstone




## OPSITE Flexifix

Smith+Nephew

Transparent, Non-Sterile, Adhesive Film Roll

- Waterproof and bacteria proof
- Breathable
- Highly conformable: easy application to awkward areas of the body
- Ideal over primary dressings, fixation of tubes

<b>5cm x 10m</b>	<a href="#">SN66000040</a>	Roll	\$70.75*	3+	\$65.15*
<b>10cm x 10m</b>	<a href="#">SN66000041</a>	Roll	\$127.15*	3+	\$117.10*



The unique and novel roll format enables rapid and economical use of the OpSite film to create comfortable, waterproof fixation for dressings.

## OPSITE spray

Smith+Nephew

Spray Film Dressing

- Transparent and quick drying film
- Waterproof
- Easy and convenient to apply
- Reduces the risk of secondary infection

Indications: The film provides protection for dry minor surgical and surface wounds.

<a href="#">SN66004978</a>	110ml	\$46.40*
	3+	\$42.70*




## OPSITE Flexigrid

Smith+Nephew


Moisture Vapour Permeable Adhesive Film Dressing

- Transparent, adhesive film
- Unique wound measurement grid
- Very conformable
- Clinically proven to provide moist wound healing for superficial wounds

<b>6 x 7cm</b>	<a href="#">SN4628</a>	Box/100	\$93.45*	3+	\$86.10*
<b>10 x 12cm</b>	<a href="#">SN4629</a>	Box/10	\$36.25*	3+	\$33.95*
<b>10 x 12cm</b>	<a href="#">SN4630</a>	Box/50	\$102.45*	3+	\$94.40*
<b>12 x 25cm</b>	<a href="#">SN4632</a>	Box/20	\$128.85*	3+	\$119.50*



Indications: Minor wounds, cuts, grazes, superficial minor burns and skin tears. Ideal for use as a secondary dressing and to help prevent pressure injuries




## IV3000

Smith+Nephew

Moisture Responsive IV Catheter Dressing

IV3000\* is a moisture responsive, transparent film dressing, specifically designed to meet the needs of IV catheter fixation.

	<b>7 x 9cm</b>	Peripheral	<a href="#">SN4006</a>	Box/100	\$130.55*	3+	\$122.40*
	<b>6 x 7cm</b>	Peripheral/Ported	<a href="#">SN4007</a>	Box/100	\$120.90*	3+	\$111.40*
	<b>6 x 8.5cm</b>	Peripheral	<a href="#">SN4924</a>	Box/100	\$136.10*	3+	\$125.40*
	<b>10 x 14cm</b>	Central/Jugular	<a href="#">SN4925</a>	Box/10	\$48.55*	3+	\$45.50*



Orange-handle Incorporates tape for security of fixation



## OPSITE Flexifix Gentle

Smith+Nephew

Transparent silicone adhesive film roll

OPSITE Flexifix Gentle combines breathability, transparency and up to a seven day wear time to offer a versatile silicone film roll for patients with fragile or at-risk skin. It removes cleanly, without disrupting fragile skin layers or causing patients undue pain.

Indication: to retain primary dressings, fixation of tubing, protect skin from friction and external contamination, for prophylaxis against pressure sores, relieve pain in painful diabetic neuropathy, protect the skin around stoma and under leg bags.

<b>2.5cm x 5m</b>	<a href="#">SN66801195</a>	Roll	\$52.55*
<b>5cm x 5m</b>	<a href="#">SN66801196</a>	Roll	\$92.40*
<b>10cm x 5m</b>	<a href="#">SN66801197</a>	Roll	\$144.35*





## ActivHeal® AQUAFIBER™

ActivHeal® AquaFiber is a soft, conformable, highly absorbent dressing that converts to a soft clear gel when in contact with wound exudate.



- ✓ High absorption of exudate
- ✓ Excellent wet and dry strength
- ✓ Hidden integrity layer
- ✓ Easy to use
- ✓ Reduced lateral wicking
- ✓ Clear gelling
- ✓ Nylon mesh reinforced layer
- ✓ Promotes healing through a moist wound environment



**Indications for use:**

ActivHeal® AquaFiber is indicated for moderately to heavily exuding chronic and acute wounds. The dressing can also be used to control minor bleeding in superficial wounds.

- ✓ Post Operative Surgical Wounds
- ✓ Superficial and Partial Thickness Burns
- ✓ Pressure Ulcers
- ✓ Venous Leg Ulcers
- ✓ Arterial Leg Ulcers
- ✓ Diabetic Ulcers
- ✓ Cavity Wounds
- ✓ Graft Donor Sites
- ✓ Skin Graft Sites

Moderate > Heavy Exudate



 HIGHLY ABSORBENT	 HIGH WET STRENGTH	 LOW RISK OF MACERATION	 LOW LATERAL WICKING
---	--	---	--

5 x 5cm	ADM10009145	Box/10	\$33.60*	2+	\$31.00*
10 x 10cm	ADM10009146	Box/10	\$78.40*	2+	\$72.35*
15 x 15cm	ADM10009147	Box/5	\$88.40*	2+	\$80.15*
2 x 42cm Rope	ADM10009148	Box/5	\$47.15*	2+	\$42.75*

Reference: 1. Lansdown, A. (2006) Bioactive Dressings: old ideas, new technology. British Journal of Nursing. Tissue Viability supplement, Vol 16, no13



## ActivHeal® AQUAFIBER™ EXTRA

ActivHeal® AquaFiber™ Extra creates a moist wound healing environment for moderate to heavily exuding wounds.



- ✓ Highly absorbent<sup>1</sup>
- ✓ High wet tensile strength<sup>1</sup>
- ✓ Increased vertical wicking<sup>1</sup>
- ✓ High gelling
- ✓ Can remain in situ for up to 7 days
- ✓ Reduces the risk of maceration
- ✓ Aids autolytic debridement
- ✓ Maintains a moist wound environment





**Indications for use:**

ActivHeal® AquaFiber™ Extra is indicated for the management of the following moderately to heavily exuding acute and chronic wounds:

- ✓ Post Operative Surgical Wounds
- ✓ Superficial and Partial Thickness Burns
- ✓ Pressure Ulcers
- ✓ Venous Leg Ulcers
- ✓ Arterial Leg Ulcers
- ✓ Diabetic Ulcers
- ✓ Cavity Wounds
- ✓ Graft Donor Sites
- ✓ Skin Graft Sites

Moderate > Heavy Exudate




 HIGHLY ABSORBENT	 HIGH WET STRENGTH	 LOW RISK OF MACERATION	 INCREASED VERTICAL WICKING
---	--	---	---

5 x 5cm	ADM9023625	Box/10	\$48.80*	2+	\$45.05*
10 x 10cm	ADM9023632	Box/10	\$90.05*	2+	\$83.05*
15 x 15cm	ADM9023649	Box/5	\$94.25*	2+	\$85.50*
20 x 30cm	ADM9023656	Box/5	\$179.10*	2+	\$162.45*
2.5 x 30.5cm Rope	ADM9023663	Box/5	\$53.05*	2+	\$48.10*
2 x 46cm Rope	ADM9023670	Box/5	\$53.05*	2+	\$48.10*

Reference: 1. Advanced Medical Solutions Data on file LD024-16, LD037-16. ActivHeal and AquaFiber are registered Trademarks of Advanced Medical Solutions Limited.

# ActivHeal® Academy

Helping you provide the best possible care.



## SIMPLIFIED GUIDES

Bitesize modules from Wound Assessment to Dressing Selection designed to support your day-to-day practice.

## ACADEMY+

Advanced training modules endorsed by the Society of Tissue Viability to further your development.

## ASSESSMENTS

Interactive tests and reflective learning to measure your progress and support your Continuing Professional Development.



Stay informed and confident in your practice with the ActivHeal® Academy

**academy.activheal.com**



## ActivHeal® ALGINATE

ActivHeal® Alginate is a sodium calcium alginate dressing indicated for the treatment of moderate to heavily exuding wounds as a primary dressing.

- ✓ Excellent absorption of exudate
- ✓ Versatile
- ✓ Haemostatic properties - controls minor bleeding
- ✓ Reduces the risk of maceration
- ✓ Aids autolytic debridement
- ✓ Encourages granulation within the wound
- ✓ Promotes healing through a moist wound environment

### Indications for use:

ActivHeal® Alginate is indicated for moderately to heavily exuding wounds that are granulating or with areas of slough<sup>1</sup> including:

- ✓ Post Operative Surgical Wounds
- ✓ Superficial and Partial Thickness Burns
- ✓ Pressure Ulcers
- ✓ Leg and Foot Ulcers
- ✓ Cavity Wounds
- ✓ Graft Donor Sites
- ✓ Skin Graft Sites

Moderate > Heavy Exudate



5 x 5cm	ADM10007432	Box/10	\$35.60*	3+	\$32.85*
10 x 10cm	ADM10007431	Box/10	\$57.60*	3+	\$53.15*
10 x 20cm	ADM10007430	Box/5	\$63.30*	3+	\$58.35*
2 x 30cm Rope	ADM10007428	Box/5	\$46.40*	3+	\$42.80*



Reference: 1. Lansdown AB, Payne MJ. (1994) An evaluation of the local reaction and biodegradation of calcium sodium alginate (Kaltostat) following subcutaneous implantation in the rat. J R Coll Surg Edinb. Oct; 39 (5) : 284-8  
 ActivHeal is a registered Trademark of Advanced Medical Solutions Limited.

## WE ARE A BUSINESS FOR GOOD

Did you know that for every 10 ActivHeal products sold, we donate digital skills training to Indigenous Communities?

Skills range from music making and multi-media to essential computer-based skills like email and internet banking. These skills stay in the community and are then passed on to others, creating a ripple effect.

Thank you for choosing quality products and being part of creating **1748 training days** so far.

Learn more about our giving impacts at [www.teammed.com.au/esg-giving-back](http://www.teammed.com.au/esg-giving-back)





### Algisite M Calcium Alginate Dressing

Algisite M is a calcium-alginate dressing which forms a soft, integral gel when it comes into contact with wound exudate. Excess wound exudate is absorbed in this process. Algisite M utilises the proven benefits of moist wound healing, which results in faster healing of even difficult wounds.

Smith+Nephew



	2 x 30cm	SN66000522	Box/5	\$93.55*	5+	\$86.20*
	5 x 5cm	SN66000519	Box/10	\$72.50*	5+	\$66.80*
	10 x 10cm	SN66000520	Box/10	\$148.70*	5+	\$137.00*

### Kaltostat Alginate Calcium Sodium Dressing

ConvaTec



KALTOSTAT® Alginate Calcium Sodium Dressing is an alginate dressing indicated as a primary dressing for moderately to highly exuding chronic and acute wounds, and for wounds with minor bleeding. On contact with exudates the dressing gradually forms a firm, moist gel which is highly absorbent.

5 x 5cm	CON168210	Box/10	\$73.30*	2+	\$66.20*
7.5 x 12cm	CON168212	Box/10	\$144.05*	2+	\$129.00*
10 x 20cm	CON168214	Box/10	\$280.05*	2+	\$263.35*
2g rope	CON1681F1	Box/5	\$92.60*	2+	\$82.95*

KALTOSTAT®

### Durafiber Highly absorbent non woven, gelling fibre dressing

Smith+Nephew

DURAFIBER is a highly absorbent, non woven gelling fibre filler dressing. Composed of a unique blend of cellulose based fibres, DURAFIBER provides a versatile solution for the management of medium to heavily exuding chronic and acute full thickness, partial thickness or shallow granulating wounds.

- Forms a soft cohesive gel sheet on contact with wound fluid
- Locks in fluid and bacteria
- Controlled lateral wicking
- Lasting integrity and gelled strength
- Highly absorbent
- Minimal dressing shrinkage



4 x 10cm	SN66800546	Box/5	\$43.60*
4 x 20cm	SN66800547	Box/5	\$76.70*
5 x 5cm	SN66800559	Box/10	\$66.35*
10 x 10cm	SN66800560	Box/10	\$176.80*
15 x 15cm	SN66800561	Box/5	\$185.70*
2 x 45cm	SN66800563	Box/5	\$108.75*

Indications: As an absorbent dressing, gelling dressing for the management of chronic and acute, full thickness, partial thickness, or shallow granulating exuding wounds

### Aquacel Extra Hydrofiber Dressing

ConvaTec

The unique gelling action of Hydrofiber® fibres convert to a soft gel when in contact with wound exudate. It absorbs and locks in the exudate, preventing spreading. Aquacel Extra provides added strength and absorption compared to AQUACEL, for when greater efficiencies are desired.

5 x 5cm	CON420671	Box/10	\$73.55*	2+	\$67.15*
10 x 10cm	CON420672	Box/10	\$173.60*	2+	\$153.40*
2cm x 45cm Rope	CON403770	Box/5	\$111.35*	2+	\$98.40*



## Solugel

### Wound Care Gel

Provides a soothing, moist healing environment for the effective treatment of minor burns, scalds and sunburn and relief of the pain and itchiness of chickenpox and shingles rash. Solugel is also excellent for the treatment of chronic ulcerative wounds, cuts and grazes.


50g	JJ10340				
Tube		\$11.95*			
12+		\$10.95*			



## Solosite gel

### Hydrogel Wound Dressing

- Preserved for multi-dose dispensing
- Effective debrider and deslougher
- Promotes moist wound healing
- Cost-effective

	20g	SN36100614	Tube	\$11.10*	3+	\$10.25*
	50g	SN36361354	Tube	\$13.20*	3+	\$12.15*




Smith+Nephew

Indications: The treatment of both shallow and deep wounds, particularly chronic lesions such as pressure injuries and leg ulcers which may be experiencing delayed healing due to the presence of necrotic and sloughy tissue. Also ideal for maintaining a moist wound healing environment in granulating wounds.

## Intrasite gel

### Premium Desloughing and Debriding Hydrogel

	8g	SN7308	Tube	\$15.00*	3+	\$13.80*
	15g	SN7311	Tube	\$17.15*	3+	\$15.80*
	25g	SN7313	Tube	\$24.30*	3+	\$22.35*



Smith+Nephew


INTRASITE Gel is an amorphous hydrogel which promotes rapid but gentle debridement of necrotic tissue, whilst being able to loosen and absorb slough and exudate.

It can also be used to provide the optimum moist wound healing environment during the later stages of wound healing.

## HydroTac® Transparent

### Transparent hydrogel sheet

HydroTac Transparent Comfort offers an additional acrylic adhesive border, providing an all-in-one dressing for water-proof protection of the wound.



## HydroTac Transparent

HydroTac Transparent will require a retention tape or bandage. Best for fragile skin.

5 x 7.5cm	HTM685905	Box/10	\$105.65*	2+	\$94.60*
10 x 10cm	HTM685906	Box/10	\$96.55*	2+	\$86.45*



## HydroTac®

### Unique foam with hydrogel layer

HydroTac is a foam dressing with a unique hydrogel wound contact layer. It effectively absorbs moderate exudate while at the same time providing the optimal environment for wound healing. HydroTac is gentle on the skin and can remain in place for up to 7 days.

- Accelerates wound healing
- Hydro-responsive dressing
- Optimal Absorption



### HydroTac Foam

HydroTac will require a retention tape or bandage. Best for fragile skin.

<b>6cm round</b>	<b>HTM685849</b>	Box/10	<b>\$86.40*</b>	2+	<b>\$80.70*</b>
<b>10 x 10cm</b>	<b>HTM685832</b>	Box/10	<b>\$154.55*</b>	2+	<b>\$145.15*</b>
<b>10 x 20cm</b>	<b>HTM685841</b>	Box/10	<b>\$236.70*</b>	2+	<b>\$221.90*</b>

### HydroTac Comfort

HydroTac Comfort features its own adhesive film border.

<b>6.5 x 10cm</b>	<b>HTM685811</b>	Box/10	<b>\$77.00*</b>	2+	<b>\$72.20*</b>
<b>8 x 8cm</b>	<b>HTM685810</b>	Box/10	<b>\$90.30*</b>	2+	<b>\$86.20*</b>
<b>12.5 x 12.5cm</b>	<b>HTM685815</b>	Box/10	<b>\$154.55*</b>	2+	<b>\$145.15*</b>



SALES - NSW

## Meet Our Team: Dianne

Dianne and her brother Harry were Team's first Sydney sales team 13 years ago. Dianne has over 25 years of experience in the medical supplies industry. Dianne loves to look after her customers and is always willing to go above and beyond.

Outside work, Dianne loves spending time in her garden, knitting many different things and singing her heart out at Karaoke. Dianne now also loves spending time with her grandson.

Dianne's role model is her dad, who fled the North of Vietnam to the South when he was only 24 and worked hard to provide for their family (wife and 7 children including Harry and herself).

**What sports do you like playing?** Badminton  
**Favourite holiday destination?** Hawaii



SALES - NSW

## Meet Our Team: Harry


Harry has been looking after his customers for Team Medical for over 13 years and is never shy to offer his expert opinion. He aims to help his customers have a smoother-running practice and is focused on always offering better value.

Harry recently went back to Vietnam to spend time with his family.

Harry enjoys staying active by playing tennis; he also enjoys watching swimming and soccer.

**What is your favourite animal?** Cat

**Which series do you watch?** SVU (Law & Order), Bridgerton and Mom (American comedy)



## Duoderm CGF

DuoDERM® CGF™ Dressing is a hydrocolloid, moisture-retentive wound dressing used for partial and full-thickness wounds with exudate. It incorporates a unique ConvaTec hydrocolloid formulation that distinguishes it from other hydrocolloid dressings.

Indicated for:

Chronic wounds - pressure ulcers (Stage I-IV), leg ulcers.

Acute wounds - traumatic wounds (minor abrasions, lacerations) partial thickness burns and donor sites

<b>10 x 10cm</b>	<b>CONH7660</b>	Box/5	<b>\$63.80*</b>	2+	<b>\$57.15*</b>
<b>15 x 15cm</b>	<b>CONH7661</b>	Box/5	<b>\$132.85*</b>	2+	<b>\$122.15*</b>



## Duoderm Extra Thin

A hydrocolloid dressing indicated for the management of lightly exuding wounds.

It combines a unique ConvaTec hydrocolloid formulation that distinguishes it from other hydrocolloid dressings and a vapour-permeable outer film to provide an occlusive moist environment.

<b>7.5 x 7.5cm</b>	<b>CONH7901</b>	Box/20	<b>\$135.20*</b>	2+	<b>\$119.50*</b>
<b>10 x 10cm</b>	<b>CONH7955</b>	Box/10	<b>\$124.85*</b>	2+	<b>\$111.80*</b>
<b>15 x 15cm</b>	<b>CONH7957</b>	Box/10	<b>\$142.50*</b>	2+	<b>\$131.00*</b>



## Meet our Team: Emma

Emma can be found in our purchasing team organising the seamless flow of stock from suppliers. She is a Lego fan, with her other hobbies including photography, hiking (with her dog Koby), go karting and following the F1.

What is your pet peeve? Loud chewers

What app do you use daily? Spotify





## Comfeel Plus Transparent

Comfeel Plus Transparent is a thin and flexible hydrocolloid dressing

It can be used for treatment of no-to-low exuding chronic wounds and superficial acute wounds in the final stages of wound healing. Comfeel Plus Transparent may be used for superficial burns, superficial partial-thickness burns, donor sites, post-operative wounds and skin abrasions. The dressing can stay on the wound for up to 7 days depending on the specific wound characteristics.

Comfeel Plus Transparent hydrocolloid dressing is the thin and flexible solution for sealing and healing no-to-low exuding wounds.



<b>5 x 7cm</b>	<b>COL3530</b>	Box/10	<b>\$81.10*</b>		
<b>5 x 25cm</b>	<b>COL3548</b>	Box/5	<b>\$82.90*</b>		
<b>10 x 10cm</b>	<b>COL3533</b>	Box/10	<b>\$140.15*</b>		
<b>9 x 14cm</b>	<b>COL3536</b>	Box/10	<b>\$179.10*</b>		

### HydroTac Comfort

HydroTac Transparent Comfort offers an additional acrylic adhesive border

<b>6.5 x 10cm</b>	<b>HTM685811</b>	Box/10	<b>\$77.00*</b>	2+	<b>\$72.20*</b>
<b>8 x 8cm</b>	<b>HTM685810</b>	Box/10	<b>\$90.30*</b>	2+	<b>\$86.20*</b>
<b>12.5 x 12.5cm</b>	<b>HTM685815</b>	Box/10	<b>\$154.55*</b>	2+	<b>\$145.15*</b>

