

## HEINE BETA® Series

The top-of-the-line full-size "BETA" Series offers the best performance, function and durability.

The BETA 400 LED F.O. Otoscope is the only otoscope with 4.2 magnification and LED Fibre Optic illumination.

The BETA 200 LED Ophthalmoscope is the only ophthalmoscope with an aspherical optical system and LED illumination.

The handles offer next-generation energy management and high-quality metal construction.



LED NOW IN HEINE QUALITY.

Heine's exclusive LED<sup>HQ</sup> illumination offers virtually unlimited hours of working life (up to 100,000 hours). No need to ever replace bulbs. Maintenance-free for a lifetime.



### HEINE BETA® 400 LED F.O. Otoscope

4.2x Magnification with LED<sup>HQ</sup> illumination

- 4.2x Magnification
- LED in HQ with continuous brightness control
- Patented, unique stepless dimming
- All-metal housing



For more information download the product brochure from our website



### HEINE BETA® 200 LED Ophthalmoscope

Superior aspherical optics with LED<sup>HQ</sup> illumination

- Aspherical optical system
- LED in HQ with continuous brightness control
- Exclusive dustproof design
- Aluminium frame system



For more information download the product brochure from our website



#### BETA 400 LED F.O. Otoscope

BETA 4 3.5V USB	BETA 400 LED F.O. Otoscope	BETA 4 USB rechargeable handle, USB cord and plug-in power supply, 1 set of reusable tips, 10 disposable tips	<a href="#">HEB14328388</a>	RRP \$1523.00	<b>\$1450.00</b>
-----------------	----------------------------	---	-----------------------------	---------------	------------------

#### BETA 400 LED F.O. Diagnostic sets

BETA 4 3.5V USB	BETA 400 LED F.O. Otoscope BETA 200 LED Ophthalmoscope	BETA 4 USB rechargeable handle, USB cord and plug-in power supply, 1 set of reusable tips, 10 disposable tips	<a href="#">HEB15328388</a>	RRP \$2543.00	<b>\$2385.00</b>
BETA 4 3.5V NT	BETA 400 LED F.O. Otoscope BETA 200 LED Ophthalmoscope	BETA 4 NT rechargeable handle, NT 4 Charger, 1 set of reusable tips, 10 disposable tips	<a href="#">HEB15324420</a>	RRP \$3028.00	<b>\$2845.00</b>

Visit our website to view the entire range, search "Heine BETA"

# Enabling health care professionals to make the earliest and most accurate diagnosis.

HEINE Diagnostics: Better illumination.  
Higher colour fidelity. More magnification.





## HEINE K Series

The HEINE-Standard full-size instrument for professionals, far beyond usual market standards. The line is characterised by unbreakable polycarbonate construction and full-feature performance to ensure the highest level of diagnostic accuracy.



If a HEINE instrument is used properly, it will still work like it did on day one even after many years of use. There's a very good reason why HEINE offers a five-year materials and workmanship guarantee on all of the products in this catalogue.

## HEINE K180® F.O. Otoscope

The standard, all around Fiber Optic Otoscope

- Fiber Optic Illumination
- Polycarbonate construction
- XHL Xenon Halogen Technology
- Acrylic viewing window with 3 x magnification
- Swivelling viewing window built into instrument



## HEINE K180® Ophthalmoscope

A full-featured instrument at an economical price

- Precision HEINE optics. Brilliant retinal image
- Bright, white light. Clear retinal image
- Dust-resistant housing
- Acrylic viewing window with 3 x magnification
- Swivelling viewing window built into instrument



For the ultimate convenience, upgrade your set with an additional handle

See page 30 for handles and accessories

K180 F.O. Otoscope						
BETA 2.5V		K180 F.O. Otoscope	BETA 2.5V handle (requires 2 x C size batteries) 1 set of reusable tips, 10 disposable tips	<a href="#">HEB18110118</a>	RRP \$579.00	\$550.00
K180 Ophthalmoscope						
BETA 2.5V		K180 Ophthalmoscope	BETA 2.5V handle (requires 2 x C size batteries)	<a href="#">HEC18210118</a>	RRP \$768.00	\$720.00
K180 Diagnostic sets						
BETA 2.5V		K180 Ophthalmoscope K180 F.O. Otoscope	BETA 2.5V handle (requires 2 x C size batteries) 1 set of reusable tips, 10 disposable tips	<a href="#">HEA27910118B</a>	RRP \$1238.00	\$1165.00
BETA 4 3.5V USB		K180 Ophthalmoscope K180 F.O. Otoscope	BETA 4 USB rechargeable handle, USB cord and plug-in power supply, 1 set of reusable tips, 10 disposable tips	<a href="#">HEA27927388</a>	RRP \$1741.00	\$1635.00
BETA 4 3.5V NT		K180 Ophthalmoscope K180 F.O. Otoscope	BETA 4 NT rechargeable handle, NT 4 Charger, 1 set of reusable tips, 10 disposable tips	<a href="#">HEA27923420</a>	RRP \$2226.00	\$2090.00

# HEINE EN200 Diagnostic Centre



**Diagnostic instruments located at your fingertips in a space-saving design**  
**Everything in its right place, ready to use. No searching for instruments. Always charged.**

Due to cutting-edge technology in the HEINE EN200 the handles switch on and off automatically when taken off or returned to the unit.

Set your brightness level once and it stays there every time you use your instrument.

When not in use and when not returned to the unit, the illumination switches off after 15 minutes.



For more information download the product brochure from our website

EN200 Diagnostic Centre Kits				
<b>BETA KIT</b>	BETA 400 LED F.O. Otoloscope BETA 200 LED Ophthalmoscope EN200 Wall Transformer with 2 Handles and Coiled Tubing	<a href="#">HEA09512208</a>	RRP \$2788.00	<b>\$2615.00</b>
<b>K180 KIT</b>	K180 F.O. Otoloscope K180 Ophthalmoscope EN200 Wall Transformer with 2 Handles and Coiled Tubing	<a href="#">HEA09512220</a>	RRP \$1963.00	<b>\$1845.00</b>



## Meet our Team: Gillian



With over 30 years of experience in Medical Sales, Gillian knows almost every product inside and out. She is the perfect matchmaker when it comes to providing personalised product solutions to meet our customers' needs.

Gillian is an animal lover and a part-time goat breeder; it's no surprise her next holiday destination involves heading to her favourite farm. Gillian can also be found in her downtime doing charity collections for Baby Give Back.

**What's a favourite app you use?**  
Wake up app

**Any hobbies?** Photography  
**Coffee or Tea?** Only coffee





## HEINE mini 3000® Series

Maximum quality with minimum dimensions.

The HEINE mini 3000 series is the successful combination of the high quality features of full-size HEINE instruments and the convenience of pocket sized design.

- **Handy and ergonomic** - compact and ergonomic design guarantees a comfortable operation for the most reliable diagnosis in any examination situation.
- **A quality product with a modern design** - special plastics provide maximum strength and durability in a lightweight and modern design
- **Quality optics and illumination** - high quality optics produce distortion-free and brilliant images of the whole area under examination. The illumination system renders structures and colours accurately and naturally.



### Two illumination options. One compact series of instruments without equal.

**HEINE LED<sup>HQ</sup>** - Heine's exclusive LED<sup>HQ</sup> illumination offers virtually unlimited hours of working life. No need to ever replace bulbs. Maintenance-free for a lifetime. Twice as bright as XHL instruments


**HEINE XHL Technology** - Heine's XHL Xenon-Halogen Technology offers the perfect colour rendering that health care professionals have examined the eye with and depended on for decades.



For more information on the HEINE Mini 3000 Series download the product brochure from our website

### HEINE mini 3000® Otoscope

- Xenon Halogen Technology
- Direct illumination
- 3x magnification
- Swivelling viewing window




### HEINE mini 3000® LED F.O. Otoscope

- Maintenance-free pocket instrument
- Bright LED illumination
- Fiber Optic (F.O.) Illumination
- Twice as bright as Xenon Halogen instruments
- Operation time up to 10 hours
- Exclusive battery performance indicator
- 3x magnification
- Swivelling viewing window



### HEINE mini 3000® F.O. Otoscope

- Xenon Halogen Technology
- Fiber Optic (F.O.) Illumination
- 3x magnification
- Swivelling viewing window



**Maintenance-free pocket instrument**  
**Bright LED illumination**  
**Up to 10 hours operation time**

mini 3000 Otoscope						
XHL 2.5V		mini 3000 XHL Otoscope	mini 3000 battery handle, 10 disposable tips, 2 x AA batteries	<a href="#">HED00170210</a>	RRP \$233.00	<b>\$230.00</b>
XHL F.O. 2.5V		mini 3000 XHL F.O. Otoscope	mini 3000 battery handle, 10 disposable tips, 2 x AA batteries	<a href="#">HED00170110</a>	RRP \$357.00	<b>\$350.00</b>
LED F.O. 2.5V		mini 3000 LED F.O. Otoscope	mini 3000 battery handle, 10 disposable tips, 2 x AA batteries	<a href="#">HED00870110</a>	RRP \$470.00	<b>\$460.00</b>

# QUALITY is in our DNA and in EVERY instrument.



## HEINE mini 3000® Ophthalmoscope

- Xenon Halogen Technology
- Crisp, clear images of the whole area under examination and minimised reflections
- 5 different apertures
- Range of lenses: -20 D to +20 D



## HEINE mini 3000® LED Ophthalmoscope

- Maintenance-free, no need to ever exchange the LED
- Optimal colour rendering
- High-quality optics provides sharp, clear images
- Up to 10 hours operation time
- 5 different apertures
- Range of lenses: -20 D to +20 D

**LED<sup>HO</sup>**  
LED NOW IN HEINE QUALITY

**With LED<sup>HO</sup> illumination; never replace your bulbs again**

mini 3000 Ophthalmoscope						
XHL 2.5V		mini 3000 XHL Ophthalmoscope	mini 3000 battery handle, 2 x AA batteries	<a href="#">HED00171120</a>	RRP \$431.00	<b>\$420.00</b>
LED 2.5V		mini 3000 LED Ophthalmoscope	mini 3000 battery handle, 2 x AA batteries	<a href="#">HED00871120</a>	RRP \$569.00	<b>\$555.00</b>

## HEINE mini 3000® Diagnostic Sets

Modern pocket otoscopes with bright illumination and accurate colour rendering. Unique design, compact and powerful.

Great 5-year warranty

mini 3000 Diagnostic Sets						
XHL 2.5V		mini 3000 XHL Oscope mini 3000 XHL Ophthalmoscope	2 mini 3000 battery handles, 5 reusable tips, 10 disposable tips	<a href="#">HED87321021</a>	RRP \$765.00	<b>\$750.00</b>
XHL F.O. 2.5V		mini 3000 XHL F.O Oscope mini 3000 XHL Ophthalmoscope	2 mini 3000 battery handles, 5 reusable tips, 10 disposable tips	<a href="#">HED87311021</a>	RRP \$889.00	<b>\$870.00</b>
LED F.O. 2.5V		mini 3000 LED F.O Oscope mini 3000 LED Ophthalmoscope	2 mini 3000 battery handles, 5 reusable tips, 10 disposable tips	<a href="#">HED88611021</a>	RRP \$1141.00	<b>\$1099.00</b>
LED F.O. 2.5V		mini 3000 LED F.O Oscope mini 3000 LED Ophthalmoscope	2 mini 3000 rechargeable battery handles, 10 disposable tips	<a href="#">HED86011022</a>	RRP \$1617.00	<b>\$1580.00</b>

# Original handles. Original power sources. Original bulbs. 100% HEINE.

5 Year Guarantee - built to last.



## HEINE NT4 Table charger

- Two independent charging compartments
- Automatic charging
- Handles cannot be overcharged
- Charge status display





### NT4 Table Charger

3.5V NT4	NT4 charger, with adaptor. For BETA 4 NT rechargeable handles, Standard F.O. 4 NT and F.O. 4 SHORT NT laryngoscope.	HEX00299494	RRP \$683.00	<b>\$640.00</b>
----------	---	-------------	--------------	-----------------

## HEINE BETA® 4 Rechargeable Handles

### Next generation energy management

- High-quality metal construction for strength and durability
- Exclusive automatic connector
- Li-ion technology
- Charge Status Indicator




[1]



[2]



[4]



BETA 4 NT rechargeable handle for NT4 table charger				
Rechargeable handle complete with Li-ion rechargeable battery and BETA 4 NT bottom insert	HEX00799396	RRP \$525.00	<b>\$490.00</b>	

BETA 4 USB rechargeable handle				
Rechargeable handle complete with Li-ion rechargeable battery, BETA 4 USB bottom insert, USB cord with medical approved plug-in power supply and handle rest	HEX00799388	RRP \$634.25	<b>\$615.00</b>	
USB Cord with E4-USB plug-in power supply [1] + [2]	HEX00099303		<b>\$105.00</b>	
E4-USB Medical approved plug-in power supply [2]	HEX00099305		<b>\$83.80</b>	

BETA 4 Spare Parts (Suits BETA 4 NT and BETA 4 USB rechargeable handles)				
BETA 4 NT bottom insert [3]	HEX00299394		<b>\$188.05</b>	
Li-ion rechargeable battery for BETA 4 rechargeable handle [4]	HEX00799383		<b>\$210.00</b>	




For more information on BETA® 4 Rechargeable Handles download the product brochure from our website



## HEINE BETA Battery Handle

BETA 2.5V handle  
Requires 2xC sized batteries



HEX00199118  
**\$190.00**

## HEINE Rechargeable Batteries


NiMH Rechargeable Battery  
for BETA NT300 rechargeable handle  
HEX00299382

**\$190.00**



Li-ion rechargeable battery  
for BETA 4 rechargeable handle  
HEX00799383

**\$210.00**



## HEINE Insufflation Bulb

- For pneumatic testing
- Suits Heine BETA400, BETA200, K180, BETA100 and K100 Otoscopes



HEB00011240  
**\$25.90**

## HEINE mini NT Charger Set for 2 mini 3000 Handles

With this set, existing mini 3000 handles can be updated to a 3.5V rechargeable system.



mini NT Charger Set			
Mini NT charger, 2 mini NT rechargeable batteries and 2 bottom inserts for the mini 3000 rechargeable handles, 10 disposable tips	HEX00199485	RRP \$626.00	<b>\$615.00</b>

## HEINE AllSpec® Disposable tips

- Secure fit
- Best light transmission and hygiene
- No sharp edges
- Grey coloured tips

2.5mm Child  
4mm Adult



### AllSpec Tips Packets of 50

2.5mm [HEB00011128](#)  
Pkt/50 \$12.25 20+ **\$11.05**



4mm [HEB00011127](#)  
Pkt/50 \$12.25 20+ **\$11.05**



### AllSpec Tips Desktop Box of 250

2.5mm [HEB00011151](#)  
4mm [HEB00011150](#)

Box/250 \$108.20  
3+ **\$87.65**



## HEINE AllSpec Tip Dispenser

- The right size, visible and instantly available
- Wall mount storage for up to 200 AllSpec tips

[HEB00011149](#)  
**\$78.85**



## HEINE UniSpec® Disposable Ear Specula (black)

For BETA100 & K100 Operating Otoscope

2.5mm [HEB00011242](#)  
Pkt/50 \$15.50 20+ **\$11.75**

4mm [HEB00011241](#)  
Pkt/50 \$15.50 20+ **\$11.75**



## HEINE XHL® Xenon Halogen Bulbs

- Exceptional performance
- Optimum light output
- Extreme long life

Otosopes		2.5V	3.5V
BETA 200 #077 - #078		<a href="#">HEX00188077</a>	<a href="#">HEX00288078</a>
K180 #077 - #078		<a href="#">HEX00188077</a>	<a href="#">HEX00288078</a>
Mini 3000 #110		<a href="#">HEX00188110</a>	
Mini 3000 F.O. #105		<a href="#">HEX00188105</a>	
Mini 2000 #056		<a href="#">HEX00188056</a>	
Mini 2000 F.O. #037		<a href="#">HEX00188037</a>	



Ophthalmoscopes		2.5V	3.5V
BETA 200 #069 - #070		<a href="#">HEX00188069</a>	<a href="#">HEX00288070</a>
K180 #084 - #086		<a href="#">HEX00188084</a>	<a href="#">HEX00288086</a>
Mini 3000 #106		<a href="#">HEX00188106</a>	
Mini 2000 #042		<a href="#">HEX00188042</a>	



SALES -NSW

## Meet our Team: Sharon

From working in the Team Medical warehouse to customer service to sales, Sharon has a great understanding of our products and processes and, most importantly, her customers' requirements. Sharon has an eye for detail and considers our customer relationships more like friendships (in the same way we think of our team as a family). After 11 years with Team Medical, Sharon still takes great joy in servicing her customers in the Hills Districts, the Northern Beaches, and the Central West.

In her downtime, Sharon enjoys cycling, hiking and painting by numbers.

What is your dream holiday destination? New York

What is your favourite animal? Alpacas





## Welch Allyn

### MACROVIEW OTOSCOPIES

- 3x larger view vs. standard otoscopes
- Long-lasting SureColor™ LED illumination
- Focus-Free Design
- Built tough with DuraShock™ technology
- Integrated locking collar

#### MacroView Basic LED Otoscope

- 3-year warranty



Otoscope Head Only  
**WA2382** RRP \$518.70  
**\$480.00**

#### MacroView Plus LED Otoscope

- Compatible with iExaminer® to digitally capture and share images
- Pair the MacroView otoscope with Welch Allyn LumiView™ Clear Single-Use Ear Specula for brighter views.
- 5-year warranty

Otoscope Head Only  
**WA2382** RRP \$518.70  
**\$480.00**



See bigger and better than ever before with PanOptic Plus Ophthalmoscopes and MacroView™ Plus Otoscope.

A new standard in physical assessment.

SCAN TO WATCH

## Welch Allyn

### PANOPTIC OPHTHALMOSCOPES

- 20x larger view vs. standard ophthalmoscopes
- Long-lasting LED illumination
- -20 to +25 diopters
- Built tough with DuraShock™ technology
- Integrated locking collar

#### PanOptic Basic LED Ophthalmoscope

- 3-year warranty



Ophthalmoscope Head Only  
**WA1182** RRP \$1094.40  
**\$1000.85**

#### PanOptic Plus LED Ophthalmoscope

- Compatible with iExaminer® to digitally capture and share images
- QuickEye alignment feature; Two indicator lights help you direct the patient's gaze for better optic disc alignment
- 5-year warranty

Ophthalmoscope Head Only  
**WA1183** RRP \$1362.30  
**\$1245.00**



**Baxter**

## Welch Allyn

GREEN SERIES 777

INTEGRATED WALL SYSTEMS



The diagnostic tools you need to maximise productivity as you capture patient in-take and vital signs — within arm's reach.

With the Green Series™ 777 Integrated Wall System, all your essential medical diagnostic tools are within reach.



### Green Series 777 Wall Transformer with PanOptic Plus LED Ophthalmoscope and MacroView Plus LED Otoscope

Includes:

- PanOptic Plus LED Ophthalmoscope
- MacroView™ Plus LED Otoscope
- iExaminer Smart Bracket
- Green Series™ 777 Wall Transformer



The Plus devices let you move from an optical exam to digital in a few clicks, expanding your capabilities. Through the free iExaminer Pro App, save and share eye and ENT images and videos for tracking and trending patient exams and consultations with specialists.

**\$2610**

RRP \$2508.00

WA777PM3XXX6

### Green Series 777 Wall Transformer with PanOptic Basic LED Ophthalmoscope and MacroView Basic LED Otoscope

Includes:

- PanOptic Plus LED Ophthalmoscope
- MacroView™ Plus LED Otoscope
- Li-Ion Basic Handle x2
- Universal Desk Charger



**\$2320**

RRP \$2257.20

WA777PM2XXX6

Choose the configuration that best suits your practice. All of your essential exam tools can be included.

Visit us online to view all available options.  
[www.teammed.com.au](http://www.teammed.com.au)



**Baxter**

## Welch Allyn

### Pocket LED Diagnostic Set

- Pocket LED Ophthalmoscope and LED F.O. Oscope
- Features SureColor™ LED illumination
- Lightest and most durable Welch Allyn Pocket Set
- Uses 2x AA batteries (not included)
- Comes in a soft, zippered case
- 3-year warranty

Black WA92871BLK

Plum WA92871PUR

Blueberry WA92871BLU **\$365**



## Welch Allyn

### Pocket LED Otoscope

- SureColor™ LED illumination
- Lightweight and durable
- Fiber Optic illumination
- Uses 2x AA batteries (not included)
- 3-year warranty

Black WA22870

Plum WA22870PUR

Blueberry WA22870BLU

**\$170**





**Change your view.**  
 Take a closer look at the next-generation diagnostic tools from Welch Allyn.



**PanOptic MacroView Plus**  
 Diagnostic Set

Effortlessly Switch Between Optical and Digital Examinations with the Plus Diagnostic Set.

Includes:

- PanOptic Plus LED Ophthalmoscope
- MacroView™ Plus LED Otoscope
- Li-Ion Plus USB-C Handle with Charging Cable
- iExaminer Smart Bracket
- Hard Case



**\$2535**

RRP \$2451

WA71PM3LXE



Enhance eye exams with the PanOptic Plus, featuring QuickEye alignment technology.

Switch between optical and digital



**PanOptic MacroView Basic**  
 Diagnostic Set

Includes:

- PanOptic Basic LED Ophthalmoscope
- MacroView™ Basic LED Otoscope
- Li-Ion Basic Handle with Charging Pod
- Hard Case



**\$1780**

RRP \$1767

WA71PM2LXU

**Coaxial MacroView Basic**  
 Portable Diagnostic Set

Includes:

- Coaxial LED Ophthalmoscope
- MacroView™ Plus LED Otoscope
- Li-Ion Basic Handle with Charging Pod
- Hard Case



**\$1350**

RRP \$1276.80

WA71SM2LXU

**PanOptic MacroView Plus**  
 Desk Diagnostic Set

Includes:

- PanOptic Plus LED Ophthalmoscope
- MacroView™ Plus LED Otoscope
- Li-Ion Plus USB-C Handle x2
- Universal Desk Charger
- iExaminer Smart Bracket



**\$2790**

RRP \$2679.00

WA71PM3LDE6



**PanOptic MacroView Basic**  
 Desk Diagnostic Set

Includes:

- PanOptic Plus LED Ophthalmoscope
- MacroView™ Plus LED Otoscope
- Li-Ion Basic Handle x2
- Universal Desk Charger



**\$2370**

RRP \$2223.00

WA71PM2LDX6


## Welch Allyn

ACCESSORIES

### LumiView Clear Single-Use Ear Specula

The new LumiView clear ear specula are one of a kind, with 4x more brightness through an adult specula, compared to a traditional black specula. Plus, it's backwards compatible to work universally on all Welch Allyn otoscopes.

2.75mm [WA52432CLR1](#) Pkt/34 \$7.00  
4.25mm [WA52434CLR1](#) 5+ \$6.65



### KleenSpec Single-Use Ear Specula

2.75mm [WA52432](#)  
4.25mm [WA52434](#)

Pkt/34 \$7.00  
5+ \$6.65



### KleenSpec Specula Dispenser


Suits 524-Series Ear Specula

[WA52400](#)  
\$90.00



### Welch Allyn 3.5V HPX Bulbs


**03100-U** [WA03100](#) 3.5 V Halogen Lamp for Otoscope (#20000)  
**06500-U** [WA06500](#) 3.5 V Halogen Lamp for MacroView™ Otosopes (#238XX)  
**03000-U** [WA03000](#) 3.5 V Halogen Lamp for Standard Ophthalmoscopes (#11710) and EpiScope (#47300)  
**04900-U** [WA04900](#) 3.5 V Halogen Lamp for Coaxial Ophthalmoscopes (#11720)  
**03800-U** [WA03800](#) 3.5 V Halogen Lamp for PanOptic Ophthalmoscopes (#11800, #11810, #11820)



### Welch Allyn SureColor LED Bulbs

- SureColor™ LED technology enabling reduction of light intensity to alleviate glare or enhance patient comfort while preserving the colour and clarity of the view
- High energy efficiency resulting in 3X the examination time with each charge cycle of a Welch Allyn rechargeable handle
- 7-year warranty

**03100-LED** [WA03100LED](#) for Traditional Otoscope (#20000, #20200, #21700, #25020)  
**06500-LED** [WA06500LED](#) for MacroView™ Diagnostic Otosopes (#23810, #23820, #23920)  
**04900-LED** [WA04900LED](#) for Traditional Coaxial Ophthalmoscopes (#11720, #11730, #11735)  
**03800-LED** [WA03800LED](#) for PanOptic Ophthalmoscopes (#11810, #11820)



### Welch Allyn Batteries

**72200** [WA72200](#)  
3.5 V Nickel-Cadmium Rechargeable Battery (Black lettering) for power handles #71670, #60835, #71000  
**72300** [WA72300](#)  
3.5 V Nickel-Cadmium Rechargeable Battery (Orange lettering) for power handles #71000-A, #71000-C, #23300  
**71960** [WA71960](#)  
3.5 V Lithium-Ion Rechargeable Battery




## Welch Allyn

HANDLES & CHARGERS

### iExaminer SmartBracket

Enhance ear and eye exams with the iExaminer SmartBracket. Simply combine your Welch Allyn PanOptic Plus Ophthalmoscope or MacroView Plus Otoscope with the iExaminer SmartBracket and iExaminer Pro App to capture and share images and video.

- Allows you to easily switch to a digital exam, simply attach the iExaminer to your device and launch the iExaminer Pro App
- Designed to be compatible with most Smart Phones
- Can be connected to your device in four different directions for increased comfort during examination
- Includes 1 SmartBracket Clip



\$150.05

RRP \$142.50

WA11840IEX

SmartBracket Clip  
Pkt/2

WA11840SC  
\$36.20

How to Set Up iExaminer

SCAN TO WATCH



### Universal Desk Charger

Compact, portable unit provides convenient drop-in recharging of Welch Allyn desk/well-charger instrument power handles.

- Integrated light indicators highlight that power is flowing to the unit and that handles are charging
- Now featuring iExaminer® Pro SmartBracket™ storage
- Designed for Welch Allyn Lithium-Ion, Lithium-Ion Plus battery handles as well as 3.5 V or 2.5 V NiCad, Lithium Ion and PocketScope™ handles



Desk Charger Set WA71XXLDE6

Universal Desk Charger with 2 x Li-Ion Plus USB C Handles

\$1890

Desk Charger WA719DSK6

Universal Desk Charger

\$700

### Welch Allyn Handles and Chargers

- 71900-USB WA71900USB 3.5V Li-Ion Handle with USB Charging Module
- 71906 WA71906 3.5V Li-Ion Handle & Charging Pod
- 71946 WA71946 Transformer & Charging Pod for Li-Ion Power Handle
- 71036 WA71036 Charging Transformer for 3.5V Power Handles



### Li-Ion Plus Handle

Optimises diagnostic tools, staying bright over the life of a battery charge with no colour change when dimmed.

- 5.5 hours on time when fully charged
- Fast Charge to 75% capacity in 20 minutes
- Auto-off feature
- Precise light control 5-stop settings
- Low battery power indicator
- USB-C Charging port
- Suitable for Welch Allyn 3.5V Otoscope and Ophthalmoscope Heads



\$845.00

WA7193

For the best performance optimise your Welch Allyn Diagnostic Exam Tool with the latest accessories.

**Baxter**

**Discover why more Doctors are choosing Riester.  
German Quality at an affordable price.**



Riester is able to offer customers the ultimate in product value in the form of high-tech functionality, ease of operation, reliability and durability at affordable prices.

**Riester e-scope®**

**Pocket Otoscope 2.7V**

- Vacuum direct illumination (2500K)
- 3x magnifying pivoting lens
- Approx. 3 hours use per set of AA batteries
- Includes reusable ear specula and soft pouch

ECOR12101200 Black  
ECOR12100200 White

**\$125**



**Riester e-scope®**

**Pocket Otoscope 3.7V LED F.O.**

- LED Fibre Optic illumination (5500K)
- 3x magnifying pivoting lens
- Approx. 50 hours use per set of AA batteries
- Includes reusable ear specula and hard case

ECOR12110203 Black  
ECOR12111203 White

**\$245**



**Riester e-scope®**

**Pocket Diagnostic Set 2.7V**

- Includes:
- e-scope® Ophthalmoscope
  - e-scope® Otoscope
  - Reusable Ear Specula
  - Soft storage bag



**\$320**

ECOR12131200 Black  
ECOR12130200 White



**Riester® Ri-Former LED Diagnostic stations**

A modular system offering up to 5 instrument-handles, with specula-dispenser, blood-pressure devices and thermometers. The ri-former diagnostic station is available as wall-mountable and is connected to your power.

- Wall Mounted
- Maintenance free LED technology
- Versatile combination of option
- Otoscope and Ophthalmoscopes are anti-theft proof
- Modules can be combined individually to suit the requirements and location
- 3-year warranty

**Option 1**

Includes:

- 1x ri-former® 2 handles 3.5V / 230V
- 1x ri-scope® L2 otoscope head LED 3.5 V
- 1x ri-scope® L2 ophthalmoscope head LED 3.5 V

ECOR10863

**\$1025**



**Riester® Otoscope Tips**

2.5mm **ECOR10531**  
Pkt/100 **\$28.00**


4.0mm **ECOR10534**  
Pkt/100 **\$28.00**



**Dispenser Ri-Spec Ear Specula Dispenser**

ECOR13654

**\$95.00**



**Customer Service Champions**

Our friendly and dedicated Customer Service Team is one of the best things about Team Medical. They are renowned for their commitment to customer service excellence.

Our team operates from Sydney, Brisbane, Melbourne, Cebu and Hobart. Wherever you are we are here to help. We cover multiple time zones and are available from 8.30am to 5pm (Friday till 4.30pm), Perth 6.30am to 3pm (Friday till 2.30pm).

**For excellent service, choose Team Medical.**

