

Soffban® Natural



Orthopaedic Padding

- Conforms well for easy shaping
- Natural fleece absorbs sweat and exudate
- Ideal for plaster casts
- Undercast padding
- 2.7m



| | | | | | |
|-------|-----------|--------|---------|----|---------|
| 5cm | SN7147200 | Pkt/12 | \$38.40 | 2+ | \$36.00 |
| 7.5cm | SN7147300 | Pkt/12 | \$41.75 | 2+ | \$38.45 |
| 10cm | SN7146900 | Pkt/12 | \$53.20 | 2+ | \$49.00 |
| 15cm | SN7147000 | Pkt/12 | \$62.10 | 2+ | \$57.20 |

Soffban® Synthetic



Water-Resistant Orthopaedic Padding

- Water-resistant
- Ideal for synthetic casts
- 2.7m



| | | | | | |
|-------|------------|--------|---------|----|---------|
| 5cm | SN36000025 | Pkt/12 | \$37.90 | 2+ | \$35.25 |
| 7.5cm | SN36000026 | Pkt/12 | \$41.75 | 2+ | \$38.45 |
| 10cm | SN36000027 | Pkt/12 | \$53.20 | 2+ | \$49.00 |
| 15cm | SN36000028 | Pkt/12 | \$62.10 | 2+ | \$57.20 |



Delta-Dry®



Waterproof Cast Liner

- Water resistant cast padding (2.4m)
- Allows patients to shower, swim and bathe like normal
- Better long term water resistance
- Offers more under cast protection
- Eliminates the need for stockinette



| | | | |
|-------|-----------|--------|----------|
| 5cm | SN7344300 | Box/12 | \$330.80 |
| 7.5cm | SN7344301 | Box/12 | \$398.65 |
| 10cm | SN7344302 | Box/12 | \$458.10 |

Velband



Absorbent Orthopaedic Padding (2.7m)

Velband is a hypoallergenic underbandage with a unique three-layer structure for easy application and long-term comfort.

It provides excellent cushioning and moisture management under plaster of Paris or synthetic casts, avoiding irritation and pressure sores by keeping the skin dry and comfortable.



| | | | | | |
|-------|-----------|--------|---------|----|---------|
| 7.5cm | SN7231501 | Pkt/12 | \$40.40 | 2+ | \$37.20 |
| 10cm | SN7231502 | Pkt/12 | \$48.55 | 2+ | \$44.70 |
| 15cm | SN7231504 | Pkt/12 | \$69.80 | 2+ | \$65.45 |

Pure Cotton



Orthopaedic under-bandages

Made from high quality, long staple cotton fibres.

Designed and made to have some self-adhesion property and appropriate tensile strength, resistant to abrasion and plaster strike through.

2.7m



| | | | | | |
|-------|----------|--------|---------|----|---------|
| 5cm | MGT30700 | Pkt/12 | \$16.45 | 3+ | \$15.00 |
| 7.5cm | MGT30701 | Pkt/12 | \$19.80 | 3+ | \$17.75 |
| 10cm | MGT30702 | Pkt/12 | \$32.35 | 3+ | \$28.95 |
| 15cm | MGT30703 | Pkt/12 | \$44.30 | 3+ | \$39.75 |

Protouch®

Synthetic and Cotton Stockinette

- Soft and comfortable providing good cushioning and comfort to the patient.
- Stretches to three times its normal width to conform to an arm, leg or body.

Protouch Synthetic Stockinette

- Made of synthetic fibres to react moisture away from the skin for rapid evaporation

| | | | |
|-------|----------|------|---------|
| 5cm | SN301002 | Roll | \$46.05 |
| 7.5cm | SN301003 | Roll | \$54.80 |
| 10cm | SN301004 | Roll | \$91.40 |



Protouch Cotton Stockinette

- Ideal material for an initial undercast layer
- Made of a high-grade knitted cotton

| | | | |
|-------|----------|------|---------|
| 5cm | SN307002 | Roll | \$46.05 |
| 7.5cm | SN307003 | Roll | \$41.90 |
| 10cm | SN307004 | Roll | \$49.55 |



Cutimed® Sorbact®

Safe and effective wound management



Safe

Binds
bacteria
and fungi

Easy

No known
contra-
indications

Due to the **purely physical mode of action**, Cutimed® Sorbact® binds and removes bacteria and fungi without the release of any antimicrobial agents.

Safe and effective wound infection management;

- ◆ For all contaminated or infected hard to heal wounds eg. Venous and diabetic ulcers.
- ◆ Reduces the bacterial load and supports the natural wound healing process.
- ◆ No undesirable side effects, no cytotoxicity or development of bacterial resistance.
- ◆ Binds common microbes in the wound bed, including multi-resistant bacteria.
- ◆ No mechanism of antimicrobial resistance has been described with Sorbact.
- ◆ No active agents are released into the wound.



For more information,
contact your local Essity Representative
or Customer Service on:
1300 276 633, customerservice.AU@essity.com



Cutimed®,
an Essity brand

ESS0552 (02/24)

Always read the label and
follow the directions for use.



Actimove®

Feel good. Move better.



Actimove
Walker

Designed to stabilise and protect the lower leg, ankle, and foot

- High patient comfort due to the lightweight design of the walker –at least 20% lighter than previous version
- Slip resistant, well-balanced rocker sole ensures walking confidence and promotes natural gait
- Improved comfort through minimised impact thanks to the shock absorbing inner sole
- Washable foam lining improves comfort and promotes good hygiene
- Available in High or Low and Standard or AIR



Always read the label and follow the directions for use.



Actimove®,
an Essity brand ESS0847 02/2024

BSN medical Australia & New Zealand an Essity company

PO Box 337,
Mount Waverley,
Victoria, Australia 3149
www.essity.com

Customer Service Australia
T 1300 276 633 F 1300 998 830
Customer Service New Zealand
T 0508 276 111 F 0508 998 830

Gypsona®

Plaster of Paris Bandage

Gypsona bandages combine fast setting, high quality plaster of Paris with the conformable and stable interlocked weave of leno cloth. The bandage demonstrates minimum plaster loss whilst the rich creamy texture allows easy handling.

| | | | | | |
|-------|------------|--------|----------|----|----------|
| 5cm | SN87060527 | Box/12 | \$85.95 | 2+ | \$79.15 |
| 7.5cm | SN87067535 | Box/24 | \$184.55 | 2+ | \$169.90 |
| 10cm | SN87061035 | Box/24 | \$238.90 | 2+ | \$219.95 |
| 15cm | SN87061535 | Box/12 | \$168.50 | 2+ | \$155.20 |



Delta-Lite Plus®

Fibreglass Casting Bandage

- Easy to apply
- Greater conformability
- High lamination: smooth finished cast which is strong, long lasting and lightweight

| | | | | | |
|--------|-----------|------|---------|-----|---------|
| 5cm | SN7345801 | Roll | \$16.85 | 10+ | \$15.75 |
| 7.5cm | SN7345802 | Roll | \$19.00 | 10+ | \$17.80 |
| 10cm | SN7345803 | Roll | \$24.30 | 10+ | \$22.55 |
| 12.5cm | SN7345804 | Roll | \$27.40 | 10+ | \$25.70 |

Also available in following colours: red, blue, purple, black, pink, green, mixed. Note that we only sell them in full cartons of 10. Price is equivalent to 10+ price of standard white. Colours are not a standard stock line.



Delta-Cast Elite®

Polyester Casting Tape

Delta-Cast Elite has been designed with ease of application and patient comfort in mind. Our patented fibreglass-free casting material combines strength and rigidity with soft edges for enhanced patient comfort.

- Smooth cast with soft edges
- Fibreglass free
- Easy to apply and remove
- Highly conformable

| | | | | | |
|--------|-----------|------|---------|-----|---------|
| 5cm | SN7227900 | Roll | \$18.65 | 10+ | \$17.35 |
| 7.5cm | SN7227901 | Roll | \$21.10 | 10+ | \$19.65 |
| 10cm | SN7227902 | Roll | \$25.75 | 10+ | \$23.95 |
| 12.5cm | SN7227903 | Roll | \$31.30 | 10+ | \$29.15 |

Also available in following colours: red, blue, purple, black, pink, green, mixed. Note that we only sell them in full cartons of 10. Price is equivalent to 10+ price of standard white. Colours are not a standard stock line.



3M™ Scotchcast™ Plus

Synthetic Casting Tape


Consists of a knitted fiberglass fabric impregnated with polyurethane resin. Exposure to moisture or water initiates a chemical reaction which causes the tape to become rigid. Finished casts produced are lightweight, strong, radiolucent, and porous. Scotchcast Plus casting tape contains a lubricating agent in the resin to provide a slippery feel when exposed to water.



| | | | | | | | | | | | |
|-------|---------|------|---------|-----|---------|--------|---------|------|---------|-----|---------|
| 5cm | 3M82002 | Roll | \$15.25 | 10+ | \$13.80 | 10cm | 3M82004 | Roll | \$20.10 | 10+ | \$18.50 |
| 7.5cm | 3M82003 | Roll | \$18.00 | 10+ | \$16.30 | 12.5cm | 3M82005 | Roll | \$23.00 | 10+ | \$21.15 |

Collar & Cuff Arm Sling

6m ROWCC06W
\$25.05



Slings

Sling White Disposable

AXP10106312
90 x 90 x 120CM

Each \$0.80



Cotton Slings

AXP10106211

110 x 110 x 155cm

Each \$1.65

12+ \$1.50



Safety Pins

Assorted Sizes

LIVSP12L

Pkt/12 \$0.60



Actimove® Sling

Collar and Cuff Support

The padded strap and Y-tab closure system provide effective vertical support to the collarbone and upper limb area.

- Extra patient comfort with padded strap
- Y-tab hook system facilitates secure clinical positioning
- Convenient dispenser box presentation
- 28 Y-tabs are included in each box



essity

Actimove®
Feel good. Move better.

5cm x 12m

SN7285918

Box/2rolls \$81.25



Actimove® Mitella

Comfort Arm Sling

The envelope sling and neck strap of the Actimove Mitella provide effective vertical support to the arm.

- Extra patient comfort with neck pad
- Adjustable neck strap, with Y-tab hook system facilitates secure clinical positioning
- Soft thumb loop maintains sling placement

Actimove®
Feel good. Move better.

Small SN7281961

Medium SN7281962

Large SN7281963

X-Large SN7281964

\$52.25

essity

JOBST® Travel Socks

Medical Compression Stockings



See page 412



Actimove® Lower Leg Walker

This walker is expertly designed for stability and protection of the lower leg, ankle, and foot. With a slip-resistant profile sole, it provides confidence on various surfaces and is at least 20% lighter than its predecessor for added comfort. The conforming plastic uprights restrict leg motion, ensuring stability, while non-slip strapping and a shock-absorbing inner sole enhance overall comfort. The well-balanced rocker sole promotes a natural gait, and dense foam lining ensures reliable walking support.

x-sml SN7627035

sml SN7627036

med SN7627037

lrg SN7627038

x-lrg SN7627039

\$83.95

Actimove®
Feel good. Move better.



| Standard | | | UK | US | EU |
|----------|----|---|---------|-----------|-----------|
| 76270-35 | XS | M | <3.5 | <4.5 | <36.5 |
| | | F | <3 | <6 | <36 |
| 76270-36 | S | M | 4-6.5 | 4.5-7 | 36-39.5 |
| | | F | 3.5-5.5 | 6-8 | 36.5-38.5 |
| 76270-37 | M | M | 7-9.5 | 7.5-10 | 40-44 |
| | | F | 6-8.5 | 8.5-11 | 39-43 |
| 76270-38 | L | M | 10-12 | 10.5-12.5 | 44.5-45.5 |
| | | F | 9-11 | 11.5-13.5 | 44-46 |
| 76270-39 | XL | M | >12 | >12.5 | >46.5 |
| | | F | >11 | >13.5 | >46 |



Actimove®

Tri-Panel Knee Immobilizer

Three Panel Knee Immobilizer provides effective immobilisation of the knee and helps to reduce pain following trauma or surgery. The tri-panel design allows splint adjustments for a broad spectrum of patient leg circumferences and the hook and loop sensitive closure system permits easy application and removal.

50cm sml SN7281914 \$123.55

50cm lrg SN7627038 \$133.55

Actimove®
Feel good. Move better.



Cervical Collar Adjustable Extraction Collar

Moulded from high density polyethylene with foam padding and a patented chin support. This collar restricts cervical spine flexion, extension and rotation to promote patient recovery.

The one-piece design enables efficient storage where space is limited. The multi-height adjustment fits almost all adults (four positions) and paediatrics (three positions).

The one-piece design can easily be used; simply measure patient, adjust and lock selected height and apply. For either a supine or seated position, application is quick and easy.



ADTMSEX700 \$48.05*

Cervical Collar

All 48cm length

Small 7.5cm SUTSM400
Medium 9cm SUTSM401



\$13.85*

Dynacast® Prelude

Prelude Single Layer Ortho Glass

Dynacast Prelude is a synthetic splinting system specifically incorporating fibreglass, covered by a polypropylene padding, to offer easy application for the support and immobilisation of fresh fractures and soft tissue injuries. It is presented in a roll format for versatility in use and application.

- Versatile splinting system for individual patient application
- Clean modern alternative to plaster of Paris slabbing
- Strong and light with no plaster mess
- All-in-one roll is easy and quick to apply

| | | | |
|---------------|----------------|------|----------|
| 2.5cm x 4.6m | SN7141800 | Roll | \$203.75 |
| 5cm x 4.6m | SN7142000 | Roll | \$266.95 |
| 7.5cm x 4.6m | SN71421001 | Roll | \$370.00 |
| 10cm x 4.6m | SN71415001 | Roll | \$573.85 |
| 12.5cm x 4.6m | SN714160000001 | Roll | \$550.60 |
| 15cm x 4.6m | SN714170000001 | Roll | \$619.15 |



Finger Stax splints



| | | | |
|-----|--------|------|--------|
| | | Each | \$5.90 |
| Sz1 | MIAS01 | Sz5 | MIAS05 |
| Sz2 | MIAS02 | Sz6 | MIAS06 |
| Sz3 | MIAS03 | Sz7 | MIAS07 |
| Sz4 | MIAS04 | | |

Finger splints

Aluminium



| | | | | |
|--------------|---------|--------|----|--------|
| 13mm x 480mm | ROW4813 | \$7.25 | 2+ | \$5.70 |
| 19mm x 480mm | ROW4819 | \$7.50 | 2+ | \$5.95 |
| 24mm x 480mm | ROW4824 | \$8.35 | 2+ | \$6.55 |
| 32mm x 480mm | ROW4832 | \$9.10 | 2+ | \$7.30 |

Sam Splint

Lightweight, reusable and very portable

The Sam Splint Original is an extremely pliable multi-purpose splint and is built from a thin core of aluminium alloy, sandwiched between two layers of closed cell foam. Bent into various shapes and curves, it becomes extremely strong and supportive for any fractured or injured limb.

9.5cm x 4cm
FERSVPSSF \$5.20

91cm x 10mm
FERSVPSS36 \$30.05

