

HEINE Gamma G7 Handheld Sphygmomanometer

The model with aluminium alloy construction and precision valve

The GAMMA G7 is unique in its advanced construction, utilising a housing made of aluminium alloy and advanced thermoplastics, producing a shockproof and modern design.

- Shockproof
- Double-walled housing with three-point manometer mounting
- Pressure-proof manometer ensures durability and reliability
- Easy to read 56mm diameter black-on-white scale is
- Large insufflation bulb for fast inflation of the cuff
- Single-tube design with exclusive pressure-distribution system
- Light alloy, adjustable spoon for both right- or left-handed operation
- Latex-Free
- 10-year warranty

RRP \$272.00

HEM00009232

\$259



Gamma G5

The HEINE standard sphygmomanometer

Combines all the design and quality features of the GAMMA G7 except:

- Housing is moulded of a 2 component thermo plastic
- Adjustable spoon (Polyamid Thermo Plastic)
- 5-year warranty

RRP \$186.00

HEM00009230

\$177



HEINE Gamma XXL LF®

The versatile professional model

The Classic GAMMA XXL LF is latex-free (LF) version with a range of mounting options that makes it the ideal sphygmomanometer in any clinical setting.

- Latex-free
- Durable polycarbonate construction will hold up in the busiest practice environments
- Oversized 130mm diameter scale can be pivoted for clear, easy reading from any angle
- Long spiral inflation tube (3m). For increased working distance
- High quality indicator is precisely calibrated and durable
- Precision air-release valve. Fine and precise air release
- Pressure-proof manometer. Ensures durability and reliability
- 10-year warranty

Desk model	HEM00009322	RRP \$331.00	\$315
Wall model	HEM00009323	RRP \$342.00	\$335
Mobile model	HEM00009326	RRP \$704.00	\$669



GAMMA G7 / G5 / GP Cuff & bladder		Size Circumference cm	
Child	HEM00009702	10x34 13-20	\$80.20
Adult small	HEM00009704	14x53 20-29	\$79.20
Adult	HEM00009619	14x58 29-41	\$82.20
Adult large	HEM00009706	17x70 35-47	\$105.95
Thigh	HEM00009615	20x86 41-59	\$106.95

GAMMA XXL LF Cuff & bladder		Size Circumference cm	
Child	HEM00009703	10x34 13-20	\$85.15
Adult small	HEM00009705	14x53 20-29	\$85.15
Adult	HEM00009629	14x58 29-41	\$87.20
Adult large	HEM00009631	17x70 35-47	\$104.20
Thigh	HEM00009625	20x86 41-59	\$112.25



Mercurial
Sphygmomanometer

An accurate, high quality
desktop Sphygmomanometer

- Measuring range : 0-300mmHg
- Graduation: 2mmHg
- Accuracy: ±3mmHg
- PVC bladder and bulb
- Aluminium case
- Nylon cuff with metal ring
- Latex free
- CE mark and TGA approved



ECOSPH0485
Light Grey

\$245

AND Mercury-free
Sphygmomanometer UM-101

The UM-101 was designed to be used as an alternative to the mercury column, while at the same time preserving the local environment and creating an easy-to-use product for hospital personnel.

- Auscultatory blood pressure measurement without mercury
- Digital read out for pulse rate
- Adjustable display column for easier analysis of results
- Portable or mobile stand options to suit diverse requirements
- Includes medium size cuff (23-33cm)
- 12 month warranty

UM-101A
Desk type
ANDUM101AW
\$215

UM-101B
Mobile stand type
ANDUM101BW
\$215

Mobile Base
For use with UM-101B
ANDUMST001
\$255



Riester

Big Ben® Sphygmomanometer

- Generously dimensioned scale 147mm diameter with high-contrast scale layout for optimum readability
- Spacious basket on the reverse side for tubes and cuff
- Spiral tube can be extended to 3m
- Non-wearing air-release valve with fine adjustment.
- Specially hardened copper-beryllium membrane, resistant to aging, can withstand pressures of up to maximum 600mm Hg
- Micro-filter to protect the release valve and measurement system
- Latex free
- 2-year warranty

Desk model ECORI1453LF \$239
Wall model ECORI1459LF \$195
Mobile model ECORI1468LF \$319



Riester

Riester R1 shock-proof®
Sphygmomanometer

- Fast, ergonomic inflation due to optimised bulb design with integrated spoon
- Maximum error tolerance: ± 2 mmH
- Includes desk storage caddy and 3 cuffs (adult, large adult and child)
- Latex free



ECORI1255LF
\$250.40


Riester

Riester Babyphon
Blood pressure monitor designed for children and babies

Great set all in a soft blue! This palm-style sphygmomanometer is ideal for paediatrician's offices and hospitals.

- Choice of either chrome or sturdy plastic/metal casing
- Ergonomic spoon-shaped handle made of stainless steel
- Precision pressure-release valve, permits fine regulation and is non-wearing
- Micro-filter protects pressure-release valve and measurement system
- Membrane can withstand overpressures of up to 600 mm Hg
- Linear scale, up to 300 mm Hg
- Includes newborn, infant and child cuff
- 2-year warranty

ECORI1440LF \$185





Aneroid Regal Sphygmomanometers

- High quality precise manometer
- Large 150mm faceplate
- Easy to read blue background with white writing
- Non-reflective glass
- Cuff, bladder, bulb and valve included
- Wall unit with 180° swivel bracket
- Mobile model with two-piece adjustable stand
- 12 month warranty

Desk model	BYDI580	\$159
Wall model	BYDI568	\$159
Mobile model	BYDI579	\$349



DESK



WALL



MOBILE



SPECTRUM - Professional Aneroid

- 300mmHg manometer with no pin stop
- Complete with heavy duty 2 tube inflation bladder, bulb and valve
- Available with or without stethoscope included
- Precision crafted premium sphygmomanometer for professional use
- Packed in oversized zippered nylon storage case
- 12 month warranty

without stethoscope	BYDI405	\$56
with stethoscope	BYDI406	\$75.90



ECONOM - Economy Aneroid

- 300mmHg manometer with no pin stop
- Complete with heavy duty 2 tube inflation bladder, bulb and valve
- Available with or without stethoscope included
- Economically priced precision sphygmomanometer
- Supplied in leatherette storage case
- 12 month warranty

without stethoscope	BYDI569	\$34.10
with stethoscope	BYDI570	\$49



Aneroid Palm held

- Economical
- Palm style manometer
- Heavy duty 1 tube latex inflation bladder
- Chrome plated inflation valve
- Supplied in nylon zipper case
- 12 month warranty

BYDI401	\$55.20
---------	---------



The Palm Lite

Shock Proof Aneroid Sphygmomanometer with enhanced performance

- Precision Crafted 300mmHg palm style manometer with enhanced shock resistant housing
- Expeditious control trigger release valve
- Complete with nylon cuff and single tube latex inflation bladder
- Quick release luer slip connector
- Supplied in large deluxe nylon zipper case
- 2-year warranty

BYDI591	\$114.20
---------	----------



Welch Allyn®

767 Series Sphygmomanometers

The 767 Series Sphygmomanometer is a smart choice for accurate and reliable readings throughout your facility. Connects to the FlexiPort Reusable Blood Pressure Cuff Range.

- Improved ergonomics and a new modern design
- Lifetime calibration warranty, reducing the cost of ownership
- Latex- and Mercury-free
- Certified accuracy to +/- 3mmHg

Desk Model WA767016 \$260 Adult Cuff included.

Mobile Model WA767010 \$479



10YR warranty



Welch Allyn®

DS-66 Gold Series Trigger Aneroid Sphygmomanometer

The robust design of the DS-66 with gear-free DuraShock™ technology ensures calibration remains longer than traditional aneroid sphygmomanometers, providing more accurate readings over time and a lower lifetime cost, helping to improve practice efficiencies and your bottom line.

- Laser-engraved dial face for accuracy
- Soft comfort grip and trigger deflation valve for smooth, one-finger operation
- Latex-free for safety
- Certified accuracy to +/- 3mmHg

WA509830 \$365

Kit Includes: Large Adult Cuff, Adult Cuff, Small Adult Cuff and Child Cuff, Sphygmomanometer and Zipper Case.



10YR warranty



Welch Allyn®

FlexiPort® Reusable Blood Pressure Cuffs

- FlexiPort technology enables every cuff to be used as a one- or two-tube cuff
- Port can be rotated for improved patient comfort
- Materials do not contain natural rubber latex
- Meet all latest clinical guidelines for proper fit from AAMI and AHA
- Primary sizes are color-coded to make it easy to find the proper cuff size



Save time when you standardise with a cuff solution that works across your blood pressure devices and patient monitors.




Small Infant	7-10cm	WAREUSE06	\$46.15
Infant	9-13cm	WAREUSE07	\$51.15
Small Child	12-16cm	WAREUSE08	\$53.10
Child	15-21cm	WAREUSE09	\$55.30
Small Adult	20-26cm	WAREUSE10	\$57.65
Adult	25-34cm	WAREUSE11	\$60.55
Adult Long	25-34cm	WAREUSE11L	\$63.50
Adult Large	32-43cm	WAREUSE12	\$85.00
Adult Large Long	32-43cm	WAREUSE12L	\$89.70
Thigh	40-55cm	WAREUSE13	\$105.40



Over the past six months our order accuracy rate has been 99.97%.
We are constantly advancing our warehouse processes to ensure quick turnaround times and accurate processing of orders.
Our Warehouse Teams are committed to going the extra mile.

Get what you need when you need it from Team Medical.



WE ARE A BUSINESS FOR GOOD



Our procurement manager, Matt reached the milestone of 11 years with Team Medical earlier this year. To celebrate, Matt was given the opportunity to choose a B1G1 project to support. Matt choose to give **5,475 days of life-saving clean water** to a community in Tanzania. This donation has allowed additional water pumps with sanitary seals to be built. Making clean water more accessible in this poverty-stricken area.

Learn more about our giving impacts at www.teammed.com.au/teamgiving



Welch Allyn®

ProBP™ 2000 Digital Blood Pressure Device

The new Welch Allyn® ProBP™ 2000 Digital Blood Pressure device was designed with advanced technology features for results you can feel confident in. Combining Welch Allyn's SureBP® technology and FlexiPort® snap-fit cuff connectors, one touch is all you need to get a fast, reliable blood pressure reading in 20 seconds.

- Designed with patient comfort in mind - SureBP® technology captures readings quickly while the cuff inflates
- Meets AAMI Standards for blood pressure accuracy
- Compatible with Welch Allyn FlexiPort® Reusable Blood Pressure Cuffs
- Adult cuff included



WA2000

\$402.55



Welch Allyn®

Connex® ProBP™ 3400 Digital Blood Pressure Device

Capture fast, accurate readings in the palm of your hand with the Connex® ProBP™ 3400 Digital Blood Pressure Device. Enhanced with SureBP® NIBP technology and compatible with FlexiPort® cuffs, this device captures readings in 15 seconds, ensuring greater patient comfort. Get more confident blood pressure readings the easier way with this practical digital upgrade, a smart choice for clinicians in any setting.

- Small, handheld design makes device easy to carry and capture readings
- SureBP® technology captures readings quickly and comfortably while the cuff inflates
- Lithium ion battery technology can take up to 100 readings on a full charge
- AAMI/BHS-AA certified blood pressure accuracy
- Compatible with Welch Allyn FlexiPort® Reusable Blood Pressure Cuffs
- Adult and Large Adult cuff included
- 3-year warranty



Handheld

WA34XFHT6

\$1055

Mobile Model

WA34XFST6

\$1495



BP Medical (WA) | www.teammed.com.au
PHONE 1300 33 55 60 | info@bpmedical.com.au

All prices are subject to 10% GST unless marked with * (GST Free)



Team Medical | www.teammed.com.au
PHONE 1300 22 44 50 | info@teammed.com.au



X5 BT "PARR" Automatic BP Monitor

- Bluetooth® low energy technology
- Arrhythmia Detection (ARR)
- Atrial Fibrillation Detection (AFib)
- Premature Contraction Detection (PC)
- Movement Detection
- Cuff Detection
- Hypertension Risk Indication
- Average of the last three readings
- 7 day AM/PM Averaging
- Patented Latex-Free Universal Cone Cuff
- 24-36cm Cuff included
- Lifetime Calibration

BYDRMX5BT

\$139



Rossmax Cuffs

Small 18-26cm	BYDRMBPCF003	\$26.95
Medium 24-36cm	BYDRMBPCF002	\$26.95
Large 34-46cm	BYDRMBPCF001	\$26.95

Professional Automatic Blood Pressure Monitor

Model HBP-1320

With a rugged design and a practical carrying handle, the HBP-1320 is designed for use in professional settings and easily carried in everyday use. The device has been clinically proven to deliver fast and reliable results.

- Dual measurement modes - oscillometric and auscultation
- Motion stop function, device stops deflation for 5 seconds when body movement is detected
- Irregular heartbeat detection
- Long battery life - 300 measurements per charge
- Rechargeable battery
- Power saving auto power off
- Highly durable pump/sensor
- Compact and light weight with carry handle
- Includes M-cuff (22 - 32 cm), L-cuff (32 - 42 cm), rechargeable battery and power supply
- 3-year warranty

DAVHBP1320

\$515



Cuffs to suit HBP-1320

Small	17-22cm	DAVHBP1320CUFFS	\$137.60
Medium	22-32cm	DAVHBP1320CUFFM	\$154.30
Large	32-42cm	DAVHBP1320CUFFL	\$175.20
X-Large	42-50cm	DAVHBP1320CUFFXL	\$203.70

MEET OUR TEAM: TIM
MANAGEMENT - VIC

As National Sales Manager, Tim is a master of relationship building and a secret Excel whizz. Tim enjoys working alongside "fantastic staff" and driving social impact within Team Medical. He thrives on constantly supporting our community of customers and partners through their ongoing challenges — consistently focusing on positive outcomes.

Outside work Tim is a proud grandfather and loves to spend time outdoors flying his drone around his local beaches or looking after his gladiolus garden.

What sport do you like watching? AFL - Go Collingwood!

What is your favourite takeaway food? Chicken and Chips



HEM-907 Professional Blood Pressure Monitor



★★★★★ PREMIUM PRODUCT



DAVHEM907KIT
\$995

Kit Includes:
Blood Pressure Monitor
Small, Medium and Large Adult Cuffs and Bladders,
Cuff tube 1-metre length
Rechargeable battery pack and AC adaptor

2YR
warranty

Accurate and reliable

Complies with AAMI* standards for accuracy in automatic blood pressure measurement.

Easy for operation and training

With advanced fuzzy logic technology, the HEM-907 determines ideal cuff inflation level based on patient arm size and blood pressure fluctuation. Simple push button operation and automatic deflation speed makes training clinical staff fast and easy.

Dual measurement mode

AUTO Mode: the oscillometric measurement method
Cuff automatically inflates to the optimal pressure according to patient’s condition.

MANU Mode: the auscultation measurement method, using a stethoscope
The unit measures blood pressure with the use of a stethoscope.

Quiet, ideal for clinical environment

The HEM-907 operates quietly reducing stress on patients and clinical staff; perfect for night time use by your patient’s bedside.



HEM- 907 Accessories

Cuffs			
Small	17 - 22cm	DAVHEM907CUFFS	\$21.25
Medium	22 - 32cm	DAVHEM907CUFFM	\$22.10
Large	32 - 42cm	DAVHEM907CUFFL	\$35.75
Bladder			
Small	17 - 22cm	DAVHEM907BLDRS	\$24.20
Medium	22 - 32cm	DAVHEM907BLDRM	\$25.90
Large	32 - 42cm	DAVHEM907BLDRL	\$37.20
Cuff and Bladder Sets Includes 1m tubing			
X-Large	42 - 50cm	DAVHEM907CUFFXL	\$214.20



Spares		
Connector Plastic	DAVHEM90704	\$3.00
Air Tube 1 Mtr	DAVHEMTUBE100	\$14.40
AC Adaptor	DAVHEM907ADAPT	\$84.05
Battery Pack	DAVHEM907PBAT	\$52.90
Optional installation options		
Pole Mounting Kit [1]	DAVHEM907PKIT	\$399.00
Wall Mounting Kit [2]	DAVHEM907WKIT	\$245.00
Mobile Stand [3]	DAVHEM907STAND	\$365.00

Smart Elite+ Upper Arm Automatic Blood Pressure Monitor

Model HEM-7600T

Introducing the OMRON Smart Elite+ HEM-7600T Blood Pressure Monitor with Intellisense Technology and IntelliWrap Cuff. With the latest tubeless technology, 360° Accuracy and simple to use one-touch measurement.

- Cuff wrapping guide lamp
- One touch use
- Irregular heartbeat detection
- Body movement detection
- OMRON connect App makes it easy to view, save, and manage data
- Bluetooth® connectivity
- 360° IntelliWrap Cuff (Med-Lge 22-42cm)



OMRON

Intelli sense

Bluetooth®



DAVHEM7600T **\$304.20**

Bluetooth Wrist Blood Pressure Monitor

Model HEM-6232T

- One touch use
- Cuff wrapping guide
- Flexible wrist cuff 13.5 - 21.5cm
- Irregular heartbeat detection
- Morning average indicator
- Average for last 2-3 readings within 10 minutes
- Omron connect app makes it easy to view, save, and manage data
- Bluetooth® connectivity
- Measurement positioning indicator



OMRON

Intelli sense

Bluetooth®



DAVHEM6232T **\$245.05**

Standard Blood Pressure Monitor

Model HEM-7121

- One touch use
- Memory for 30 sets of readings
- Hypertension indicator
- Irregular heartbeat detection
- Body movement detection
- Accurate measurements with enhanced Intellisense Technology
- Portable and compact design
- Cuff wrapping guide
- 4-year warranty (when registered online)
- Soft Cuff (Med - Lrg)



OMRON

Intelli sense



DAVHEM7121 **\$109**



HEM7361T Advanced with AFIB Indicator Automatic Blood Pressure Monitor

The OMRON HEM-7361T detects and measures AFIB and High Blood Pressure simultaneously. It also features an irregular heartbeat indicator. Early diagnosis and detection of Atrial Fibrillation may help to prevent stroke or heart failure. You can download the free OMRON connect App to easily track readings. Dual check for dual users allows users to have a comprehensive check for their risk of stroke with one device, in one sitting.

Great for at home or at the clinic.

- AFIB Indicator
- Irregular Heartbeat Indicator
- Bluetooth Connectivity
- 360° IntelliWrap Cuff
- Intellisense Technology
- Dual Users x 100 Measurement memory
- Hypertension Indicator
- Cuff Wrapping Guide
- Body Movement Indicator
- Monitor Pulse
- Five year Standard Warranty

DAVHEM7361T

\$245



Automatic Blood Pressure Monitor

Model HEM-7156T

- 360° IntelliWrap™ Cuff
- Intellisense Technology
- Irregular Heartbeat Indicator
- Bluetooth Connectivity
- 60 Memory
- Hypertension Indicator
- Cuff wrapping guide
- Body movement Indicator
- Monitor Pulse
- 4x AA Batteries
- 5 Year Standard Warranty

DAVHEM7156T

\$175



Automatic Blood Pressure Monitor

Model HEM-7144T1

- Intellisense Technology
- Irregular Heartbeat Indicator
- Bluetooth Connectivity
- 30 Memory
- Hypertension Indicator
- Cuff Wrapping guide
- Body Movement Indicator
- Monitor Pulse
- 5 Year Standard Warranty



BP Medical (WA) | www.teammed.com.au
PHONE 1300 33 55 60 | info@bpmedical.com.au



DAVHEM7144T1

\$145



Team Medical | www.teammed.com.au
PHONE 1300 22 44 50 | info@teammed.com.au

Cuffs

Newborn	18.5L x 5.5W cm	BYDI453	\$10.85
Infant	26L x 7.5W cm	BYDI440	\$12.40
Child	34.5L x 11W cm	BYDI474	\$16.25
Adult	54L x 15W cm	BYDI486	\$17.15
Large Adult	61L x 17.5W cm	BYDI485	\$29.65
Thigh/Obese	82L x 21.5W cm	BYDI524	\$38.45



Bladders (2 Tube)

Newborn	8L x 4Wcm	BYDI454	\$15.70
Infant	11L x 5.5Wcm	BYDI444	\$16.45
Child	17L x 8.5Wcm	BYDI475	\$19.20
Adult	22.5L x 12.5Wcm	BYDI455	\$22.30
Large Adult	32L x 14.5Wcm	BYDI487	\$25.60
Thigh/Obese	38L x 18.5Wcm	BYDI523	\$40.55



Complete cuff & bladder sets (2 Tube)

Includes cuff, bladder, bulb, valve and connector

Child	34.5L x 11Wcm	BYDI469	\$41.25
Adult	54L x 15Wcm	BYDI470	\$57.25
Large Adult	61L x 17.5Wcm	BYDI471	\$62.35
Thigh/Obese	82L x 21.5Wcm	BYDI472	\$68.35



Complete cuff & bladder set Latex-Free (2 Tube)

Includes cuff, bladder, bulb, valve and connector

Child	34L x 11Wcm	BYDI506	\$35.75
Adult	53L x 15Wcm	BYDI507	\$42.85
Large Adult	62L x 18Wcm	BYDI508	\$50.35

LATEX-FREE

Full range of latex free spares is available.
Check our website or ask customer service



Bladderless cuffs (2 Tube)

- Fully washable
- Latex-free
- No need for separate cuff and bladder

Child	34L x 11Wcm	BYDI464	\$20.35
Adult	53L x 15Wcm	BYDI465	\$22.55
Large Adult	62L x 18Wcm	BYDI466	\$31.35



Spares



1	Coiled Tubing 240cm	BYDI578	\$18.15
2	Double ended connector plastic	BYDI525	\$2.40
3	Screw connector, Male, brass	BYDI526	\$5.45
4	Screw connector, Female, brass	BYDI527	\$5.45
5	Screw connector, Set, brass	BYDI528	\$10.95
6	Air release valve standard	BYDI574	\$6.00
7	Regular bulb with intake valve	BYDI577	\$6.90
8	Large bulb with intake valve	BYDI576573	\$8.00



Riester Spares

Riester Cuff/bladder LF child 1 tube	ECORIBPC13
Riester Cuff/bladder LF small adult 1 tube	ECORIBPC14
Riester Cuff/bladder LF adult 1 tube	ECORIBPC15

Riester Cuff/bladder LF child 2 tube	ECORIBPC23
Riester Cuff/bladder LF small adult 2 tube	ECORIBPC24
Riester Cuff/bladder LF adult 2 tube	ECORIBPC25
Riester Cuff/bladder LF obese 2 tube	ECORIBPC26

Riester inflation bulb latex free	ECORISPH2LF
Riester air release valve	ECORISPH3
Riester tube connector metal 'm/f set	ECORISPH4
Riester tube connector metal 'male'	ECORISPH5
Riester tube connector metal 'female'	ECORISPH6
Riester tube connector plastic 'male'	ECORI10387
Riester tube connector plastic 'female'	ECORI10388
Riester sphyg coil tubing latex free	ECORISPH7



Accessories for Omron Digital B.P. Monitors



Standard flat cuffs - A				
Child	17-22cm	DAVCUFFCHILD	\$54.45	to suit models IA2, T9P, HEM7211/7121/7130/7/7322
Medium	22-32cm	DAVCUFFMEDIUM	\$59.55	
Large	32-42cm	DAVCUFFLARGE	\$54.45	
Easy Fit Cuffs - B				
M to L	22-42cm	DAVCUFFML	\$54.45	to suit HEM7121/7130/7322/7280
Cuff to suit HEM7320 - C				
M to L	22-42cm	DAVCUFF7320	\$54.45	to suit HEM7320
Cuff to suit HBP1320 - D				
Small	17-22cm	DAVHBP1320CUFFS	\$137.60	
Medium	22-32cm	DAVHBP1320CUFFM	\$154.30	
Large	32-42cm	DAVHBP1320CUFFI	\$175.20	
X-Large	42-50cm	DAVHBP1320CUFFXL	\$203.70	



Extra-Large Blood Pressure Monitor Cuff (BPM-CUFFXL)

- Upper Arm Circumference 32 - 50cm
- Compatible with most automatic upper-arm blood pressure monitors when using the existing air connector provided with the device.
- 1 x Cuff & 1 x Air Tube

DAVBPMCUFFXL \$24.95

For more blood pressure spare parts — visit our website or call our team, if we don't have it we can source it for you www.teammed.com.au



MEET OUR TEAM: NAOMI
WAREHOUSE - BRISBANE

Naomi has a steady hand, having helped run our Brisbane warehouse for nearly 6 years. She loves playing cricket and collects NRL stuff (preferably Brisbane Broncos or Queensland Maroons, that is!) — and has quite a collection of Footy Cards.

What is one of your pet peeves?
People snorting



MEET OUR TEAM: ANTHONY
WAREHOUSE - SYDNEY

Anthony leads the team in our Sydney Warehouse, ensuring your deliveries are out on time. He loves being in planes and watching them; he hopes his next travel destination is Hawaii. When he's grounded on land, he enjoys watching the NRL, the Amazing Race, or Go Karting.

Which live show/concert have you attended lately? Pitbull. It was unreal!

MEET OUR TEAM: NATHAN
WAREHOUSE - PERTH

When not running the WA warehouse or spending time with 3-year-old labrador Meika, Nathan is a keen golfer. He also enjoys playing video games and watching AFL.

What is your favourite takeaway food? Varsity (the Perth burger joint)



MEET OUR TEAM: RAFAEL
WAREHOUSE - MELBOURNE

Rafael has spent over 15 years in the medical supplies industry and currently leads the operations at our busy Melbourne Warehouse. He is big on winter sports and tries to get amongst the action at the snow whenever he visits.

Where were you born? Brazil