


## Terumo Syringes

latex-free

The clear barrel with easy to read measurements allows for quick recognition of air bubbles and accurate dosing.

1ml	Slip Tuberculin Tip	TER6SS01T	Box/100	\$21.65*	10+	\$21.25*
3ml	Slip Tip	TER6SS03S	Box/100	\$12.95*	10+	\$12.75*
5ml	Slip Tip	TER6SS05S	Box/100	\$19.60*	10+	\$19.10*
10ml	Slip Tip	TER6SS10S	Box/100	\$24.45*	5+	\$22.55*
3ml	Lock Tip	TER6SS03L	Box/100	\$15.50*	10+	\$15.20*
5ml	Lock Tip	TER6SS05L	Box/100	\$23.00*	10+	\$22.45*
10ml	Lock Tip	TER6SS10L	Box/100	\$29.20*	5+	\$26.90*
20ml	Lock Tip	TER6SS20L	Box/50	\$33.45*	5+	\$30.85*
30ml	Lock Tip	TER1SS30LE1	Box/50	\$54.80*	5+	\$50.45*
50ml	Lock Tip	TER1SS50LE	Box/20	\$33.05*	5+	\$30.45*
20ml	Slip Eccentric Tip	TER6SS20ES	Box/50	\$24.10*	5+	\$22.20*
30ml	Slip Eccentric Tip	TER1SS30ESE1	Box/50	\$57.80*	5+	\$53.25*



**TERUMO**

## Syringes

1ml	Slip Tuberculin Tip	BD302100	Box/100	\$27.10*	10+	\$26.95*
2ml	Slip Tip	BD302204	Box/100	\$17.75*	10+	\$17.65*
5ml	Slip Tip	BD302130	Box/100	\$24.55*	10+	\$24.10*
10ml	Slip Tip	BD302143	Box/100	\$28.45*	4+	\$26.40*
1ml	Lock Tip	BD309628	Box/100	\$91.75*	10+	\$30.10*
3ml	Lock Tip	BD302113	Box/100	\$21.95*	10+	\$21.89*
5ml	Lock Tip	BD302135	Box/100	\$30.25*	10+	\$30.10*
10ml	Lock Tip	BD302149	Box/100	\$40.95*	4+	\$37.85*
10ml	Slip Eccentric Tip	BD302146	Box/100	\$27.75*	4+	\$25.30*
20ml	Slip Eccentric Tip	BD300613	Box/120	\$83.75*	4+	\$77.35*
30ml	Slip Eccentric Tip	BD301231	Box/60	\$60.50*	4+	\$56.70*
50ml	Slip Eccentric Tip	BD300145	Box/40	\$61.85*	4+	\$57.95*
50ml	Catheter Tip	BD300867	Box/60	\$104.40*	4+	\$97.80*



**BD**  
Advancing the world of health

## Ultra-Fine Insulin syringes

**BD**  
Advancing the world of health



0.3ml	29G x 1/2"	BD326103	Box/100	\$49.15*	5+	\$39.85*
0.3ml	31G x 5/16"	BD328822	Box/100	\$40.45*	5+	\$38.95*
0.5ml	29G x 1/2"	BD326769	Box/100	\$45.65*	5+	\$41.85*
0.5ml	30G x 1/2"	BD326725	Box/100	\$52.30*	5+	\$50.80*
0.5ml	31G x 5/16"	BD328821	Box/100	\$40.45*	5+	\$37.70*
1ml	29G x 1/2"	BD326719	Box/100	\$43.20*	5+	\$40.25*
1ml	31G x 5/16"	BD328820	Box/100	\$53.30*	5+	\$42.70*

## Insulin Syringes with Needles

latex-free

Terumo insulin syringes are designed to maximise patient comfort. Needle sharpness combined with a smooth plunger movement from the lubricated barrel, ensures only gentle injection pressure is required.




**TERUMO**

0.5ml	27G x 1/2"	TER3SS05M2713A	Box/100	\$28.50*	5+	\$26.25*
0.5ml	29G x 1/2"	TER3SS05M2913A	Box/100	\$29.25*	5+	\$26.95*
1ml	27G x 1/2"	TER3SS10M2713A	Box/100	\$28.50*	5+	\$26.25*
1ml	29G x 1/2"	TER3SS10M2913A	Box/100	\$29.25*	5+	\$26.95*

## Unisharp

### 0.3ml Ultra-fine Fixed Needle Syringe

- Ultra-sharp double bevel for easy skin penetration
- Double graduations in both ml and units
- Clear marks every 10 units for precise control when used for injection
- Low dead space design minimises waste
- Latex and DEHP free



Markings in ml and units

8mm needle

**TERUMO** **Latex-free**

31G x 8mm **TERU31PK**  
Box/100 **\$33.10\***

## Agani™ Hypodermic Needles

Agani™ needles feature an ultra sharp lancet point.

The precisely honed edges and unique bevel design easily penetrate tissue with an incision that minimises trauma. Silicone lubrication allows the needle to advance smoothly.

- |                                  |                                  |                                  |
|----------------------------------|----------------------------------|----------------------------------|
| ● 18G x 38mm <b>TER8AN1838R1</b> | ● 22G x 32mm <b>TER8AN2232R1</b> | ● 26G x 13mm <b>TER8AN2613R1</b> |
| ● 19G x 25mm <b>TER8AN1925R1</b> | ● 22G x 38mm <b>TER8AN2238R1</b> | ● 26G x 23mm <b>TER8AN2623R1</b> |
| ● 19G x 38mm <b>TER8AN1938R1</b> | ● 23G x 16mm <b>TER8AN2316R1</b> | ● 27G x 16mm <b>TER8AN2716R1</b> |
| ● 20G x 25mm <b>TER8AN2025R1</b> | ● 23G x 25mm <b>TER8AN2325R1</b> | ● 27G x 19mm <b>TER8AN2719R1</b> |
| ● 20G x 38mm <b>TER8AN2038R1</b> | ● 23G x 32mm <b>TER8AN2332R1</b> | ● 30G x 19mm <b>TER8AN3013R1</b> |
| ● 21G x 16mm <b>TER8AN2116R1</b> | ● 24G x 25mm <b>TER8AN2425R1</b> |                                  |
| ● 21G x 25mm <b>TER8AN2125R1</b> | ● 25G x 16mm <b>TER8AN2516R1</b> |                                  |
| ● 21G x 38mm <b>TER8AN2138R1</b> | ● 25G x 25mm <b>TER8AN2525R1</b> |                                  |

Box/100 **\$6.60\***  
10+ **\$6.20\***

### Drawing Up Needle (blunt needle)

- 18G x 1.5" **TER1XXNN1838** Box/100

**\$6.75\***



**TERUMO** **Latex-free**

## Needles

- |                                     |                                     |
|-------------------------------------|-------------------------------------|
| ● 18G x 1.5" (38mm) <b>BD302032</b> | ● 25G x 1" (25mm) <b>BD300400</b>   |
| ● 20G x 1" (25mm) <b>BD302025</b>   | ● 25G x 1.5" (38mm) <b>BD301808</b> |
| ● 21G x 1.5" (38mm) <b>BD302017</b> | ● 27G x 1/2" (13mm) <b>BD301801</b> |
| ● 23G x 1" (25mm) <b>BD300800</b>   | ● 30G x 1/2" (13mm) <b>BD304000</b> |
| ● 25G x 5/8" (16mm) <b>BD300600</b> |                                     |

Box/100 **\$8.35\*** 10+ **\$7.85\***



### Drawing Up Needle (blunt needle)

- 18G x 1.5" **BD300204**

Box/100 **\$8.75\***

## Agani™ Safety Needles

- |                                |                                |
|--------------------------------|--------------------------------|
| 25G x 16mm <b>TERSAN2516R1</b> | 22G x 38mm <b>TERSAN2238R1</b> |
| 26G x 13mm <b>TERSAN2613R1</b> | 23G x 25mm <b>TERSAN2325R1</b> |
| 18G x 38mm <b>TERSAN1838R1</b> | 23G x 38mm <b>TERSAN2338R1</b> |
| 19G x 38mm <b>TERSAN1938R1</b> | 25G x 25mm <b>TERSAN2525R1</b> |
| 21G x 38mm <b>TERSAN2138R1</b> | 25G x 38mm <b>TERSAN2538R1</b> |



Box/100 **\$20.05\*** 8+ **\$18.85\***



## Blunt Fill Needle 18G x 1.5"

- Specifically designed for safe injection preparation
- Minimises the risk of a needlestick injury
- 10 x the force to penetrate the skin, compared to a conventional needle

Box/100 **\$39.75\***  
10+ **\$37.35\***

**BD305180**



## Spinal Needles

**Latex-free**

- |                                |
|--------------------------------|
| ● 18G x 90mm <b>TER1SN1890</b> |
| ● 19G x 90mm <b>TER1SN1990</b> |
| ● 20G x 90mm <b>TER1SN2090</b> |
| ● 21G x 70mm <b>TER1SN2170</b> |
| ● 21G x 90mm <b>TER1SN2190</b> |
| ● 22G x 38mm <b>TER1SN2238</b> |
| ● 22G x 70mm <b>TER1SN2270</b> |
| ● 22G x 90mm <b>TER1SN2290</b> |
| ● 23G x 70mm <b>TER1SN2370</b> |
| ● 23G x 90mm <b>TER1SN2390</b> |
| ● 25G x 90mm <b>TER1SN2590</b> |

Box/25 **\$65.50\***  
5+ **\$58.85\***

38mm = 1 1/2"  
70mm = 2 3/4"  
90mm = 3 1/2"







## Terumo solutions for patient care

### AGANI™ Safety Hypodermic Needles

#### Reliable sharps protection with easy one-hand activation

- Agani safety hypodermic needles fit all luer lock syringes.
- Easily recognise the needle gauge by the colour-coded hub and safety shield.
- After withdrawing the needle from the patient, immediately activate the safety shield to minimise the risk of an accidental needlestick injury.
- 3 easy ways to activate the safety shield with one-hand:
  - Press the shield against a hard surface
  - Press the shield with your finger
  - Press the shield with your thumb



#### Terumo Australia Pty Ltd

Macquarie Park NSW 2113 Australia  
Mangere Auckland 2022 New Zealand  
E: cs\_australia@terumo.co.jp

T: 1800 837 866  
T: 0800 66 77 57

AGANI™ trademark of Terumo Corporation. © 2024 Terumo Australia Pty Ltd. MCS240220D1577V2 | PBD00413



## Surshield™ Safety Winged Infusion Sets

Reduce the risk of needlestick injuries with Surshield™

- Easy to use one-handed technique
- Large colour coded wings
- Flexible, memory-less tubing and transparent hub

- 30cm 21G ● TER5SVS21BL
- 30cm 23G ● TER5SVS23BL
- 30cm 25G ● TER5SVS25BL Box/50



\$103.05\*

## Surflo® Winged Infusion Sets Scalpvein Sets (butterfly needles)

latex-free

The precisely honed edges and double bevel design facilitate smooth penetration with minimal tissue trauma. The thin walled cannula allows for higher flow rates being ideal for rapid infusions and dosages with high viscosity. The tubing is highly flexible and virtually memory-less which greatly reduces risk of kinking and disruption to the flow.

### Short Tube 9cm

- 9cm 19G TER1SV19BLS

### Long Tube 30cm

- 30cm 19G TER1SV19BLK
- 30cm 21G TER1SV21BLK
- 30cm 23G TER1SV23BLK
- 30cm 25G TER1SV25BLK

Each \$1.20\*  
100+ \$1.00\*



## Surflo® IV Catheters

latex-free

The point of the introducing needle is honed to an ultrasharp triple-bevel designed to minimise patient discomfort and tissue trauma. The gently tapered catheter tip aligns with the needle bevel to provide smooth transition from needle to catheter during insertion.

- 14G x 2" (51mm) TER6SROX1451C
- 16G x 2" (51mm) TER6SROX1651C
- 18G x 1.25" (32mm) TER6SROX1832C
- 18G x 2" (51mm) TER6SROX1851C
- 20G x 1" (25mm) TER6SROX2025C
- 20G x 1.25" (32mm) TER6SROX2032C
- 20G x 2" (51mm) TER6SROX2051C
- 22G x 1" (25mm) TER6SROX2225C
- 24G x 3/4" (19mm) TER6SROX2419C

Each \$2.50\* 50+ \$1.85\* 100+ \$1.60\*



## WE ARE A BUSINESS FOR GOOD.



When our Customer Service processes 800 flu orders for the quarter, we give clothing to orphaned children in India.

Thanks to the ongoing support of you, our loyal customers, we have provided 270 orphaned children in India with clothing. A clean set of clothes not only provides warmth during cold weather in India, but also protects from the sun through the year.

Lack of proper clothing may also make a child outcasted and face discrimination. But thanks to you we are able to help combat this ongoing problem.

Learn more about our giving impacts at [www.teammed.com.au/esg-giving-back/](http://www.teammed.com.au/esg-giving-back/)



## Insyte IV Catheters



- 14G x 45mm BD381267
- 16G x 30mm BD381254
- 16G x 45mm BD381257
- 18G x 30mm BD381244
- 18G x 48mm BD381247
- 20G x 30mm BD381234
- 20G x 48mm BD381237
- 22G x 25mm BD381223
- 24G x 19mm BD381212

Each \$3.85\*  
50+ \$2.80\*

## Insyte Autoguard Shielded IV Catheters



Designed to reduce needlestick injuries



The unique push-button shielding mechanism releases the spring and allows the needle and flash chamber to quickly retract into the safety barrel. The clinician maintains control of the process by deciding when to activate the push-button shielding mechanism.

Each \$6.45\*  
50+ \$4.80\*

- 16G x 30mm BD381854
- 18G x 30mm BD381844
- 18G x 48mm BD381847
- 20G x 25mm BD381833
- 20G x 30mm BD381834
- 22G x 25mm BD381823
- 24G x 19mm BD381812

## BD Insyte™ Autoguard™ BC Pro Shielded IV Catheter

With Blood Control technology



When you access a vein, you are at risk of a needlestick injury and blood exposure. If you have to apply pressure after removing the needle, the risk of blood exposure increases. Minimize your worries with BD Insyte™ Autoguard™ BC Pro shielded IV catheter with blood control technology.



- 16G x 30mm [BD381054](#)    ● 18G x 48mm [BD381047](#)    ● 20G x 48mm [BD381037](#)
- 16G x 45mm [BD381057](#)    ● 20G x 25mm [BD381033](#)    ● 22G x 25mm [BD381023](#)
- 18G x 30mm [BD381044](#)    ● 20G x 30mm [BD381034](#)    ● 24G x 19mm [BD381012](#)

Box/50 **\$275.00\***  
 2+ **\$245.30\***    Some more popular sizes are sold as eaches - please check unit when ordering

BD Blood Control Technology - stops blood so it doesn't escape the hub. Controlling blood escape means a cleaner area helping to increase patient confidence.

## Jelco® Optiva® IV Catheters

- 14G x 32mm [SMI5048INT](#)    ● 20G x 32mm [SMI5066INT](#)
- 16G x 32mm [SMI5061INT](#)    ● 20G x 50mm [SMI5069INT](#)
- 16G x 50mm [SMI5062INT](#)    ● 22G x 25mm [SMI5060INT](#)
- 18G x 32mm [SMI5065INT](#)    ● 24G x 19mm [SMI5063INT](#)
- 18G x 44mm [SMI5064INT](#)

Each **\$3.60\***    50+ **\$2.80\***



## Q-Syte™ Luer Access Split Septum

Ultimate performance in closed luer access

A simple split-septum penetrated by a standard male luer (syringe tip or administration set) to become the fluid path

[BD385100](#) Each **\$4.65\***    50+ **\$4.30\***



## SmartSite Needle-free Valve

The SmartSite needle-free valve helps protect you against the estimated 600,000 to 1 million needle-stick injuries that occur yearly. Widely used and award winning, the valve has earned its reputation as a proven, reliable system.

[CAR2000E](#) Each **\$3.35\***    10+ **\$3.15\***



## Injection site

Interlink injection site with luer lock adapter



[BXT2N3399K](#) Each **\$3.40**    10+ **\$2.90**

## Safsite

Capped safety connector for infusion systems

The "needle-free" alternative to injection ports; Luer activated valve that is developed as needle-free injection ports in IV applications

[BN415068](#) Each **\$4.15**    50+ **\$3.85**



# Join The Team that Delivers

Australia's fastest growing family-owned healthcare supplies provider.



## Stopcock 3 Way with Luer Lock




[BD394600](#)  
 Each **\$1.80**    50+ **\$1.50**

### Blood bag system

For the collection of 450ml of blood  
CompoFlex® CPDA-1, 450 ml,

Each \$13.20  
6+ \$12.55

AUT111150S



### Blood bag system

For the collection of 450ml of blood  
Blood bag with needle injury protector, CPDA-1

TER8PB1CD456U0Y  
Pkt/8 \$83.40  
2+ \$81.85




### Giving Sets



#### Straight Type Blood/Solution Set with 200 Micron Filter, Interlink

Blood/solution set, single non-vented drip chamber, 200 micron filter, roller clamp, 1 Interlink injection site 15cm from luer lock adapter (capped) with a retractable collar.




BXTFNC2211

Each \$10.60\*  
48+ \$9.65\*

#### Continue-Flo Solution Set, Clearlink

Solution Continue-Flo set, non-vented drip chamber with 15 micron filter, non-return valve, roller clamp, 2 Clearlink luer activated valves 15cm and 200cm from luer lock adapter (capped) with a retractable collar



BXTFNC3120

Each \$14.40\*  
48+ \$13.45\*

### Giving sets

Transfusion set, 220cm length, vented 20 drops/ml chamber, roller regulator, rotating male luer lock collar with needleless access site

MMD020001

Each \$11.65\*  
50+ \$11.20\*




### Infusion set

Infusion Set with 15um Filter Vented Chamber and NeutralSite™

Includes a vented spike attached to a drip chamber with 15µm filter roller regulator. This range provides the option of a NeutralSite™ or injection site. The 220cm kink resistant standard tubing with male luer lock and rotating collar is made from medical grade PVC.

MMDIV010003

Each \$5.95\*  
50+ \$5.55\*



### Burette Set In-Line Port

Infusion Set with 15um Filter Vented Chamber and NeutralSite™

The m|devices Burette Set In-Line Port has a 150mL capacity vented chamber with a floating safety shut off valve. There is also a vented spike, roller regulator and a Swabable Needleless Access Site at one end and In-line Port with click clamp at the other.

MMDIV040001

Each \$5.55\*  
40+ \$5.30\*



Looking for another set?  
Talk to us and we will help you out.