

Sterile Catheter Pack

#06-688P

- 1 x Sterile Field Blue 60cm x 60cm
- 2 x Dressing Tray 24cm x 9.5cm x 3.5cm
- 8 x Non-woven balls
- 2 x Non-woven swabs 7.5 x 7.5cm
- 2 x Forceps blue
- 1 x Draper 50 x 80cm, Fenestration 8cm x 8cm Blue
- 1 x Dressing Towel 40 x 38cm
- 1 x pair gloves blue (large)

MGT06688P

Each \$4.10
5+ \$3.85



Nelaton Catheters

For intermittent self-catheterisation

Nelaton Catheters are intended for intermittent bladder catheterisation. They are made from medical grade thermo-sensitive plastic and offer a perfectly finished distal end for smooth insertion.



Male Nelaton Catheter

FG10	Black	40cm	MMDUR042001	Each	\$1.15*
FG12	White	40cm	MMDUR042002	Each	\$1.15*
FG14	Green	40cm	MMDUR042003	Each	\$1.40*
FG16	Orange	40cm	MMDUR042004	Each	\$1.15*

Female Nelaton Catheter

FG08	Light Blue	20cm	MMDUR041000	Each	\$1.15*
FG10	Black	20cm	MMDUR041001	Each	\$1.15*
FG12	White	20cm	MMDUR041002	Each	\$1.15*
FG14	Green	20cm	MMDUR041003	Each	\$1.15*




2-way Foley Catheters

Approved for both urethral and suprapubic catheterisation, our 2-Way Foley is made from 100% medical grade translucent and flexible silicone. A symmetrical balloon design supports reliable positioning and minimises the potential for ridging.



FG12	5-10CC	White	MMDUR011000	Each	\$13.95*
FG14	5-10CC	Green	MMDUR011001	Each	\$13.95*
FG16	5-10CC	Orange	MMDUR011002	Each	\$13.95*
FG18	5-10CC	Red	MMDUR011003	Each	\$11.05*
FG20	5-10CC	Yellow	MMDUR011004	Each	\$11.05*



Drainage Bag 2000mL T-Tap NRV Sample Port 120cm Clear Tubing Sterile

Closed system drainage to help reduce the risk of infection

This drainage bag is intended to collect urine via an indwelling catheter or leg bag. It is made from durable material and comes with a shelf life of 5 years. It can store up to 2000ml.

- Clear volume graduation
- With T-Tap
- Varying tube lengths to support patient's needs
- Tubing length: 120cm

MMDUR062003

Each \$4.40*
50+ \$4.15*



2000 ml A4 closed System Bag with Kombikon



NO SET-UP IS TOO BIG OR TOO SMALL

Choose our experienced team for your Practice Set-Up. Call us to discuss your set-up requirements. We can help throughout the entire process.





CARBOMER FREE

Surgi-Gel Lubricating Gel

SURGIgel is a sterile, greaseless, transparent, water-based, non-staining, non-irritating lubricating gel.

Recommended for: Urinary catheterisation, pap smear, Endoscopy and cystoscopy.

- Carbomer Free
- High Lubricity
- Low Viscosity



3ml Sachets (sterile)	Box/144	ORISURO3037F	\$19.10*	3+	\$16.90*
50g Tube (non-sterile)	Each	ORISUR01190F	\$14.50*	4+	\$11.80*

Ingredients known as "carbomers" or "carbopol polymers" may be prone to interfere with cervical screening tests. When doing pap smears ensure your lubricant is Carbomer free.

SURGIgel Plus

Water based lubricating gel for general clinical and lubricating needs. SURGIgel PLUS offers high viscosity and adherence for surgical instruments. Recommended for: Ultrasound, ENT procedures, Suppository insertion.

- Sterile
- Low Lubricity
- High Viscosity (thicker consistency)



50g Tube (sterile)	ORISUR3029	\$11.95*	20+	\$10.95*
--------------------	------------	----------	-----	----------



Medical Lubricating Gel

CARBOMER FREE

Water-based lubricant intended for facilitating the insertion of diagnostic therapeutic devices

- Water Soluble
- Non greasy
- Unscented & Dye free
- Extra Gentle with low tonicity
- High viscosity for strong adherence



WHI132182	
500ml Pump	\$24.40

WE ARE A BUSINESS FOR GOOD



When a customer places their first order, we provide **1 day of food for farmers** affected by natural disaster. So far we have provided **11,090 giving impacts** to date.



Learn more about our giving impacts at <https://www.teammed.com.au/esg-giving-back/>



K-Y Lubricating Jelly

Sterile
CARBOMER FREE



JJ1067500	
42g tube	\$18.40*
12+	\$22.50*

L-gel

Lubricating gel non-sterile



MIA283663	
100g tube	\$16.45*
3+	\$13.85*

MIALG1004	
500g pump	\$41.95*
12+	\$40.45*

Obstetric Care Lotion

Chlorhexidine Gluconate 1%

For disinfection or lubrication during gynaecological and obstetric procedures or childbirth.

200ml	XXC022	\$18.90
50g tube	XXC026	\$4.80



Hysterometer

Disposable Uterine Sound

For Dilatation of the endocervix, and measurements of the uterine cavity. Visible graduations offer clear measurements. Soft distal part on a rigid handle.

Hysterometer CH10	3.3mm dia	CTN1510	Box/25	\$145.65
Hysterometer CH12	4.3mm dia	CTN1512	Box/25	\$160.25
Hysterometer CH14	3.3mm - 4.67mm dia	CTN1514	Box/25	\$160.25

Diameter 4.67 mm with a nozzle of progressive dilatation, from CH10 to CH14 over 25 mm. Indicated for the fitting in of the progestogen IUD (intra-uterine contraceptive device). Dilates the endocervix and measures the depth of the uterine cavity.



WE ARE A BUSINESS FOR GOOD



Marking 13 years of hard work with Team Medical this year, our founder Marijke was allocated donation tokens for B1G1 as her last anniversary rolled over. With this she was able to provide **43 days of access to a safe space and stimulating activities for refugee children in Greece.**



This project provides child residents a welcome respite from the dreary camp life.

Learn more about our giving impacts at www.teammed.com.au/esg-giving-back



IUD insertion kit

This procedure Pack is for an IUD insertion procedure. The kit is sterile and disposable, ensuring all components are available for every insertion.

- 2 x Swab Sticks 20cm Paper Shaft
 - 2 x Gallipots 150ml
 - 1 x Kidney Dish 1000ml Blue
 - 1 x Vaginal Speculum Large
 - 1 x Scissor Mayo 23cm Straight
 - 1 x Forcep Tenaculum 25cm
 - 1 x Sponge Holder 19cm Plastic
 - 1 x Uterine Sound (Hysterometer CH14) 4.67mm with Nozzle
- Wrapped in a Blue Plastic Field 75cm x 75cm



DEF3256

Each \$46.55
10+ \$43.00

IUD Littauer scissor

Disposable and sterile long suture scissors with Littauer blades.

This pair of long, suture scissors has been specially designed to cut the threads of an IUD after inserting it.



CTN1268

Box/25 \$242.10

IUD Hook

- Disposable and sterile IUD hook
- Length 260 mm



CTN1299

Box/20 \$229.60

Saunders IUD Removing Hook

- Stainless Steel



LRKB599

\$66.90

Pelican IUD Thread Retriever

The PM-CARE® RETRIEVAL IUD hook is used to retrieve IUD strings that have been lost in the endocervical canal or the uterine cavity, as well to remove the IUD, if necessary.

263mm - sterile



CTNPS3830026

Each \$15.15

Pipelle

For Endometrial Biopsy

MIS250303 Each \$24.80



Gyn & Push

Cervix Biopsy Forceps



CTN01390

Box/25 \$470.00

Baxter
Welch Allyn

Specialty Disposable Instruments

KleenSpec Cordless Illumination System

The Welch Allyn Cordless Illumination System is a unique design with enhanced features including bright, white LED light to help improve patient exams. The compact, efficient system offers on-time of 80 minutes that will get you through the day, and recharges in 6 hours. Standardise your facility with KleenSpec Cordless Illumination Systems and KleenSpec Disposable Specula to help please both your staff and your patients.

- For use with the KleenSpec Disposable Vaginal Specula
- No bulb replacement necessary, for reduced cost of ownership
- Compact LED light source with homogenous spot improves exam area visualisation
- No cord issues and on-time of 80 minutes, to help improve exam efficiency
- Automatic illumination when specula is attached / Once removed light goes off



WA80015 **\$450**

Complete KleenSpec Cordless Illuminator System with Charging Station

KleenSpec Vaginal Specula

For use with illumination system

Small	WA59000	Box/24	\$185.70	4+	\$169.55
Medium	WA59001	Box/24	\$185.70	4+	\$169.55
Large	WA59004	Box/18	\$169.40	4+	\$156.30

With built-in light source

Small	WA59000LEDB	Box/24	\$315.15
Medium	WA59001LEDB	Box/24	\$315.15
Large	WA59004LEDB	Box/18	\$289.65

Sanitary Pads

Stayfree Regular	no wings	XXS044	Pkt/20	\$12.00
Stayfree Super	no wings	XXS051	Pkt/18	\$13.05
Libra Extra	with wings	XXL021	Pkt/14	\$10.15



Pessary Rings

53mm	SMI700300053	90mm	SMI700300090
56mm	SMI700300056	95mm	SMI700300095
59mm	SMI700300059	100mm	SMI700300100
62mm	SMI700300062	110mm	SMI700300110
65mm	SMI700300065		
68mm	SMI700300068		
71mm	SMI700300071		
74mm	SMI700300074		
77mm	SMI700300077		
80mm	SMI700300080		
85mm	SMI700300085		



\$62.75

Coccyx Wedge cushion with cut out

When doing a Speculum Examination, it may be easier to perform a pap smear if the pelvis is raised off the couch or bed. The coccyx wedge will lift the pelvis up to make the examination easier.



ALLCEXC

\$57.30

Condoms

LIFESTYLE Lubricated Condoms AN6369 Pkt/20 \$30.15

LIFESTYLE Lubricated Condoms XXC148

Box/144 \$65.35




SureStep™ hCG

One Step Pregnancy test

- Rapid results in 3-4 minutes
- Detects hCG at a concentration of 25 mIU/mL or greater
- Proven accuracy of >99% - reliable
- Confirmation of correctly performed test

ALE102425

Box/25 \$45.00



TestPack™ +Plus with OBC hCG Urine

Complete confidence for the rapid detection of early pregnancy

- Easy to use
- An indicator tells you when it is ready to read
- Detects hCG in urine to levels of 25mIU/ml and over
- Highly accurate
- Rapid results, in around 5 minutes, allows immediate follow up
- Unique in-built controls give complete confidence throughout testing

ALEIV505798


Box/20 \$176.45



Bioline™ Influenza Antigen Test

Rapid flu test

- Differential detection of influenza virus type A & B
- Kit contains all key materials required to perform a test
- Results in 10-15 minutes
- Simple procedure
- Throat and nasal swab compatible
- Highly sensitive and specific (refer to IFU for details)



ALESD19FK11

Box/10 \$168.35

Panbio COVID-19 Antigen Self-Test

- Test individuals with or without symptoms or suspected of COVID-19 exposure
- Proven performance (98.1% sensitivity, 99.8% specificity v ersus nasal PCR)
- Can be used for all ages
- Patient-friendly nasal swab
- Results in only 15 minutes
- Confirmed performance against known variants
- Pre-filled vials
- Break-off swab minimizes staff exposure

ALE41FK91

Box/4 \$49.25*



Professional Kits also available - contact us for more information.