

TAT-2000C Temporal Artery Thermometer

- Fast non-invasive measurement
- Less than one second response time
- Range: 15.5°C - 42.0°C
- Display Resolution: 0.1°C
- Measurement technique: Scan across forehead
- Highly accurate
- SmartGlow illuminated display
- Store up to 8 temperature readings
- Optional Probe Covers available
- Can be wiped down with alcohol wipes
- 3-year warranty



ECOTAT2000C
\$139

Forehead Thermometer

- Fast non-invasive measurement
- Instant one second response time
- Range: 32.0°C-43.0°C
- Display Resolution: 0.1°C
- Measurement place: Forehead
- Highly accurate
- Display screen: LCD screen with 4-digit and special icons
- 1-year warranty



PACYT1
\$102



Forehead Thermometer

Model MC-720

- Instant, one second reading
- 3 measurement modes (Forehead, surface and room temperature)
- Backlight provides visibility in the dark
- Last reading displayed on the same screen as the current reading for easy comparison
- Last 25 readings memory
- 2-year warranty



DAVMC720
\$99

WELCARE Standard Digital Thermometer WDT404

- 10 second reading time
- Rigid tip



DAVWDT404 \$12.50

WELCARE Deluxe Digital Thermometer WDT505

- 10 second reading time
- Flexible tip



DAVWDT505 \$13.50

Terumo Digital Thermometer

- Oral/Rectal TER1ETC402S \$28.50
- Under Arm TER1ETC202S \$29.85



Welch Allyn® SureTemp® Plus 690 Electronic Thermometer

The SureTemp 690 Thermometer is fast—Readings in 4 to 6 seconds for oral, 10 to 13 seconds for rectal/paediatric axillary and 12 to 15 seconds for adult axillary temperatures. Featuring the technology that sets the standard for accuracy. You can have confidence the SureTemp 690 will deliver accurate, repeatable readings on patients of all age groups.

- Large LCD displays temperatures in Fahrenheit or Celsius
- Convenient storage for 25 probe covers
- Last temperature recall button
- Takes approximately 6,000 readings on 3 AA batteries
- Interchangeable, removable wells to help reduce the risk of cross-contamination

\$325
WA1690410



Probe Covers
WA1454828 Box/250 \$44.30 5+ \$40.50 10+ \$39.80



Welch Allyn®
Braun ThermoScan PRO 6000

Tympanic Infrared Thermometer

The Welch Allyn Braun ThermoScan PRO 6000 ear thermometer featuring advanced measurement technology, has the speed and simplicity you need to capture accurate and reliable temperature readings in a comfortable way for patients

- Measurement time of 2 – 3 seconds
- Pre-heated probe tip and sensor promotes accuracy
- PerfectTemp™ technology compensates for changes in placement of probe
- ExacTemp™ technology ensures stability of probe during measurement
- Automatic power down energy saving
- 60 second pulse timer measures pulse rate and respiration
- Good for 1000 measurements on each set of batteries
- 3-year warranty

3YR
warranty



BRAUN



Thermoscan Probe Covers
Suitable for PRO 6000 and PRO 4000

Box/200 \$31.50
WAPC200

Small cradle - WA6000200

\$325



Genius™ 3 Tympanic Thermometer

Durable, Dependable, and Accurate

- Peak Select technology takes multiple readings and selects the highest temperature for display in 1 to 2 seconds
- No-touch probe cover application and cover ejector for enhanced infection control
- Storage compartment for up to 32 tips in the base of the thermometer
- Highly accurate - +/- 0.3°C
- Low ongoing costs - economically priced probe covers
- Long battery life - 15,000 temperatures (3 x AAA batteries)
- Widely used throughout hospitals
- 3-year warranty

3YR
warranty



COVIDIEN

\$325

CAR303013

Genius Probe Covers
Suitable for Genius 2 and Genius 3

Box/96 COV303030 \$16.85 4+ \$14.85

Hospital Grade Thermometer

OMRON

WET100 Ear Thermometer

2in1 Ear and Object Temperature Detection

- 1 second quick measurement
- Fever alarm with beeper
- Colour Fever Indicator:
Green: Result is Normal
Yellow: Result is Elevated
Red: Result is High
- 30 readings memory
- Real time clock
- Probe cover detection & probe ejection function
- 2-year warranty



DAVWET100

\$89

Welcare Probe Covers

Pack/60 DAVWET10060 \$8.95

Thermometer Probe Covers & Sheaths

Genius 2 Probe Covers Suitable for Genius 2 and Genius 3	COV303030	Box/96	\$16.85	4+	\$14.85
Thermoscan Probe Covers Suitable for PRO4000 & PRO6000	WAPC200	Box/200	\$31.50		
SureTemp Probe Covers	WA1454828	Box/250	\$44.30	5+	\$40.50
Omron TH839S Probe Covers	DAVTH839S40	Box/40	\$24.15		
ABN Probe Covers/Sheaths	BYDI429	Box/100	\$12.05	2+	\$10.95

