

# Power your practice

## Efficiency



Practice management consultancy



Practice administration



Practice supplies



Practice applications

## Protection



Practice advice



Life insurance advice



Practice insurance

## Growth



Practice legal services



Practice finance



Power your practice with a full suite of products and services to optimise your medical practice

[avant.org.au/practicesolutions](http://avant.org.au/practicesolutions) | 1800 010 236

 **Avant**  
Practice  
Solutions

## HEINE HR® 2.5x Binocular Loupe

**High resolution. Crisp image. Various magnifications.**

Your optimal loupe is as individual as your work. Experience the difference of getting less tired - even during long procedures.

- Crisp, high-resolution image
- Large field of view
- Lightweight
- Multi-coated Achromatic Optics
- i-View loupe mount provides any angle of view and flips up the optics independently of the optional LED LoupeLight



HR 2.5x Binocular Loupes on S-Frame	
340mm working distance	<a href="#">HEC00032355</a>
420mm working distance	<a href="#">HEC00032356</a>
520mm working distance	<a href="#">HEC00032567</a>

Contains: Optics HR with i-View on S-Frame, protective Lenses, 2 sterilisable swivel levers, retaining cord, cleaning fluid, cleaning cloth, carrying case



**\$2925**  
RRP \$3075

HR 2.5x Binocular Loupes on S-Frame with LoupeLight 2	
340mm working distance	<a href="#">HEC00832450</a>
420mm working distance	<a href="#">HEC00832451</a>
520mm working distance	<a href="#">HEC00832452</a>

Contains: LoupeLight 2, Detachable Yellow Filter, mPack mini, mPack mini Belt Clip, E4-USB plug-in transformer, protective lenses, 2 sterilisable levers, retaining cord, cleaning fluid, cleaning cloth, carrying case



**\$5820**  
RRP \$5995



For more information download the product brochure from our website



## HEINE LoupeLight 2

**Brighter. Lighter. Improved functionality.**

Most trusted by dentists and ENT specialists: The combination of light and loupe delivers the ultimate in performance

- Shadow free
- Aluminium housing
- i-View loupe mount with flip-up function
- Bright, homogeneous and coaxial illumination (typ. 55,000 lux)
- Optional Polarisation Filter System



LED HQ  
LED NOW IN HEINE QUALITY.



LoupeLight 2 for i-view loupe mount (HR/HRP)			
LoupeLight 2 (light only)	<a href="#">HEC00832226</a>	RRP \$1995.00	<b>\$1840.00</b>
LoupeLight 2 with Mpack mini and E4-USB plug-in transformer	<a href="#">HEC00832265</a>	RRP \$3425.00	<b>\$3255.00</b>
mPack mini Belt Clip	<a href="#">HEX00099007</a>		<b>\$25.00</b>

For more information download the product brochure from our website



For more information download the product brochure from our website

## HEINE ML4 LED Headlight

**Perfect fit. Perfect view. During long procedures.**

Top seller in all surgical specialties: The HEINE ML4 LED HeadLight is specially designed for healthcare professionals who need a perfect view during long surgical operations or treatments. The headband of the ML4 LED HeadLight offers multiple adjustments points and soft padding ensuring highest comfort and a firm fit.

- Perfect fit
- Shadow free
- Adjustable spot size




ML4 LED Headlight			
with mPack and plug-in transformer	<a href="#">HEJ00831413</a>	RRP \$4025.00	<b>\$3820</b>
with mPack UNPLUGGED and plug-in transformer UNPLUGGED	<a href="#">HEJ00831415</a>	RRP \$4225.00	<b>\$4025</b>

**For HEINE Loupe and Headlight Power Options Visit our Website**



### Illuco optics - Professional and Precise

Introducing the Flip-up loupes and wireless headlights from the Optical and LED specialist engineers at Illuco. Renowned for their dermatoscope range, they now bring this sophisticated designed loupe range.



### Illuco Flip-up Loupes 3.5x

- 3 years warranty for oculars, 2 years for the frame and hinge
- Working distance 42cm
- Easily adjustable pupil distance to your preference
- Possible to add prescription lens
- Mount directly on your own frame

Black Frame/Black Loupe  
Silver Frame/Black Loupe  
Pink Frame/White Loupe

ILLIFLB1000  
ILLIFLS1000  
ILLIFLP100035X

**\$1330**



### IHL-2000 LED Wireless Headlight

- Cable Free
- Quick battery change
- Uniform, Vivid Round Spot over run time
- Max 35,000 lux
- 7.5+ Hours Run Time (2.5+ hrs per battery)
- 37g of Lightweight with battery
- Touch sensor - on/off, 3 step bright control
- Fast USB C charging
- Compatible Adapters to fit most manufacturers' loupes



The system works without cables as it is directly coupled to a compact battery unit. For improved hygiene, the power and intensity of the light can be controlled via a touch sensor using a hand gesture. The system is very lightweight so as to reduce pressure on the nose, and it is attached to the loupes using a easy slide-on dovetail connector. The battery casing is available in two elegant colours; its finely ribbed surface structure adds sophistication to the design.

### Wireless Headlight, 2 Battery Packs

Silver [ILLIHL2000](#) **\$1270**  
Gold [ILLIHL2000G](#) **\$1300**



**Ask us about  
other products,  
configurations and  
connectivity  
options in the  
Illuco loupe range**

## MaggyScan Portable LED Magnifier Lamp ML219

Cordless operation, user-tested intuitive design, and 2.5-hour run time, this is the perfect inspection tool to supplement your everyday medical procedures.

- 6,000 lux @ 25cm
- 2x Magnification
- Large 125mm glass lens
- Automatically recharges when placed in wall dock
- 2.5-hour continuous run time

**Ideal for skin inspection**  
**Popular choice for**  
**Dermatologists and Skin Clinics**

RBKML219

**\$810**



Visit our website to watch the product video



**MAGGYLAMP**  
medical

## Illuco Magnum IDS-3100 UV Light w/ Charger

The IDS-3100 is a new concept of UV Light Magnifier. Its high-power UV LED and specially designed lens coating present an excellent UV fluorescence image. It has 3 options in UV light and 2 options for white light examination, both polarised and non-polarised.

- 2.5 x Magnification
- 100mm large field of view
- 300g light weight is easy to handle
- Nano Hybrid Hard Coating enhances UV transmittance and lens protection
- Aspherical lens minimizes distortion and improves visual acuity.
- Automatic shut down after 5 minutes idle run

ILLIDS3100

**\$1635**



**ILLUCO**

## Mag-Eyes Hands Free Magnifier

With 2 Lenses

- Ultra light weight 50 grams
- Use with or without glasses
- Easy to use/clean

1.6 / 2.0 x magnification **RBKML400** **\$85**

2.25 / 2.75 x magnification **RBKML410** **\$105**



## Guthrie Loupe

- Adjustable head band padded for user comfort
- Plastic lens 2.7 x magnification
- Retractable visor
- Precision lens

**LRM2001** **\$125**



### Standard penlight

- High intensity light beam
- Pocket clip on/off press switch
- Battery powered



BYDS231  
\$11.20

### Deluxe penlight

- High intensity light beam
- Pocket clip doubles as on/off switch
- Useful 5cm rule and dilation guide
- Battery powered



BYDS233  
\$30

### Disposable Penlight Torch

- High intensity light beam
- Pocket clip doubles as on/off switch
- Useful 5cm rule and dilation guide
- Battery powered



LRDDT

Each \$7.15  
6+ \$6.30

### Ultra Violet Light

- Convenient hand size
- Takes 4 x AA size batteries (not included)

GPIDL01 \$50



### Mini Maglite® LED



- Powerful projecting LED beam
- 4 Selectable Modes: High Power, Low Power (25%), Blink Mode, SOS
- Powered by two (2) AA alkaline batteries
- 127 Lumens
- Battery powered



\$65  
Black SHE89210

### Wood's Lamp

- 4 x high power long life UV tubes (4 watts)
- Bulb protection cover included
- 5 Diopter magnification (2.25x)
- Large magnifying lens (5 x 10cm)
- Includes light blocking drape

RBKML215  
\$590



### ACTIK®



- Several lighting modes suitable for different situations, from proximity lighting to rapid movement
- Boost mode for temporary access to maximum brightness
- Lighting colours: white and red
- Shines up to 80 metres
- 350 Lumens (max mode)
- Watertightness: IPX4 (weather resistant)
- Battery powered
- Lightweight only 86g

Ideal for the doctor's bag or in any emergency situation



\$65  
SPEE99AAC

### TIKKINA®



- Three lighting modes: proximity, movement and distance vision
- Shines up to 55 metres (maximum mode)
- Lightweight only 81g
- 250 Lumen brightness
- Battery powered
- Watertightness: IPX4 (weather resistant)



Hybrid



Hybrid

\$45  
Blue SPE91ABC  
Red SPE91ABB